



QUALITY ASSURANCE PROCEDURE SOLUTION TEST REPORT

BATCH REPORT: 14005

CUSTOMER INFORMATION

Washington State Patrol – Breath Test Program
811 East Roanoke SEATTLE, WA 98102

TESTING PROCEDURE USED: TLD Technical Manual, Chapter 4.0 Certification of Simulator Solutions; Headspace-Gas Chromatography.

TESTING ITEM INFORMATION

TARGET VAPOR CONCENTRATION: 0.08 g/210L
DATE PREPARED: 02/12/2014
BATCH UNITS: g/100mL

IDENTITY: QAP Solution
PREPARED BY: Amanda Chandler

	AC	AG	JLK
1	0.101	0.098	0.101
2	0.099	0.099	0.099
3	0.100	0.099	0.099
4	0.100	0.099	0.099
5	0.099	0.099	0.099
C	0.100	0.100	0.100

ETHANOL CONTROL INFORMATION

LOT NUMBER: A095230 EXPIRATION: 09/2017 CONCENTRATION: 0.10 g/100mL

RESULTS OF TESTING

AVERAGE SOLUTION CONCENTRATION: 0.0993 g/100mL PRECISION CV (%): 0.82
STANDARD DEVIATION: 0.00082 NUMBER OF TESTS: 15

EQUIVALENT VAPOR CONCENTRATION: 0.0808 g/210L
COMBINED STANDARD UNCERTAINTY: ± 0.0009 (k=1, 68% confidence interval)

WASHINGTON STATE PATROL – TOXICOLOGY LABORATORY DIVISION




Lisa Noble Forensic Scientist Supervisor

3/5/14

DATE REPORT ISSUED

THIS TESTING WAS PERFORMED BY:

ANALYST	NAME	SIGNATURE	DATE TESTED
AC	Amanda Chandler		02/12/2014
AG	Andrew Gingras		02/13/2014
JLK	Justin L. Knoy		02/18/2014

This report applies only to the item being tested and shall not be reproduced except in full, without the written approval of the WSP Toxicology Laboratory Division. Page 1 of 1

Washington State Patrol - Toxicology Laboratory Division
QAP Test Report Calculation Record

QAP Solution Batch #: 14005

Date Prepared: 2/12/2014

Analyst:	AC	AG	JLK
Date Tested:	2/12/2014	2/13/2014	2/18/2014
Instrument:	HSGC #3	HSGC #3	HSGC #3
1	0.101	0.098	0.101
2	0.099	0.099	0.099
3	0.100	0.099	0.099
4	0.100	0.099	0.099
5	0.099	0.099	0.099
C	0.100	0.100	0.100

CV^2_{COA}	$CV^2_{QAP\ Solution}$	$CV^2_{Control}$	$CV^2_{Part\ Coef}$
0.0000084100	0.0000045043	0.0000000000	0.0001016326

Ethanol Control Lot #: A095230

Control Uncertainty (%): 0.29

Average Solution Concentration: 0.0993 g/100mL
 Standard Deviation: 0.00082 g/100mL
 Precision CV (%): 0.82
 Equivalent Vapor Concentration: 0.0808 g/210L
 Combined Standard Uncertainty (\pm): 0.0009 g/210L

Calculations performed by: Lisa Noble [Signature] 2/26/14
 Name Signature Date

Calculations verified by: Amanda M. Black [Signature] 3-4-2014
 Name Signature Date

Method: Hand calculation

Tech. review performed by: Lisa Noble [Signature] 2/26/14
 Name Signature Date

SIMULATOR SOLUTION DATA ENTRY REVIEW

Reviewer/s: Amanda M. Black Date: 3-4-2014

Location: WSP-FLSB Seattle, WA Solution Batch Number: 14005

	YES	NO	N/A
Analysis dates do not precede preparation date:	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Declarations signed and properly dated:	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Data entry corresponds to all chromatograms:	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
All signatures present on Test Report:	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Average solution concentration correct:	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Standard deviation correct:	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
CV (%) correct:	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Equivalent vapor concentration correct:	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
All chromatograms and sequences included in file:	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Ethanol control information present: (lot # present & used within expiration)	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Complies with accuracy and precision requirements established by the State Toxicologist:	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

Comments:

Reviewer Signature:  Date: 3-4-2014

Reviewer Signature: N/A as 3-4-14 Date: _____

SOLUTION CERTIFICATE REVIEW

Please check that the data on your chromatograms is the data entered into the Test Report, that the date to the right of your name is the date that you tested the solution, and then sign the Test Report.

Please initial and date below to affirm that you have:

- 1) Checked your data
- 2) Checked the date to the right of your name on the Test Report
- 3) Signed the Test Report

	Initials	Date
Amanda Chandler	AC	2/26/14
Andrew Gingras	AG	2/26/14
Asa Louis		
Brittany Ball		
Christie Mitchell-Mata		
Christopher Johnston		
Dawn Sklerov		
Justin Knoy	JK	2/27/14
Katie Knorr		
Lyndsey Lowe		
Naziha Nuwayhid		
Rebecca Fiaherty		
Sarah Swenson		

Batch # 14005

JAY INSLEE
Governor



JOHN R. BATISTE
Chief

STATE OF WASHINGTON
WASHINGTON STATE PATROL
WASHINGTON STATE TOXICOLOGY LABORATORY

2203 Airport Way South, Suite 360 • Seattle, Washington 98134-2927 • (206) 262-6100 • FAX (206) 262-6145

**DATAMASTER 0.08 QAP SOLUTION
CERTIFICATION FOR LOT 14005**

I, Amanda Chandler, do certify under penalty of perjury that:

I am employed by the Washington State Toxicology Laboratory, and a part of my responsibilities includes preparing and testing the alcohol solutions for the DataMaster breath test instrument.

I possess the following qualifications: BS degree in Forensic Science.

The qap solution, Lot Number 14005, was prepared in the Washington State Toxicology Laboratory on 2/12/2014. I examined and tested this solution. It was found to conform to those standards established by the state toxicologist for the certification of simulator solution. It shall not be used to perform a quality assurance procedure after 2/12/2015.

Seattle, WA

Amanda Chandler *2/24/14*

Amanda Chandler
Forensic Toxicologist

Date



JAY INSLEE
Governor



JOHN R. BATISTE
Chief

STATE OF WASHINGTON
WASHINGTON STATE PATROL
WASHINGTON STATE TOXICOLOGY LABORATORY

2203 Airport Way South, Suite 360 • Seattle, Washington 98134-2927 • (206) 262-6100 • FAX (206) 262-6145

**DATAMASTER 0.08 QAP SOLUTION
CERTIFICATION FOR LOT 14005**

I, Andrew Gingras, do certify under penalty of perjury that:

I am employed by the Washington State Toxicology Laboratory, and a part of my responsibilities includes preparing and testing the alcohol solutions for the DataMaster breath test instrument.

I possess the following qualifications: BS degree in Cell and Molecular Biology and MS degree in Forensic Science.

The gap solution, Lot Number 14005, was prepared in the Washington State Toxicology Laboratory on 2/12/2014. I examined and tested this solution. It was found to conform to those standards established by the state toxicologist for the certification of simulator solution. It shall not be used to perform a quality assurance procedure after 2/12/2015.

Seattle, WA

A handwritten signature in blue ink that reads "Andrew Gingras" followed by a date "2/26/14" written in the same ink.

Andrew Gingras
Forensic Toxicologist

Date



JAY INSLEE
Governor



JOHN R. BATISTE
Chief

STATE OF WASHINGTON
WASHINGTON STATE PATROL
WASHINGTON STATE TOXICOLOGY LABORATORY

2203 Airport Way South, Suite 360 • Seattle, Washington 98134-2927 • (206) 262-6100 • FAX (206) 262-6145

**DATAMASTER 0.08 QAP SOLUTION
CERTIFICATION FOR LOT 14005**


I, Justin L. Knoy, do certify under penalty of perjury that:

I am employed by the Washington State Toxicology Laboratory, and a part of my responsibilities includes preparing and testing the alcohol solutions for the DataMaster breath test instrument.

I possess the following qualifications: BS degree in Biology, and MS degree in Forensic Science.

The gap solution, Lot Number 14005, was prepared in the Washington State Toxicology Laboratory on 2/12/2014. I examined and tested this solution. It was found to conform to those standards established by the state toxicologist for the certification of simulator solution. It shall not be used to perform a quality assurance procedure after 2/12/2015.

Seattle, WA

 2-27-14
Justin L. Knoy Date
Forensic Toxicologist

FILE A COPY IN THE BATCH FILE FOR EACH SOLUTION LISTED ON THE WORKSHEET

Preparation Date: 2/12/2014 Initials of Preparer: AC

Expiration Date: 2/12/2015

Lot # of 200-proof Ethanol used in preparation: 2CB0070

Date the 200-proof Ethanol bottle was opened: 2/12/2014

After opening, each bottle of 200-proof Ethanol is approved for use for 6 months unless an extension is approved by the State Toxicologist. This timeframe applies to the 200-proof Ethanol only, not to simulator solutions which have a 1 year expiration.

Simulator Solution	Volume of Ethanol (mL)	Volume of Deionized Water (L)		Batch Number
QAP 0.04	11.2	18	<input checked="" type="checkbox"/>	<u>2/12/14 14004</u>
QAP 0.08	22.4	18	<input checked="" type="checkbox"/>	<u>14005</u>
QAP 0.10	28.1	18	<input checked="" type="checkbox"/>	<u>14006</u>
QAP 0.15	42.0	18	<input checked="" type="checkbox"/>	<u>14007</u>
ESS	66.5	52	<input type="checkbox"/>	<u> </u>

Stir bar is rotating

Stirred for minimum 30 minutes; 2 hours for ESS

Spigot purged

Aliquot taken

Batch labeled, packaged and sealed

2/12/2014
Date

If different ethanol lot numbers are used in the preparation of solutions, record them and the corresponding solution batch numbers in the comments section.

Comments:

Amanda Chandler
Analyst Signature

2/12/2014
Date

Sequence Parameters:

Operator: Amanda Chandler
 Data File Naming: Prefix/Counter
 Signal 1 Prefix: SIG1
 Counter: 0001
 Signal 2 Prefix: SIG2
 Counter: 0001
 Data Directory: C:\HPCHEM\2\DATA\
 Data Subdirectory: 140212AC
 Part of Methods to run: According to Runtime Checklist
 Barcode Reader: not used
 Shutdown Cmd/Macro: none

Sequence Comment:

Cal 1 (0.079 g/100mL) - Lot#E1213-01 - Exp 06/03/2014
 Cal 2 (0.158 g/100mL) - Lot#E1213-02 - Exp 06/03/2014
 Cal 3 (0.316 g/100mL) - Lot#E1213-03 - Exp 06/03/2014

CTRL 1 (0.04 g/100mL) - Lot#FN120110-04 - Exp 12/2015
 CTRL 2 (0.10 g/100mL) - Lot#~~A095232~~ - Exp 09/2017 *A095236 ac 2/12/14*
 CTRL 3 (0.20 g/100mL) - Lot#FN100511-01 - Exp 10/2016

n-Propanol ISTD - Lot# P0114 - Exp 4/26/2014

Sequence Table (Front Injector):

Method and Injection Info Part:


Line	Location	SampleName	Method	Inj	SampleType	InjVolume	DataFile
1	Vial 1	BLANK	SIMALC3	1	Sample		
2	Vial 2	0.079 CAL 1	SIMALC3	1	Calib		
3	Vial 3	0.158 CAL 2	SIMALC3	1	Calib		
4	Vial 4	0.316 CAL 3	SIMALC3	1	Calib		
5	Vial 5	Neg Control-AC	SIMALC3	1	Ctrl Samp		
6	Vial 6	0.04 Control-AC	SIMALC3	1	Ctrl Samp		
7	Vial 7	0.10 Control-AC	SIMALC3	1	Ctrl Samp		
8	Vial 8	0.20 Control-AC	SIMALC3	1	Ctrl Samp		
9	Vial 9	Neg Control-AC	SIMALC3	1	Ctrl Samp		
10	Vial 10	QAP 14004 0.04 1	SIMALC3	1	Sample		
11	Vial 11	QAP 14004 0.04 2	SIMALC3	1	Sample		
12	Vial 12	QAP 14004 0.04 3	SIMALC3	1	Sample		
13	Vial 13	QAP 14004 0.04 4	SIMALC3	1	Sample		
14	Vial 14	QAP 14004 0.04 5	SIMALC3	1	Sample		
15	Vial 15	0.10 Control-AC	SIMALC3	1	Ctrl Samp		
16	Vial 16	Neg Control-AC	SIMALC3	1	Ctrl Samp		
17	Vial 17	QAP 14005 0.08 1	SIMALC3	1	Sample		
18	Vial 18	QAP 14005 0.08 2	SIMALC3	1	Sample		
19	Vial 19	QAP 14005 0.08 3	SIMALC3	1	Sample		
20	Vial 20	QAP 14005 0.08 4	SIMALC3	1	Sample		
21	Vial 21	QAP 14005 0.08 5	SIMALC3	1	Sample		
22	Vial 22	0.10 Control-AC	SIMALC3	1	Ctrl Samp		
23	Vial 23	Neg Control-AC	SIMALC3	1	Ctrl Samp		
24	Vial 24	QAP 14006 0.10 1	SIMALC3	1	Sample		
25	Vial 25	QAP 14006 0.10 2	SIMALC3	1	Sample		
26	Vial 26	QAP 14006 0.10 3	SIMALC3	1	Sample		

14007


ac

Washington State Patrol Toxicology Laboratory
2203 Airport Way S Seattle, WA 98134


Inj. Date: 2/12/2014 1:35:21 PM Sample Name: QAP 14005 0.08 1
Instrument: HSGC#3 Operator: Amanda Chandler
Column: DB-ALC2 Location: Vial 17
Method: C:\HPCHEM\2\METHODS\SIMALC3.M
Sample Info:



#	Compound	Peak Area	RT (min)
1	Ethanol	801	1.018
2	n-Propanol	1658	1.744



Ethanol 0.101 g/100mL




n-Propanol 0.012 g/100mL


ac

Washington State Patrol Toxicology Laboratory
2203 Airport Way S Seattle, WA 98134


Inj. Date: 2/12/2014 1:38:28 PM Sample Name: QAP 14005 0.08 2
Instrument: HSGC#3 Operator: Amanda Chandler
Column: DB-ALC2 Location: Vial 18
Method: C:\HPCHEM\2\METHODS\SIMALC3.M
Sample Info:



#	Compound	Peak Area	RT (min)
1	Ethanol	680	1.018
2	n-Propanol	1435	1.742



Ethanol 0.099 g/100mL




n-Propanol 0.012 g/100mL


AC

Washington State Patrol Toxicology Laboratory
 2203 Airport Way S Seattle, WA 98134


Inj. Date: 2/12/2014 1:41:35 PM Sample Name: QAP 14005 0.08 3
 Instrument: HSGC#3 Operator: Amanda Chandler
 Column: DB-ALC2 Location: Vial 19
 Method: C:\HPCHEM\2\METHODS\SIMALC3.M
 Sample Info:



#	Compound	Peak Area	RT (min)
1	Ethanol	705	1.018
2	n-Propanol	1478	1.743



Ethanol 0.100 g/100mL



n-Propanol 0.012 g/100mL

ac


Washington State Patrol Toxicology Laboratory
2203 Airport Way S Seattle, WA 98134

Inj. Date: 2/12/2014 1:44:43 PM
Instrument: HSGC#3
Column: DB-ALC2


Sample Name: QAP 14005 0.08 4
Operator: Amanda Chandler
Location: Vial 20

Method: C:\HPCHEM\2\METHODS\SIMALC3.M


Sample Info:



#	Compound	Peak Area	RT (min)
1	Ethanol	698	1.018
2	n-Propanol	1464	1.742



Ethanol 0.100 g/100mL



n-Propanol 0.012 g/100mL


Washington State Patrol Toxicology Laboratory
2203 Airport Way S Seattle, WA 98134

Inj. Date: 2/12/2014 1:47:51 PM
Instrument: HSGC#3
Column: DB-ALC2


Sample Name: QAP 14005 0.08 5
Operator: Amanda Chandler
Location: Vial 21

Method: C:\HPCHEM\2\METHODS\SIMALC3.M


Sample Info:



#	Compound	Peak Area	RT (min)
1	Ethanol	688	1.018
2	n-Propanol	1455	1.742




Ethanol 0.099 g/100mL




n-Propanol 0.012 g/100mL

Washington State Patrol Toxicology Laboratory
 2203 Airport Way S Seattle, WA 98134

Inj. Date: 2/12/2014 1:50:58 PM Sample Name: 0.10 Control-AC
 Instrument: HSGC#3 Operator: Amanda Chandler
 Column: DB-ALC2 Location: Vial 22
 Method: C:\HPCHEM\2\METHODS\SIMALC3.M
 Sample Info:




#	Compound	Peak Area	RT (min)
1	Ethanol	698	1.018
2	n-Propanol	1460	1.742



Ethanol 0.100 g/100mL

14005




n-Propanol 0.012 g/100mL


ac

Washington State Patrol Toxicology Laboratory
 2203 Airport Way S Seattle, WA 98134


Inj. Date: 2/12/2014 1:54:06 PM Sample Name: Neg Control-AC
 Instrument: HSGC#3 Operator: Amanda Chandler
 Column: DB-ALC2 Location: Vial 23
 Method: C:\HPCHEM\2\METHODS\SIMALC3.M
 Sample Info:



#	Compound	Peak Area	RT (min)
1	Ethanol	0	0.000
2	n-Propanol	1598	1.743



Ethanol 0.000 g/100mL



n-Propanol 0.012 g/100mL

14005

AK

Sequence Parameters:

Operator: Andrew Gingras
 Data File Naming: Prefix/Counter
 Signal 1 Prefix: SIG1
 Counter: 0001
 Signal 2 Prefix: SIG2
 Counter: 0001
 Data Directory: C:\HPCHEM\2\DATA\
 Data Subdirectory: 140213AG
 Part of Methods to run: According to Runtime Checklist
 Barcode Reader: not used
 Shutdown Cmd/Macro: none

Sequence Comment:

Cal 1 (0.079 g/100mL) - Lot#E1213-01 - Exp 6/3/2014
 Cal 2 (0.158 g/100mL) - Lot#E1213-02 - Exp 6/3/2014
 Cal 3 (0.316 g/100mL) - Lot#E1213-03 - Exp 6/3/2014

 CTRL 1 (0.04 g/100mL) - Lot#FN120110-04 - Exp 12/2015
 CTRL 2 (0.10 g/100mL) - Lot#A095230 - Exp 9/2017
 CTRL 3 (0.20 g/100mL) - Lot#~~FN100511-01~~ - Exp 10/2016
Cal 2/13/14
 n-Propanol ISTD - Lot# P0114 - Exp 4/26/2014

Sequence Table (Front Injector):

Method and Injection Info Part:

Line	Location	SampleName	Method	Inj	SampleType	InjVolume	DataFile
1	Vial 1	BLANK	SIMALC3	1	Sample		
2	Vial 2	0.079 CAL 1	SIMALC3	1	Calib		
3	Vial 3	0.158 CAL 2	SIMALC3	1	Calib		
4	Vial 4	0.316 CAL 3	SIMALC3	1	Calib		
5	Vial 5	Neg Control	SIMALC3	1	Ctrl Samp		
6	Vial 6	0.04 Control	SIMALC3	1	Ctrl Samp		
7	Vial 7	0.10 Control	SIMALC3	1	Ctrl Samp		
8	Vial 8	0.20 Control	SIMALC3	1	Ctrl Samp		
9	Vial 9	Neg Control	SIMALC3	1	Ctrl Samp		
10	Vial 10	QAP 14004 0.04 1	SIMALC3	1	Sample		
11	Vial 11	QAP 14004 0.04 2	SIMALC3	1	Sample		
12	Vial 12	QAP 14004 0.04 3	SIMALC3	1	Sample		
13	Vial 13	QAP 14004 0.04 4	SIMALC3	1	Sample		
14	Vial 14	QAP 14004 0.04 5	SIMALC3	1	Sample		
15	Vial 15	0.10 Control	SIMALC3	1	Ctrl Samp		
16	Vial 16	Neg Control	SIMALC3	1	Ctrl Samp		
17	Vial 17	QAP 14005 0.08 1	SIMALC3	1	Sample		
18	Vial 18	QAP 14005 0.08 2	SIMALC3	1	Sample		
19	Vial 19	QAP 14005 0.08 3	SIMALC3	1	Sample		
20	Vial 20	QAP 14005 0.08 4	SIMALC3	1	Sample		
21	Vial 21	QAP 14005 0.08 5	SIMALC3	1	Sample		
22	Vial 22	0.10 Control	SIMALC3	1	Ctrl Samp		
23	Vial 23	Neg Control	SIMALC3	1	Ctrl Samp		
24	Vial 24	QAP 14006 0.10 1	SIMALC3	1	Sample		
25	Vial 25	QAP 14006 0.10 2	SIMALC3	1	Sample		
26	Vial 26	QAP 14006 0.10 3	SIMALC3	1	Sample		

14007

Washington State Patrol Toxicology Laboratory
2203 Airport Way S Seattle, WA 98134

Inj. Date: 2/13/2014 3:43:54 PM

Sample Name: QAP 14005 0.08 1

Instrument: HSGC#3


Operator: Andrew Gingras

Column: DB-ALC2


Location: Vial 17

Method: C:\HPCHEM\2\METHODS\SIMALC3.M


Sample Info:



#	Compound	Peak Area	RT (min)
1	Ethanol	683	1.017
2	n-Propanol	1450	1.740



Ethanol 0.098 g/100mL




n-Propanol 0.012 g/100mL

Washington State Patrol Toxicology Laboratory
 2203 Airport Way S Seattle, WA 98134


Inj. Date: 2/13/2014 3:47:01 PM
 Instrument: HSGC#3
 Column: DB-ALC2
 Method: C:\HPCHEM\2\METHODS\SIMALC3.M

Sample Name: QAP 14005 0.08 2
 Operator: Andrew Gingras
 Location: Vial 18


Sample Info:



#	Compound	Peak Area	RT (min)
1	Ethanol	696	1.017
2	n-Propanol	1470	1.740




Ethanol 0.099 g/100mL




n-Propanol 0.012 g/100mL

Washington State Patrol Toxicology Laboratory
2203 Airport Way S Seattle, WA 98134


Inj. Date: 2/13/2014 3:50:09 PM Sample Name: QAP 14005 0.08 3
Instrument: HSGC#3 Operator: Andrew Gingras
Column: DB-ALC2 Location: Vial 19
Method: C:\HPCHEM\2\METHODS\SIMALC3.M
Sample Info:



#	Compound	Peak Area	RT (min)
1	Ethanol	692	1.016
2	n-Propanol	1460	1.740




Ethanol 0.099 g/100mL




n-Propanol 0.012 g/100mL

Washington State Patrol Toxicology Laboratory
 2203 Airport Way S Seattle, WA 98134


Inj. Date: 2/13/2014 3:53:16 PM Sample Name: QAP 14005 0.08 4
 Instrument: HSGC#3 Operator: Andrew Gingras
 Column: DB-ALC2 Location: Vial 20
 Method: C:\HPCHEM\2\METHODS\SIMALC3.M
 Sample Info:



#	Compound	Peak Area	RT (min)
1	Ethanol	684	1.017
2	n-Propanol	1446	1.740




Ethanol 0.099 g/100mL




n-Propanol 0.012 g/100mL

Washington State Patrol Toxicology Laboratory
2203 Airport Way S Seattle, WA 98134


Inj. Date: 2/13/2014 3:56:24 PM Sample Name: QAP 14005 0.08 5
Instrument: HSGC#3 Operator: Andrew Gingras
Column: DB-ALC2 Location: Vial 21
Method: C:\HPCHEM\2\METHODS\SIMALC3.M
Sample Info:



#	Compound	Peak Area	RT (min)
1	Ethanol	692	1.016
2	n-Propanol	1461	1.740



Ethanol 0.099 g/100mL




n-Propanol 0.012 g/100mL


Washington State Patrol Toxicology Laboratory
2203 Airport Way S Seattle, WA 98134

Inj. Date: 2/13/2014 3:59:31 PM
Instrument: HSGC#3
Column: DB-ALC2
Method: C:\HPCHEM\2\METHODS\SIMALC3.M
Sample Info:


Sample Name: 0.10 Control
Operator: Andrew Gingras
Location: Vial 22



#	Compound	Peak Area	RT (min)
1	Ethanol	706	1.017
2	n-Propanol	1472	1.740



Ethanol 0.100 g/100mL



n-Propanol 0.012 g/100mL


14005

Washington State Patrol Toxicology Laboratory
2203 Airport Way S Seattle, WA 98134


Inj. Date: 2/13/2014 4:02:38 PM
Instrument: HSGC#3
Column: DB-ALC2
Method: C:\HPCHEM\2\METHODS\SIMALC3.M

Sample Name: Neg Control
Operator: Andrew Gingras
Location: Vial 23


Sample Info:



#	Compound	Peak Area	RT (min)
1	Ethanol	0	0.000
2	n-Propanol	1475	1.741



Ethanol 0.000 g/100mL



n-Propanol 0.012 g/100mL

14005

Sequence Parameters:

Operator: Justin Knoy

Data File Naming: Prefix/Counter

Signal 1 Prefix: SIG1

Counter: 0001

Signal 2 Prefix: SIG2

Counter: 0001

Data Directory: C:\HPCHEM\2\DATA\

Data Subdirectory: 140218JU

Part of Methods to run: According to Runtime Checklist

Barcode Reader: not used

Shutdown Cmd/Macro: none

Sequence Comment:

Ethanol Calibrator 1, E1213-01 - Exp. 06/03/2014

Ethanol Calibrator 2, E1213-02 - Exp. 06/03/2014

Ethanol Calibrator 3, E1213-03 - Exp. 06/03/2014

CTRL1 (0.04g/100mL), Lot # FN120110-04 - Exp. 12/2015

CTRL2 (0.10g/100mL), Lot # A095230 - Exp. 09/2017

CTRL3 (0.20g/100mL), Lot # FN100511-01 - Exp. 10/2016

Internal Standard Lot#P0114 - Exp. 4/26/14

Sequence Table (Front Injector):

Method and Injection Info Part:

Line	Location	SampleName	Method	Inj	SampleType	InjVolume	DataFile
1	Vial 1	BLANK	SIMALC3	1	Sample		
2	Vial 2	CAL1 0.079	SIMALC3	1	Calib		
3	Vial 3	CAL2 0.158	SIMALC3	1	Calib		
4	Vial 4	CAL3 0.316	SIMALC3	1	Calib		
5	Vial 5	NEG CTRL	SIMALC3	1	Ctrl Samp		
6	Vial 6	CTRL1 (0.04)	SIMALC3	1	Ctrl Samp		
7	Vial 7	CTRL2 (0.10)	SIMALC3	1	Ctrl Samp		
8	Vial 8	CTRL3 (0.20)	SIMALC3	1	Ctrl Samp		
9	Vial 9	NEG CTRL	SIMALC3	1	Ctrl Samp		
10	Vial 10	14004-1	SIMALC3	1	Sample		
11	Vial 11	14004-2	SIMALC3	1	Sample		
12	Vial 12	14004-3	SIMALC3	1	Sample		
13	Vial 13	14004-4	SIMALC3	1	Sample		
14	Vial 14	14004-5	SIMALC3	1	Sample		
15	Vial 15	CTRL2 (0.10)	SIMALC3	1	Ctrl Samp		
16	Vial 16	NEG CTRL	SIMALC3	1	Ctrl Samp		
17	Vial 17	14005-1	SIMALC3	1	Sample		
18	Vial 18	14005-2	SIMALC3	1	Sample		
19	Vial 19	14005-3	SIMALC3	1	Sample		
20	Vial 20	14005-4	SIMALC3	1	Sample		
21	Vial 21	14005-5	SIMALC3	1	Sample		
22	Vial 22	CTRL2 (0.10)	SIMALC3	1	Ctrl Samp		
23	Vial 23	NEG CTRL	SIMALC3	1	Ctrl Samp		
24	Vial 24	14006-1	SIMALC3	1	Sample		
25	Vial 25	14006-2	SIMALC3	1	Sample		
26	Vial 26	14006-3	SIMALC3	1	Sample		
27	Vial 27	14006-4	SIMALC3	1	Sample		

14007

*Calibration data
w/14004
LK 2-19-14*

IK

Line	Location	SampleName	Method	Inj	SampleType	InjVolume	DataFile
28	Vial 28	14006-5	SIMALC3	1	Sample		
29	Vial 29	CTRL2 (0.10)	SIMALC3	1	Ctrl Samp		
30	Vial 30	NEG CTRL	SIMALC3	1	Ctrl Samp		
31	Vial 31	14007-1	SIMALC3	1	Sample		
32	Vial 32	14007-2	SIMALC3	1	Sample		
33	Vial 33	14007-3	SIMALC3	1	Sample		
34	Vial 34	14007-4	SIMALC3	1	Sample		
35	Vial 35	14007-5	SIMALC3	1	Sample		
36	Vial 36	CTRL2 (0.10)	SIMALC3	1	Ctrl Samp		
37	Vial 37	NEG CTRL	SIMALC3	1	Ctrl Samp		

Calibration Part:

Line	Location	SampleName	Method	CalLev	Update	RF Update	RT Interval
2	Vial 2	CAL1 0.079	SIMALC3	1	Replace	Replace	
3	Vial 3	CAL2 0.158	SIMALC3	2	Replace	Replace	
4	Vial 4	CAL3 0.316	SIMALC3	3	Replace	Replace	

Sequence Table (Back Injector):

No entries - empty table!

14004
14005


NK

Washington State Patrol Toxicology Laboratory
 2203 Airport Way S Seattle, WA 98134


Inj. Date: 2/18/2014 2:18:23 PM
 Instrument: HSGC#3
 Column: DB-ALC2
 Method: C:\HPCHEM\2\METHODS\SIMALC3.M

Sample Name: 14005-1
 Operator: Justin Knoy
 Location: Vial 17


Sample Info:



#	Compound	Peak Area	RT (min)
1	Ethanol	721	1.018
2	n-Propanol	1507	1.744



Ethanol 0.101 g/100mL




n-Propanol 0.012 g/100mL


Handwritten signature/initials

Washington State Patrol Toxicology Laboratory
2203 Airport Way S Seattle, WA 98134


Inj. Date: 2/18/2014 2:21:31 PM Sample Name: 14005-2
Instrument: HSGC#3 Operator: Justin Knoy
Column: DB-ALC2 Location: Vial 18
Method: C:\HPCHEM\2\METHODS\SIMALC3.M
Sample Info:



#	Compound	Peak Area	RT (min)
1	Ethanol	698	1.016
2	n-Propanol	1488	1.741



Ethanol 0.099 g/100mL




n-Propanol 0.012 g/100mL


R

Washington State Patrol Toxicology Laboratory
 2203 Airport Way S Seattle, WA 98134


Inj. Date: 2/18/2014 2:24:38 PM Sample Name: 14005-3
 Instrument: HSGC#3 Operator: Justin Knoy
 Column: DB-ALC2 Location: Vial 19
 Method: C:\HPCHEM\2\METHODS\SIMALC3.M
 Sample Info:



#	Compound	Peak Area	RT (min)
1	Ethanol	706	1.017
2	n-Propanol	1511	1.741



Ethanol 0.099 g/100mL




n-Propanol 0.012 g/100mL


1X

Washington State Patrol Toxicology Laboratory
2203 Airport Way S Seattle, WA 98134


Inj. Date: 2/18/2014 2:27:46 PM Sample Name: 14005-4
Instrument: HSGC#3 Operator: Justin Knoy
Column: DB-ALC2 Location: Vial 20
Method: C:\HPCHEM\2\METHODS\SIMALC3.M
Sample Info:



#	Compound	Peak Area	RT (min)
1	Ethanol	697	1.016
2	n-Propanol	1496	1.741



Ethanol 0.099 g/100mL




n-Propanol 0.012 g/100mL


AK

Washington State Patrol Toxicology Laboratory
 2203 Airport Way S Seattle, WA 98134


Inj. Date: 2/18/2014 2:30:53 PM Sample Name: 14005-5
 Instrument: HSGC#3 Operator: Justin Knoy
 Column: DB-ALC2 Location: Vial 21
 Method: C:\HPCHEM\2\METHODS\SIMALC3.M
 Sample Info:



#	Compound	Peak Area	RT (min)
1	Ethanol	701	1.017
2	n-Propanol	1501	1.741



Ethanol 0.099 g/100mL




n-Propanol 0.012 g/100mL


Handwritten signature or initials

Washington State Patrol Toxicology Laboratory
 2203 Airport Way S Seattle, WA 98134


Inj. Date: 2/18/2014 2:34:01 PM Sample Name: CTRL2 (0.10)
 Instrument: HSGC#3 Operator: Justin Knoy
 Column: DB-ALC2 Location: Vial 22
 Method: C:\HPCHEM\2\METHODS\SIMALC3.M
 Sample Info: 0.10g/100mL



#	Compound	Peak Area	RT (min)
1	Ethanol	699	1.016
2	n-Propanol	1483	1.741



Ethanol 0.100 g/100mL




n-Propanol 0.012 g/100mL

14005


Handwritten signature

Washington State Patrol Toxicology Laboratory
 2203 Airport Way S Seattle, WA 98134


Inj. Date: 2/18/2014 2:37:08 PM Sample Name: NEG CTRL
 Instrument: HSGC#3 Operator: Justin Knoy
 Column: DB-ALC2 Location: Vial 23
 Method: C:\HPCHEM\2\METHODS\SIMALC3.M
 Sample Info:



#	Compound	Peak Area	RT (min)
1	Ethanol	0	0.000
2	n-Propanol	1484	1.743



Ethanol 0.000 g/100mL



n-Propanol 0.012 g/100mL

14005

N