



QUALITY ASSURANCE PROCEDURE SOLUTION TEST REPORT

BATCH REPORT: 14052

CUSTOMER INFORMATION

Washington State Patrol – Breath Test Program
811 East Roanoke SEATTLE, WA 98102

TESTING PROCEDURE USED: TLD Technical Manual, Chapter 4.0 Certification of Simulator Solutions; Headspace-Gas Chromatography.

TESTING ITEM INFORMATION

TARGET VAPOR CONCENTRATION: 0.20 g/210L
DATE PREPARED: 10/01/2014
BATCH UNITS: g/100mL

IDENTITY: QAP Solution
PREPARED BY: Katie Knorr

	KK	JLK	RF
1	0.253	0.253	0.250
2	0.254	0.253	0.251
3	0.251	0.253	0.251
4	0.253	0.253	0.251
5	0.252	0.254	0.253
C	0.103	0.102	0.102

ETHANOL CONTROL INFORMATION

LOT NUMBER: FN08051301 EXPIRATION: 10/2018 CONCENTRATION: 0.10 g/100mL

RESULTS OF TESTING

AVERAGE SOLUTION CONCENTRATION: 0.2523 g/100mL PRECISION CV (%): 0.49
STANDARD DEVIATION: 0.00123 NUMBER OF TESTS: 15

EQUIVALENT VAPOR CONCENTRATION: 0.2051 g/210L
EXPANDED UNCERTAINTY: ± 0.0046 (k=2, 95.45% confidence interval)

WASHINGTON STATE PATROL – TOXICOLOGY LABORATORY DIVISION






Lisa Noble Forensic Scientist Supervisor

10/14/14

DATE REPORT ISSUED

THIS TESTING WAS PERFORMED BY:

ANALYST	NAME	SIGNATURE	DATE TESTED
KK	Katie Knorr		10/01/2014
JLK	Justin L. Knoy		10/09/2014
RF	Rebecca Flaherty		10/10/2014

Washington State Patrol - Toxicology Laboratory Division
QAP Test Report Calculation Record

QAP Solution Batch #: 14052

Date Prepared: 10/1/2014

Analyst:	KK	JLK	RF
Date Tested:	10/1/2014	10/9/2014	10/10/2014
Instrument:	HSGC #1	HSGC #1	HSGC #1
1	0.253	0.253	0.250
2	0.254	0.253	0.251
3	0.251	0.253	0.251
4	0.253	0.253	0.251
5	0.252	0.254	0.253
C	0.103	0.102	0.102

CV^2_{COA}	$CV^2_{QAP\ Solution}$	$CV^2_{Control}$	$CV^2_{Part\ Coef}$
0.0000084100	0.0000015955	0.0000106102	0.0001016326

Ethanol Control Lot #: FN08051301
Control Uncertainty (%): 0.29

Average Solution Concentration: 0.2523 g/100mL
Standard Deviation: 0.00123 g/100mL
Precision CV (%): 0.49
Equivalent Vapor Concentration: 0.2051 g/210L
Combined Standard Uncertainty (\pm): 0.0023 g/210L
Expanded Uncertainty (\pm): 0.0046 coverage factor (k) = 2 (95.45% level of confidence)

Calculations performed by: Lisa Noble [Signature] 10/13/14
Name Signature Date

Calculations verified by: Amanda M. Black [Signature] 10-13-2014 Method: Hand calculation
Name Signature Date

Tech. review performed by: Lisa Noble [Signature] 10/13/14
Name Signature Date

8

SIMULATOR SOLUTION DATA ENTRY REVIEW

Reviewer/s: Amanda H. Black

Date: 10-13-2014

Location: WSP-FLSB Seattle, WA

Solution Batch Number: 14052

	YES	NO	N/A
Analysis dates do not precede preparation date:	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Declarations signed and properly dated:	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Data entry corresponds to all chromatograms:	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
All signatures present on Test Report:	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Average solution concentration correct:	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Standard deviation correct:	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
CV (%) correct:	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Equivalent vapor concentration correct:	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
All chromatograms and sequences included in file:	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Ethanol control information present: (lot # present & used within expiration)	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Complies with accuracy and precision requirements established by the State Toxicologist:	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

Comments:

Reviewer Signature: 

Date: 10-13-2014

Reviewer Signature: AB 10-13-14

Date: _____



SOLUTION CERTIFICATE REVIEW

Please check that the data on your chromatograms is the data entered into the Test Report, that the date to the right of your name is the date that you tested the solution, and then sign the Test Report.

Please initial and date below to affirm that you have:

- 1) Checked your data
- 2) Checked the date to the right of your name on the Test Report
- 3) Signed the Test Report

	Initials	Date
Amanda Chandler		
Andrew Gingras		
Asa Louis		
Brittany Ball		
Christie Mitchell-Mata		
Christopher Johnston		
Dawn Sklerov		
Justin Knoy	JK	10.13.14
Katie Knorr	KK	10/13/14
Lyndsey Lowe		
Naziha Nuwayhid		
Rebecca Flaherty	RF	10/13/14

Batch # 14052 for 10/13/14



JAY INSLEE
Governor



JOHN R. BATISTE
Chief

STATE OF WASHINGTON
WASHINGTON STATE PATROL
WASHINGTON STATE TOXICOLOGY LABORATORY

2203 Airport Way South, Suite 360 • Seattle, Washington 98134-2927 • (206) 262-6100 • FAX (206) 262-6145

**0.20 g/210 L QUALITY ASSURANCE PROCEDURE SOLUTION
CERTIFICATION FOR LOT 14052**

I, Katie Knorr, do certify under penalty of perjury that:

I am employed by the Washington State Toxicology Laboratory, and my responsibilities include the preparation and certification of alcohol solutions for use with evidential breath test instruments.

I possess the following qualifications: BS in Biochemistry and MS degree in Forensic Science.

The quality assurance procedure (QAP) solution, Lot Number 14052, was prepared in the Washington State Toxicology Laboratory on 10/1/2014. I tested this solution and it was found to conform to those standards established by the State Toxicologist for the certification of simulator solution. It shall not be used to perform a quality assurance procedure after 10/1/2015.

Seattle, WA

Katie Knorr 10/13/14

Katie Knorr

Date

Forensic Toxicologist



JAY INSLEE
Governor



JOHN R. BATISTE
Chief

STATE OF WASHINGTON
WASHINGTON STATE PATROL
WASHINGTON STATE TOXICOLOGY LABORATORY

2203 Airport Way South, Suite 360 • Seattle, Washington 98134-2927 • (206) 262-6100 • FAX (206) 262-6145

**0.20 g/210 L QUALITY ASSURANCE PROCEDURE SOLUTION
CERTIFICATION FOR LOT 14052**


I, Justin L. Knoy, do certify under penalty of perjury that:

I am employed by the Washington State Toxicology Laboratory, and my responsibilities include the preparation and certification of alcohol solutions for use with evidential breath test instruments.

I possess the following qualifications: BS degree in Biology, and MS degree in Forensic Science.

The quality assurance procedure (QAP) solution, Lot Number 14052, was prepared in the Washington State Toxicology Laboratory on 10/1/2014. I tested this solution and it was found to conform to those standards established by the State Toxicologist for the certification of simulator solution. It shall not be used to perform a quality assurance procedure after 10/1/2015.

Seattle, WA

 Justin L. Knoy 10.13.14
Justin L. Knoy Date
Forensic Toxicologist

JAY INSLEE
Governor



JOHN R. BATISTE
Chief

STATE OF WASHINGTON
WASHINGTON STATE PATROL
WASHINGTON STATE TOXICOLOGY LABORATORY

2203 Airport Way South, Suite 360 • Seattle, Washington 98134-2927 • (206) 262-6100 • FAX (206) 262-6145

**0.20 g/210 L QUALITY ASSURANCE PROCEDURE SOLUTION
CERTIFICATION FOR LOT 14052**

I, Rebecca Flaherty, do certify under penalty of perjury that:

I am employed by the Washington State Toxicology Laboratory, and my responsibilities include the preparation and certification of alcohol solutions for use with evidential breath test instruments.

I possess the following qualifications: BS degrees in Biochemistry and Psychobiology and MS degree in Forensic Science.

The quality assurance procedure (QAP) solution, Lot Number 14052, was prepared in the Washington State Toxicology Laboratory on 10/1/2014. I tested this solution and it was found to conform to those standards established by the State Toxicologist for the certification of simulator solution. It shall not be used to perform a quality assurance procedure after 10/1/2015.

Seattle, WA

 _____ 10/13/2014

Rebecca Flaherty
Forensic Toxicologist

Date



FILE A COPY IN THE BATCH FILE FOR EACH SOLUTION LISTED ON THE WORKSHEET

Preparation Date: 10/1/14 Expiration Date: 10/1/15 Initials of Preparer: KK

Lot # of 200-proof Ethanol used in preparation: 2EK0002

Date the 200-proof Ethanol bottle was opened: 9/29/14

After opening, each bottle of 200-proof Ethanol is approved for use for 6 months unless an extension is approved by the State Toxicologist. This timeframe applies to the 200-proof Ethanol only, not to simulator solutions which have a 1 year expiration.

Environmental conditions verified as acceptable:

Simulator Solution	Volume of Ethanol (mL)	Volume of Deionized Water (L)		Batch Number
QAP 0.04	11.2	18	<input checked="" type="checkbox"/>	<u>14049</u>
QAP 0.08	22.4	18	<input checked="" type="checkbox"/>	<u>14050</u>
QAP 0.10	28.1	18	<input type="checkbox"/>	<u> </u>
QAP 0.15	42.1	18	<input checked="" type="checkbox"/>	<u>14051</u>
QAP 0.20	56.1	18	<input checked="" type="checkbox"/>	<u>14052</u>
ESS	66.5	52	<input type="checkbox"/>	<u> </u>

Stir bar is rotating

Stirred for minimum 30 minutes; 2 hours for ESS

Spigot purged

Aliquot taken

Batch labeled, packaged and sealed

10/1/14
Date

If different ethanol lot numbers are used in the preparation of solutions, record them and the corresponding solution batch numbers in the comments section.

Comments:

Katie Knowlton
Analyst Signature

10/1/14 JK
Date

Sequence Parameters:

Operator: Katie Knorr
 Data File Naming: Prefix/Counter
 Signal 1 Prefix: SIG1
 Counter: 0001
 Signal 2 Prefix: SIG2
 Counter: 0001
 Data Directory: C:\HPCHEM\1\DATA\
 Data Subdirectory: 141001KK
 Part of Methods to run: According to Runtime Checklist
 Barcode Reader: not used
 Shutdown Cmd/Macro: none

Sequence Comment:

CAL 1 (0.079g/100mL) - LOT# E0814-01 - EXP 2/19/2015
 CAL 2 (0.158g/100mL) - LOT# E0814-02 - EXP 2/19/2015
 CAL 3 (0.316g/100mL) - LOT# E0814-03 - EXP 2/19/2015
 n-Propanol ISTD - LOT# P0814 - 10/30/2014
 CTRL 1 (0.04g/100mL) - LOT# FN05011301 - EXP 05/2018
 CTRL 2 (0.10g/100mL) - LOT# FN08051301 - EXP 10/2018
 CTRL 3 (0.20g/100mL) - LOT# FN03211401 - EXP 06/2019

KK 10/13/14

Calibration filed with batch 14049 10/13/14

14052

10/13/14

Sequence Table (Front Injector):

Method and Injection Info Part:

Line	Location	SampleName	Method	Inj	SampleType	InjVolume	DataFile
1	Vial 1	BLANK	SIMALC1	1	Sample		
2	Vial 2	0.079 CAL 1	SIMALC1	1	Calib		
3	Vial 3	0.158 CAL 2	SIMALC1	1	Calib		
4	Vial 4	0.316 CAL 3	SIMALC1	1	Calib		
5	Vial 5	NEG CTRL	SIMALC1	1	Ctrl Samp		
6	Vial 6	0.04 CTRL	SIMALC1	1	Ctrl Samp		
7	Vial 7	0.10 CTRL	SIMALC1	1	Ctrl Samp		
8	Vial 8	0.20 CTRL	SIMALC1	1	Ctrl Samp		
9	Vial 9	NEG CTRL	SIMALC1	1	Ctrl Samp		
10	Vial 10	14049-1	SIMALC1	1	Sample		
11	Vial 11	14049-2	SIMALC1	1	Sample		
12	Vial 12	14049-3	SIMALC1	1	Sample		
13	Vial 13	14049-4	SIMALC1	1	Sample		
14	Vial 14	14-049-5 14049-5	SIMALC1	1	Sample		
15	Vial 15	0.10 CTRL <i>KK</i>	SIMALC1	1	Ctrl Samp		
16	Vial 16	NEG CTRL <i>10/12/14</i>	SIMALC1	1	Ctrl Samp		
17	Vial 17	14050-1	SIMALC1	1	Sample		
18	Vial 18	14050-2 <i>KK</i>	SIMALC1	1	Sample		
19	Vial 19	14050-3 <i>10/13/14</i>	SIMALC1	1	Sample		
20	Vial 20	14050-4	SIMALC1	1	Sample		
21	Vial 21	14050-5	SIMALC1	1	Sample		
22	Vial 22	0.10 CTRL	SIMALC1	1	Ctrl Samp		
23	Vial 23	NEG CTRL	SIMALC1	1	Ctrl Samp		
24	Vial 24	14051-1	SIMALC1	1	Sample		
25	Vial 25	14051-2	SIMALC1	1	Sample		
26	Vial 26	14051-3	SIMALC1	1	Sample		
27	Vial 27	14051-4	SIMALC1	1	Sample		
28	Vial 28	14051-5	SIMALC1	1	Sample		

JK

KK

Line	Location	SampleName	Method	Inj	SampleType	InjVolume	DataFile
29	Vial 29	0.10CTRL	SIMALC1	1	Ctrl Samp		
30	Vial 30	NEG CTRL	SIMALC1	1	Ctrl Samp		
31	Vial 31	14052-1	SIMALC1	1	Sample		
32	Vial 32	14052-2	SIMALC1	1	Sample		
33	Vial 33	14052-3	SIMALC1	1	Sample		
34	Vial 34	14052-4	SIMALC1	1	Sample		
35	Vial 35	14052-5	SIMALC1	1	Sample		
36	Vial 36	0.10 CTRL	SIMALC1	1	Ctrl Samp		
37	Vial 37	NEG CTRL	SIMALC1	1	Ctrl Samp		

Calibration Part:

Line	Location	SampleName	Method	CalLev	Update	RF	Update	RT	Interval
2	Vial 2	0.079 CAL 1	SIMALC1	1	Replace		Replace		
3	Vial 3	0.158 CAL 2	SIMALC1	2	Replace		Replace		
4	Vial 4	0.316 CAL 3	SIMALC1	3	Replace		Replace		

Sequence Table (Back Injector):

No entries - empty table!

14052
10/13/14

fr

KK

Washington State Patrol Toxicology Laboratory
 2203 Airport Way S Seattle, WA 98134

Inj. Date: 10/1/2014 2:13:39 PM

Sample Name: 14052-1

Instrument: HSGC#1


Operator: Katie Knorr

Column: DB-ALC1


Location: Vial 31

Method: C:\HPCHEM\1\METHODS\SIMALC1.M


Sample Info:



#	Compound	Peak Area	RT (min)
1	Ethanol	2748	1.037
2	n-Propanol	2480	1.700



Ethanol 0.253 g/100mL



n-Propanol 0.012 g/100mL

JK

JK

Washington State Patrol Toxicology Laboratory
 2203 Airport Way S Seattle, WA 98134

Inj. Date: 10/1/2014 2:16:52 PM

Sample Name: 14052-2

Instrument: HSGC#1


Operator: Katie Knorr

Column: DB-ALC1


Location: Vial 32

Method: C:\HPCHEM\1\METHODS\SIMALC1.M


Sample Info:



#	Compound	Peak Area	RT (min)
1	Ethanol	2776	1.036
2	n-Propanol	2490	1.701



Ethanol 0.254 g/100mL



n-Propanol 0.012 g/100mL

Handwritten signature

Handwritten initials 'KK'

Washington State Patrol Toxicology Laboratory
 2203 Airport Way S Seattle, WA 98134

Inj. Date: 10/1/2014 2:20:06 PM

Sample Name: 14052-3

Instrument: HSGC#1


Operator: Katie Knorr

Column: DB-ALC1


Location: Vial 33

Method: C:\HPCHEM\1\METHODS\SIMALC1.M


Sample Info:



#	Compound	Peak Area	RT (min)
1	Ethanol	2708	1.036
2	n-Propanol	2463	1.700



Ethanol 0.251 g/100mL



n-Propanol 0.012 g/100mL

fn

PK

Washington State Patrol Toxicology Laboratory
 2203 Airport Way S Seattle, WA 98134

Inj. Date: 10/1/2014 2:23:20 PM

Sample Name: 14052-4

Instrument: HSGC#1


Operator: Katie Knorr

Column: DB-ALC1


Location: Vial 34

Method: C:\HPCHEM\1\METHODS\SIMALC1.M


Sample Info:



#	Compound	Peak Area	RT (min)
1	Ethanol	2755	1.037
2	n-Propanol	2485	1.700



Ethanol 0.253 g/100mL




n-Propanol 0.012 g/100mL


fr

KK


Inj. Date: 10/1/2014 2:26:32 PM Sample Name: 14052-5
 Instrument: HSGC#1 Operator: Katie Knorr
 Column: DB-ALC1 Location: Vial 35
 Method: C:\HPCHEM\1\METHODS\SIMALC1.M
 Sample Info:



#	Compound	Peak Area	RT (min)
1	Ethanol	2752	1.037
2	n-Propanol	2488	1.701



Ethanol 0.252 g/100mL



n-Propanol 0.012 g/100mL

fu

KK

Washington State Patrol Toxicology Laboratory
2203 Airport Way S Seattle, WA 98134

Inj. Date: 10/1/2014 2:29:45 PM

Sample Name: 0.10 CTRL

Instrument: HSGC#1


Operator: Katie Knorr

Column: DB-ALC1


Location: Vial 36

Method: C:\HPCHEM\1\METHODS\SIMALC1.M

Sample Info:




#	Compound	Peak Area	RT (min)
1	Ethanol	1159	1.036
2	n-Propanol	2555	1.700



Ethanol 0.103 g/100mL

14052
2/10/13/14



n-Propanol 0.012 g/100mL

fu

KK

Washington State Patrol Toxicology Laboratory
 2203 Airport Way S Seattle, WA 98134

Inj. Date: 10/1/2014 2:33:00 PM

Sample Name: NEG CTRL

Instrument: HSGC#1


Operator: Katie Knorr

Column: DB-ALC1


Location: Vial 37

Method: C:\HPCHEM\1\METHODS\SIMALC1.M

Sample Info:




#	Compound	Peak Area	RT (min)
1	Ethanol	0	0.000
2	n-Propanol	2597	1.701



Ethanol 0.000 g/100mL

14052
 8/10/13/14



n-Propanol 0.012 g/100mL

fn

KK

Sequence Parameters:

Operator: Justin Knoy

Data File Naming: Prefix/Counter

Signal 1 Prefix: SIG1
Counter: 0001

Signal 2 Prefix: SIG2
Counter: 0001

Data Directory: C:\HPCHEM\1\DATA\

Data Subdirectory: 141009JK

Part of Methods to run: According to Runtime Checklist

Barcode Reader: not used

Shutdown Cmd/Macro: none

Sequence Comment:

Ethanol Calibrator 1, E0814-01 - Exp. 02/19/2015
 Ethanol Calibrator 2, E0814-02 - Exp. 02/19/2015
 Ethanol Calibrator 3, E0814-03 - Exp. 02/19/2015
 CTRL1 (0.04g/100mL), Lot # FN05011301 - Exp. 05/2018
 CTRL2 (0.10g/100mL), Lot # FN08051301 - Exp. 10/2018
 CTRL3 (0.20g/100mL), Lot # FN03211401 - Exp. 06/2019

Internal Standard Lot#P0914 - Exp. 12/29/14

Calibration filed with batch 14049 on 10/13/14

*14052
on 10/13/14*

Sequence Table (Front Injector):

Method and Injection Info Part:

Line	Location	SampleName	Method	Inj	SampleType	InjVolume	DataFile
1	Vial 1	BLANK	SIMALC1	1	Sample		
2	Vial 2	0.079 CAL 1	SIMALC1	1	Calib		
3	Vial 3	0.158 CAL 2	SIMALC1	1	Calib		
4	Vial 4	0.316 CAL 3	SIMALC1	1	Calib		
5	Vial 5	NEG CTRL JK	SIMALC1	1	Ctrl Samp		
6	Vial 6	0.04 CTRL - JK	SIMALC1	1	Ctrl Samp		
7	Vial 7	0.10 CTRL - JK	SIMALC1	1	Ctrl Samp		
8	Vial 8	0.20 CTRL - JK	SIMALC1	1	Ctrl Samp		
9	Vial 9	NEG CTRL JK	SIMALC1	1	Ctrl Samp		
10	Vial 10	14049-1	SIMALC1	1	Sample		
11	Vial 11	14049-2	SIMALC1	1	Sample		
12	Vial 12	14049-3	SIMALC1	1	Sample		
13	Vial 13	14049-4	SIMALC1	1	Sample		
14	Vial 14	14049-5	SIMALC1	1	Sample		
15	Vial 15	0.10 CTRL JK	SIMALC1	1	Ctrl Samp		
16	Vial 16	NEG CTRL JK	SIMALC1	1	Ctrl Samp		
17	Vial 17	14050-1	SIMALC1	1	Sample		
18	Vial 18	14050-2	SIMALC1	1	Sample		
19	Vial 19	14050-3	SIMALC1	1	Sample		
20	Vial 20	14050-4	SIMALC1	1	Sample		
21	Vial 21	14050-5	SIMALC1	1	Sample		
22	Vial 22	0.10 CTRL JK	SIMALC1	1	Ctrl Samp		
23	Vial 23	NEG CTRL JK	SIMALC1	1	Ctrl Samp		
24	Vial 24	14051-1	SIMALC1	1	Sample		
25	Vial 25	14051-2	SIMALC1	1	Sample		
26	Vial 26	14051-3	SIMALC1	1	Sample		
27	Vial 27	14051-4	SIMALC1	1	Sample		

J

JK

Line	Location	SampleName	Method	Inj	SampleType	InjVolume	DataFile
28	Vial 28	14051-5	SIMALC1	1	Sample		
29	Vial 29	0.10 CTRL JK	SIMALC1	1	Ctrl Samp		
30	Vial 30	NEG CTRL JK	SIMALC1	1	Ctrl Samp		
31	Vial 31	14052-1	SIMALC1	1	Sample		
32	Vial 32	14052-2	SIMALC1	1	Sample		
33	Vial 33	14052-3	SIMALC1	1	Sample		
34	Vial 34	14052-4	SIMALC1	1	Sample		
35	Vial 35	14052-5	SIMALC1	1	Sample		
36	Vial 36	0.10 CTRL JK	SIMALC1	1	Ctrl Samp		
37	Vial 37	NEG CTRL JK	SIMALC1	1	Ctrl Samp		

Calibration Part:

Line	Location	SampleName	Method	CalLev	Update	RF	Update	RT	Interval
2	Vial 2	0.079 CAL 1	SIMALC1	1	Replace		Replace		
3	Vial 3	0.158 CAL 2	SIMALC1	2	Replace		Replace		
4	Vial 4	0.316 CAL 3	SIMALC1	3	Replace		Replace		

Sequence Table (Back Injector):

No entries - empty table!

14052

du
10/13/14

du

JK

Inj. Date: 10/9/2014 12:59:31 PM

Sample Name: 14052-1

Instrument: HSGC#1


Operator: Justin Knoy

Column: DB-ALC1


Location: Vial 31

Method: C:\HPCHEM\1\METHODS\SIMALC1.M


Sample Info:



#	Compound	Peak Area	RT (min)
1	Ethanol	2703	1.037
2	n-Propanol	2380	1.700



Ethanol 0.253 g/100mL



n-Propanol 0.012 g/100mL

J

JK

Washington State Patrol Toxicology Laboratory
 2203 Airport Way S Seattle, WA 98134

Inj. Date: 10/9/2014 1:02:44 PM

Sample Name: 14052-2

Instrument: HSGC#1


Operator: Justin Knoy

Column: DB-ALC1


Location: Vial 32

Method: C:\HPCHEM\1\METHODS\SIMALC1.M


Sample Info:



#	Compound	Peak Area	RT (min)
1	Ethanol	2693	1.036
2	n-Propanol	2368	1.700



Ethanol 0.253 g/100mL



n-Propanol 0.012 g/100mL

Handwritten mark

Handwritten mark

Washington State Patrol Toxicology Laboratory
 2203 Airport Way S Seattle, WA 98134

Inj. Date: 10/9/2014 1:05:57 PM

Sample Name: 14052-3

Instrument: HSGC#1


Operator: Justin Knoy

Column: DB-ALC1


Location: Vial 33

Method: C:\HPCHEM\1\METHODS\SIMALC1.M


Sample Info:



#	Compound	Peak Area	RT (min)
1	Ethanol	2701	1.036
2	n-Propanol	2374	1.700



Ethanol 0.253 g/100mL



n-Propanol 0.012 g/100mL

JK

JK

Inj. Date: 10/9/2014 1:09:11 PM

Sample Name: 14052-4

Instrument: HSGC#1


Operator: Justin Knoy

Column: DB-ALC1


Location: Vial 34

Method: C:\HPCHEM\1\METHODS\SIMALC1.M


Sample Info:



#	Compound	Peak Area	RT (min)
1	Ethanol	2699	1.036
2	n-Propanol	2372	1.700



Ethanol 0.253 g/100mL



n-Propanol 0.012 g/100mL

h

JK

Washington State Patrol Toxicology Laboratory
 2203 Airport Way S Seattle, WA 98134

Inj. Date: 10/9/2014 1:12:24 PM

Sample Name: 14052-5

Instrument: HSGC#1


Operator: Justin Knoy

Column: DB-ALC1


Location: Vial 35

Method: C:\HPCHEM\1\METHODS\SIMALC1.M


Sample Info:



#	Compound	Peak Area	RT (min)
1	Ethanol	2748	1.036
2	n-Propanol	2413	1.700



Ethanol 0.254 g/100mL



n-Propanol 0.012 g/100mL

JK

JK

Washington State Patrol Toxicology Laboratory
 2203 Airport Way S Seattle, WA 98134

Inj. Date: 10/9/2014 1:15:37 PM

Sample Name: 0.10 CTRL JK

Instrument: HSGC#1


Operator: Justin Knoy

Column: DB-ALC1


Location: Vial 36

Method: C:\HPCHEM\1\METHODS\SIMALC1.M

Sample Info:




#	Compound	Peak Area	RT (min)
1	Ethanol	1116	1.036
2	n-Propanol	2415	1.700



Ethanol 0.102 g/100mL

14052
 Jh 10/13/14



n-Propanol 0.012 g/100mL

Jh

K

Washington State Patrol Toxicology Laboratory
2203 Airport Way S Seattle, WA 98134

Inj. Date: 10/9/2014 1:18:52 PM

Sample Name: NEG CTRL JK

Instrument: HSGC#1


Operator: Justin Knoy

Column: DB-ALC1


Location: Vial 37

Method: C:\HPCHEM\1\METHODS\SIMALC1.M

Sample Info:




#	Compound	Peak Area	RT (min)
1	Ethanol	0	0.000
2	n-Propanol	2459	1.700



Ethanol 0.000 g/100mL

14052
initially



n-Propanol 0.012 g/100mL

JK

K

Sequence Parameters:

Operator: Rebecca Flaherty
 Data File Naming: Prefix/Counter
 Signal 1 Prefix: SIG1
 Counter: 0001
 Signal 2 Prefix: SIG2
 Counter: 0001
 Data Directory: C:\HPCHEM\1\DATA\
 Data Subdirectory: 141010RF
 Part of Methods to run: According to Runtime Checklist
 Barcode Reader: not used
 Shutdown Cmd/Macro: none

Sequence Comment:

Ethanol Calibrator 1, E0814-01 - Exp. 02/19/2015
 Ethanol Calibrator 2, E0814-02 - Exp. 02/19/2015
 Ethanol Calibrator 3, E0814-03 - Exp. 02/19/2015
 CTRL1 (0.04g/100mL), Lot # FN05011301 - Exp. 05/2018
 CTRL2 (0.10g/100mL), Lot # FN08051301 - Exp. 10/2018
 CTRL3 (0.20g/100mL), Lot # FN03211401 - Exp. 06/2019
 Internal Standard Lot#P0914 - Exp. 12/29/14

Calibration vials 1-9 filed with 14049.

14052

Inj 1314

Sequence Table (Front Injector):

Method and Injection Info Part:

Line	Location	SampleName	Method	Inj	SampleType	InjVolume	DataFile
1	Vial 1	BLANK	SIMALC1	1	Sample		
2	Vial 2	0.079 CAL 1	SIMALC1	1	Calib		
3	Vial 3	0.158 CAL 2	SIMALC1	1	Calib		
4	Vial 4	0.316 CAL 3	SIMALC1	1	Calib		
5	Vial 5	NEG CTRL	SIMALC1	1	Ctrl Samp		
6	Vial 6	0.04 CTRL	SIMALC1	1	Ctrl Samp		
7	Vial 7	0.10 CTRL	SIMALC1	1	Ctrl Samp		
8	Vial 8	0.20 CTRL	SIMALC1	1	Ctrl Samp		
9	Vial 9	NEG CTRL	SIMALC1	1	Ctrl Samp		
10	Vial 10	14049-1	SIMALC1	1	Sample		
11	Vial 11	14049-2	SIMALC1	1	Sample		
12	Vial 12	14049-3	SIMALC1	1	Sample		
13	Vial 13	14049-4	SIMALC1	1	Sample		
14	Vial 14	14049-5	SIMALC1	1	Sample		
15	Vial 15	0.10 CTRL	SIMALC1	1	Ctrl Samp		
16	Vial 16	NEG CTRL	SIMALC1	1	Ctrl Samp		
17	Vial 17	14050-1	SIMALC1	1	Sample		
18	Vial 18	14050-2	SIMALC1	1	Sample		
19	Vial 19	14050-3	SIMALC1	1	Sample		
20	Vial 20	14050-4	SIMALC1	1	Sample		
21	Vial 21	14050-5	SIMALC1	1	Sample		
22	Vial 22	0.10 CTRL	SIMALC1	1	Ctrl Samp		
23	Vial 23	NEG CTRL	SIMALC1	1	Ctrl Samp		
24	Vial 24	14051-1	SIMALC1	1	Sample		
25	Vial 25	14051-2	SIMALC1	1	Sample		
26	Vial 26	14051-3	SIMALC1	1	Sample		

sh

RF

Line	Location	SampleName	Method	Inj	SampleType	InjVolume	DataFile
27	Vial 27	14051-4	SIMALC1	1	Sample		
28	Vial 28	14051-5	SIMALC1	1	Sample		
29	Vial 29	0.10 CTRL	SIMALC1	1	Ctrl Samp		
30	Vial 30	NEG CTRL	SIMALC1	1	Ctrl Samp		
31	Vial 31	14052-1	SIMALC1	1	Sample		
32	Vial 32	14052-2	SIMALC1	1	Sample		
33	Vial 33	14052-3	SIMALC1	1	Sample		
34	Vial 34	14052-4	SIMALC1	1	Sample		
35	Vial 35	14052-5	SIMALC1	1	Sample		
36	Vial 36	0.10 CTRL	SIMALC1	1	Ctrl Samp		
37	Vial 37	NEG CTRL	SIMALC1	1	Ctrl Samp		

Sequence Table (Back Injector):

No entries - empty table!

14052

Analysis

Handwritten signature

RF

Washington State Patrol Toxicology Laboratory
 2203 Airport Way S Seattle, WA 98134

Inj. Date: 10/10/2014 3:30:53 PM

Sample Name: 14052-1

Instrument: HSGC#1


Operator: Rebecca Flaherty

Column: DB-ALC1


Location: Vial 31

Method: C:\HPCHEM\1\METHODS\SIMALC1.M


Sample Info:



#	Compound	Peak Area	RT (min)
1	Ethanol	2694	1.036
2	n-Propanol	2383	1.700



Ethanol 0.250 g/100mL



n-Propanol 0.012 g/100mL


h

RF


Washington State Patrol Toxicology Laboratory
 2203 Airport Way S Seattle, WA 98134

Inj. Date: 10/10/2014 3:34:06 PM Sample Name: 14052-2
 Instrument: HSGC#1 Operator: Rebecca Flaherty
 Column: DB-ALC1 Location: Vial 32
 Method: C:\HPCHEM\1\METHODS\SIMALC1.M


Sample Info:



#	Compound	Peak Area	RT (min)
1	Ethanol	2696	1.036
2	n-Propanol	2382	1.700



Ethanol 0.251 g/100mL




n-Propanol 0.012 g/100mL

h


RF

Washington State Patrol Toxicology Laboratory
 2203 Airport Way S Seattle, WA 98134


Inj. Date: 10/10/2014 3:37:20 PM Sample Name: 14052-3
 Instrument: HSGC#1 Operator: Rebecca Flaherty
 Column: DB-ALC1 Location: Vial 33
 Method: C:\HPCHEM\1\METHODS\SIMALC1.M
 Sample Info:



#	Compound	Peak Area	RT (min)
1	Ethanol	2724	1.036
2	n-Propanol	2402	1.700



Ethanol 0.251 g/100mL




n-Propanol 0.012 g/100mL

RF


RF

Washington State Patrol Toxicology Laboratory
 2203 Airport Way S Seattle, WA 98134


Inj. Date: 10/10/2014 3:40:33 PM Sample Name: 14052-4
 Instrument: HSGC#1 Operator: Rebecca Flaherty
 Column: DB-ALC1 Location: Vial 34
 Method: C:\HPCHEM\1\METHODS\SIMALC1.M
 Sample Info:



#	Compound	Peak Area	RT (min)
1	Ethanol	2681	1.036
2	n-Propanol	2366	1.700



Ethanol 0.251 g/100mL



n-Propanol 0.012 g/100mL

f

RF

Washington State Patrol Toxicology Laboratory
 2203 Airport Way S Seattle, WA 98134

Inj. Date: 10/10/2014 3:43:46 PM

Sample Name: 14052-5

Instrument: HSGC#1


Operator: Rebecca Flaherty

Column: DB-ALC1


Location: Vial 35

Method: C:\HPCHEM\1\METHODS\SIMALC1.M


Sample Info:



#	Compound	Peak Area	RT (min)
1	Ethanol	2712	1.036
2	n-Propanol	2378	1.700



Ethanol 0.253 g/100mL



n-Propanol 0.012 g/100mL

RF

RF

Washington State Patrol Toxicology Laboratory
 2203 Airport Way S Seattle, WA 98134

Inj. Date: 10/10/2014 3:46:59 PM

Sample Name: 0.10 CTRL

Instrument: HSGC#1


Operator: Rebecca Flaherty

Column: DB-ALC1


Location: Vial 36

Method: C:\HPCHEM\1\METHODS\SIMALC1.M


Sample Info: 14052



#	Compound	Peak Area	RT (min)
1	Ethanol	1113	1.037
2	n-Propanol	2393	1.700



Ethanol 0.102 g/100mL



n-Propanol 0.012 g/100mL

RF

RF

Washington State Patrol Toxicology Laboratory
 2203 Airport Way S Seattle, WA 98134

Inj. Date: 10/10/2014 3:50:14 PM

Sample Name: NEG CTRL

Instrument: HSGC#1


Operator: Rebecca Flaherty

Column: DB-ALC1


Location: Vial 37

Method: C:\HPCHEM\1\METHODS\SIMALC1.M


Sample Info: 14052



#	Compound	Peak Area	RT (min)
1	Ethanol	0	0.000
2	n-Propanol	2614	1.701



Ethanol 0.000 g/100mL



n-Propanol 0.012 g/100mL

fr

RF