



QUALITY ASSURANCE PROCEDURE SOLUTION TEST REPORT

BATCH REPORT: 15006

CUSTOMER INFORMATION

Washington State Patrol – Breath Test Program
811 East Roanoke SEATTLE, WA 98102

TESTING PROCEDURE USED: TLD Technical Manual, Chapter 4.0 Certification of Simulator Solutions; Headspace-Gas Chromatography.

TESTING ITEM INFORMATION

TARGET VAPOR CONCENTRATION: 0.10 g/210L
DATE PREPARED: 01/13/2015
BATCH UNITS: g/100mL

IDENTITY: QAP Solution
PREPARED BY: Elizabeth Wehner

	EW	DN	NN
1	0.128	0.127	0.127
2	0.127	0.128	0.126
3	0.127	0.128	0.127
4	0.127	0.128	0.127
5	0.128	0.128	0.127
C	0.104	0.102	0.102

ETHANOL CONTROL INFORMATION

LOT NUMBER: FN08051301 EXPIRATION: 10/2018 CONCENTRATION: 0.10 g/100mL

RESULTS OF TESTING

AVERAGE SOLUTION CONCENTRATION: 0.1273 g/100mL PRECISION CV (%): 0.48
STANDARD DEVIATION: 0.00062 NUMBER OF TESTS: 15

EQUIVALENT VAPOR CONCENTRATION: **0.1035 g/210L**
EXPANDED UNCERTAINTY: ± 0.0026 (k=2, 95.45% confidence interval)

WASHINGTON STATE PATROL – TOXICOLOGY LABORATORY DIVISION






Lisa Noble Forensic Scientist Supervisor

2/4/15

DATE REPORT ISSUED

THIS TESTING WAS PERFORMED BY:

ANALYST	NAME	SIGNATURE	DATE TESTED
EW	Elizabeth Wehner		01/13/2015
DN	David Nguyen		01/14/2015
NN	Naziha Nuwayhid		01/14/2015

Washington State Patrol - Toxicology Laboratory Division
QAP Test Report Calculation Record

QAP Solution Batch #: 15006

Date Prepared: 1/13/2015

Analyst:	EW	DN	NN
Date Tested:	1/13/2015	1/14/2015	1/14/2015
Instrument:	HSGC #3	HSGC #3	HSGC #3
1	0.128	0.127	0.127
2	0.127	0.128	0.126
3	0.127	0.128	0.127
4	0.127	0.128	0.127
5	0.128	0.128	0.127
C	0.104	0.102	0.102

CV^2_{COA}	$CV^2_{QAP\ Solution}$	$CV^2_{Control}$	$CV^2_{Part\ Coef}$
0.0000084100	0.0000015664	0.0000421656	0.0001016326

Ethanol Control Lot #: FN08051301
Control Uncertainty (%): 0.29

Average Solution Concentration: 0.1273 g/100mL
Standard Deviation: 0.00062 g/100mL
Precision CV (%): 0.48
Equivalent Vapor Concentration: 0.1035 g/210L
Combined Standard Uncertainty (\pm): 0.0013 g/210L
Expanded Uncertainty (\pm): 0.0026 coverage factor (k) = 2 (95.45% level of confidence)

Calculations performed by: Lisa Noble [Signature] 1/27/15
Name Signature Date

Calculations verified by: Amanda M. Black [Signature] 2-3-2015
Name Signature Date

Method: Hand calculation

Tech. review performed by: Lisa Noble [Signature] 1/27/15
Name Signature Date

[Signature]

SIMULATOR SOLUTION DATA ENTRY REVIEW

Reviewer/s: Amanda M. Black Date: 2-3-2015

Location: WSP-FLSB Seattle, WA Solution Batch Number: 15004

	YES	NO	N/A
Analysis dates do not precede preparation date:	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Declarations signed and properly dated:	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Data entry corresponds to all chromatograms:	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
All signatures present on Test Report:	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Average solution concentration correct:	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Standard deviation correct:	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
CV (%) correct:	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Equivalent vapor concentration correct:	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
All chromatograms and sequences included in file:	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Ethanol control information present: (lot # present & used within expiration)	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Complies with accuracy and precision requirements established by the State Toxicologist:	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

Comments:

Reviewer Signature:  Date: 2-3-2015

for

SOLUTION CERTIFICATE REVIEW

Please check that the data on your chromatograms is the data entered into the Test Report, that the date to the right of your name is the date that you tested the solution, and then sign the Test Report.

Please initial and date below to affirm that you have:

- 1) Checked your data
- 2) Checked the date to the right of your name on the Test Report
- 3) Signed the Test Report

	Initials	Date
Amanda Chandler		
Andrew Gingras		
Asa Louis		
Brittany Thomas		
Christie Mitchell-Mata		
Christopher Johnston		
David Nguyen	DN	1/28/15
Dawn Sklerov		
Elizabeth Wehner	EW	01/28/15
Justin Knoy		
Katie Harris		
Lyndsey Lowe		
Naziha Nuwayhid	NN	1.28.15
Rebecca Flaherty		

Batch # 15006 Jan 27/15



JAY INSLEE
Governor



JOHN R. BATISTE
Chief

STATE OF WASHINGTON
WASHINGTON STATE PATROL
WASHINGTON STATE TOXICOLOGY LABORATORY

2203 Airport Way South, Suite 360 • Seattle, Washington 98134-2927 • (206) 262-6100 • FAX (206) 262-6145

**0.10 g/210 L QUALITY ASSURANCE PROCEDURE SOLUTION
CERTIFICATION FOR LOT 15006**

I, Elizabeth Wehner, do certify under penalty of perjury that:

I am employed by the Washington State Toxicology Laboratory, and my responsibilities include the preparation and certification of alcohol solutions for use with evidential breath test instruments.

I possess the following qualifications: BS degree in Biochemistry.

The quality assurance procedure (QAP) solution, Lot Number 15006, was prepared in the Washington State Toxicology Laboratory on 1/13/2015. I tested this solution and it was found to conform to those standards established by the State Toxicologist for the certification of simulator solution. It shall not be used to perform a quality assurance procedure after 1/13/2016.

Seattle, WA

Elizabeth Wehner 01/28/15

Elizabeth Wehner

Date

Forensic Scientist

JAY INSLEE
Governor



JOHN R. BATISTE
Chief

STATE OF WASHINGTON
WASHINGTON STATE PATROL
WASHINGTON STATE TOXICOLOGY LABORATORY

2203 Airport Way South, Suite 360 • Seattle, Washington 98134-2927 • (206) 262-6100 • FAX (206) 262-6145

**0.10 g/210 L QUALITY ASSURANCE PROCEDURE SOLUTION
CERTIFICATION FOR LOT 15006**


I, Naziha Nuwayhid, do certify under penalty of perjury that:

I am employed by the Washington State Toxicology Laboratory, and my responsibilities include the preparation and certification of alcohol solutions for use with evidential breath test instruments.

I possess the following qualifications: Bachelor and Masters Degrees in Biology, Ph.D. degree in Basic Medical Science, ten years experience in clinical laboratory sciences, one year in clinical toxicology and more than ten years in forensic toxicology. I am also board certified by the American Board of Clinical Chemistry.

The quality assurance procedure (QAP) solution, Lot Number 15006, was prepared in the Washington State Toxicology Laboratory on 1/13/2015. I tested this solution and it was found to conform to those standards established by the State Toxicologist for the certification of simulator solution. It shall not be used to perform a quality assurance procedure after 1/13/2016.

Seattle, WA



Naziha Nuwayhid Date
Forensic Scientist



FILE A COPY IN THE BATCH FILE FOR EACH SOLUTION LISTED ON THE WORKSHEET

Preparation Date: 01/13/15 Expiration Date: 01/13/16 Initials of Preparer: EW

Lot # of 200-proof Ethanol used in preparation: 2DG0231

Date the 200-proof Ethanol bottle was opened: 12/03/14

After opening, each bottle of 200-proof Ethanol is approved for use for 6 months unless an extension is approved by the State Toxicologist. This timeframe applies to the 200-proof Ethanol only, not to simulator solutions which have a 1 year expiration.

Environmental conditions verified as acceptable:

Simulator Solution	Volume of Ethanol (mL)	Volume of Deionized Water (L)		Batch Number
QAP 0.04	11.2	18	<input checked="" type="checkbox"/>	<u>15004</u>
QAP 0.08	22.4	18	<input checked="" type="checkbox"/>	<u>15005</u>
QAP 0.10	28.1	18	<input checked="" type="checkbox"/>	<u>15006</u>
QAP 0.15	42.1	18	<input checked="" type="checkbox"/>	<u>15007</u>
QAP 0.20	56.1	18	<input checked="" type="checkbox"/>	<u>15008</u>
ESS	66.5	52	<input type="checkbox"/>	<u> </u>

Stir bar is rotating

Stirred for minimum 30 minutes; 2 hours for ESS

Spigot purged

Aliquot taken

Batch labeled, packaged and sealed

01/13/15
Date

If different ethanol lot numbers are used in the preparation of solutions, record them and the corresponding solution batch numbers in the comments section.

Comments:

Elizabeth Wehner
Analyst Signature

01/13/15
Date

[Handwritten mark]

Sequence Parameters:

Operator: Elizabeth Wehner
 Data File Naming: Prefix/Counter
 Signal 1 Prefix: SIG1
 Counter: 0001
 Signal 2 Prefix: SIG2
 Counter: 0001
 Data Directory: C:\HPCHEM\2\DATA\
 Data Subdirectory: 150113EW
 Part of Methods to run: According to Runtime Checklist
 Barcode Reader: not used
 Shutdown Cmd/Macro: none

Sequence Comment:

CAL 1: 0.079 g/100mL - LOT# E0814-01 - EXP 02/19/2015
 CAL 2: 0.158 g/100mL - LOT# E0814-02 - EXP 02/19/2015
 CAL 3: 0.316 g/100mL - LOT# E0814-03 - EXP 02/19/2015

CTRL 1: 0.04 g/100mL - LOT# FN05011301 - EXP 05/2018
 CTRL 2: 0.10 g/100mL - LOT# FN08051301 - EXP 10/2018
 CTRL 3: 0.20 g/100mL - LOT# FN03211401 - EXP 06/2019

n-PROPANOL ISTD - LOT# P1114 - EXP 02/20/2015

Vials# 1-9 are filed with 15004

Sequence Table (Front Injector):

Method and Injection Info Part:

Line	Location	SampleName	Method	Inj	SampleType	InjVolume	DataFile
1	Vial 1	BLANK	SIMALC3	1	Sample		
2	Vial 2	CAL 1 (0.079)	SIMALC3	1	Calib		
3	Vial 3	CAL 2 (0.158)	SIMALC3	1	Calib		
4	Vial 4	CAL 3 (0.316)	SIMALC3	1	Calib		
5	Vial 5	NEG CTRL	SIMALC3	1	Ctrl Samp		
6	Vial 6	CTRL 1 (0.04)	SIMALC3	1	Ctrl Samp		
7	Vial 7	CTRL 2 (0.10)	SIMALC3	1	Ctrl Samp		
8	Vial 8	CTRL 3 (0.20)	SIMALC3	1	Ctrl Samp		
9	Vial 9	NEG CTRL	SIMALC3	1	Ctrl Samp		
10	Vial 10	15004 #1	SIMALC3	1	Sample		
11	Vial 11	15004 #2	SIMALC3	1	Sample		
12	Vial 12	15004 #3	SIMALC3	1	Sample		
13	Vial 13	15004 #4	SIMALC3	1	Sample		
14	Vial 14	15004 #5	SIMALC3	1	Sample		
15	Vial 15	POS CTRL (0.10)	SIMALC3	1	Ctrl Samp		
16	Vial 16	NEG CTRL	SIMALC3	1	Ctrl Samp		
17	Vial 17	15005 #1	SIMALC3	1	Sample		
18	Vial 18	15005 #2	SIMALC3	1	Sample		
19	Vial 19	15005 #3	SIMALC3	1	Sample		
20	Vial 20	15005 #4	SIMALC3	1	Sample		
21	Vial 21	15005 #5	SIMALC3	1	Sample		
22	Vial 22	POS CTRL (0.10)	SIMALC3	1	Ctrl Samp		
23	Vial 23	NEG CTRL	SIMALC3	1	Ctrl Samp		
24	Vial 24	15006 #1	SIMALC3	1	Sample		

15006
 7/12/15
 JW
 EW

Line	Location	SampleName	Method	Inj	SampleType	InjVolume	DataFile
25	Vial 25	15006 #2	SIMALC3	1	Sample		
26	Vial 26	15006 #3	SIMALC3	1	Sample		
27	Vial 27	15006 #4	SIMALC3	1	Sample		
28	Vial 28	15006 #5	SIMALC3	1	Sample		
29	Vial 29	POS CTRL (0.10)	SIMALC3	1	Ctrl Samp		
30	Vial 30	NEG CTRL	SIMALC3	1	Ctrl Samp		
31	Vial 31	15007 #1	SIMALC3	1	Sample		
32	Vial 32	15007 #2	SIMALC3	1	Sample		
33	Vial 33	15007 #3	SIMALC3	1	Sample		
34	Vial 34	15007 #4	SIMALC3	1	Sample		
35	Vial 35	15007 #5	SIMALC3	1	Sample		
36	Vial 36	POS CTRL (0.10)	SIMALC3	1	Ctrl Samp		
37	Vial 37	NEG CTRL	SIMALC3	1	Ctrl Samp		
38	Vial 38	15008 #1	SIMALC3	1	Sample		
39	Vial 39	15008 #2	SIMALC3	1	Sample		
40	Vial 40	15008 #3	SIMALC3	1	Sample		
41	Vial 41	15008 #4	SIMALC3	1	Sample		
42	Vial 42	15008 #5	SIMALC3	1	Sample		
43	Vial 43	POS CTRL (0.10)	SIMALC3	1	Ctrl Samp		
44	Vial 44	NEG CTRL	SIMALC3	1	Ctrl Samp		

Calibration Part:

Line	Location	SampleName	Method	CalLev	Update	RF	Update	RT	Interval
2	Vial 2	CAL 1 (0.079)	SIMALC3	1	Replace		Replace		
3	Vial 3	CAL 2 (0.158)	SIMALC3	2	Replace		Replace		
4	Vial 4	CAL 3 (0.316)	SIMALC3	3	Replace		Replace		

Sequence Table (Back Injector):

No entries - empty table!

15006
jm/27/15

EW

Washington State Patrol Toxicology Laboratory
 2203 Airport Way S Seattle, WA 98134

Inj. Date: 1/13/2015 12:18:20 PM

Sample Name: 15006 #1

Instrument: HSGC#3


Operator: Elizabeth Wehner

Column: DB-ALC2


Location: Vial 24

Method: C:\HPCHEM\2\METHODS\SIMALC3.M


Sample Info:



#	Compound	Peak Area	RT (min)
1	Ethanol	902	1.020
2	n-Propanol	1534	1.747



Ethanol 0.128 g/100mL




n-Propanol 0.012 g/100mL

fw


EW

Washington State Patrol Toxicology Laboratory
2203 Airport Way S Seattle, WA 98134


Inj. Date: 1/13/2015 12:21:34 PM Sample Name: 15006 #2
Instrument: HSGC#3 Operator: Elizabeth Wehner
Column: DB-ALC2 Location: Vial 25
Method: C:\HPCHEM\2\METHODS\SIMALC3.M
Sample Info:



#	Compound	Peak Area	RT (min)
1	Ethanol	903	1.019
2	n-Propanol	1544	1.745



Ethanol 0.127 g/100mL



n-Propanol 0.012 g/100mL

EW

Washington State Patrol Toxicology Laboratory
 2203 Airport Way S Seattle, WA 98134

Inj. Date: 1/13/2015 12:24:47 PM

Sample Name: 15006 #3

Instrument: HSGC#3


Operator: Elizabeth Wehner

Column: DB-ALC2


Location: Vial 26

Method: C:\HPCHEM\2\METHODS\SIMALC3.M


Sample Info:



#	Compound	Peak Area	RT (min)
1	Ethanol	892	1.020
2	n-Propanol	1524	1.746



Ethanol 0.127 g/100mL



n-Propanol 0.012 g/100mL

EW

Washington State Patrol Toxicology Laboratory
 2203 Airport Way S Seattle, WA 98134

Inj. Date: 1/13/2015 12:28:00 PM

Sample Name: 15006 #4

Instrument: HSGC#3


Operator: Elizabeth Wehner

Column: DB-ALC2


Location: Vial 27

Method: C:\HPCHEM\2\METHODS\SIMALC3.M


Sample Info:



#	Compound	Peak Area	RT (min)
1	Ethanol	902	1.019
2	n-Propanol	1542	1.744



Ethanol 0.127 g/100mL



n-Propanol 0.012 g/100mL

EW

EW

Washington State Patrol Toxicology Laboratory
 2203 Airport Way S Seattle, WA 98134

Inj. Date: 1/13/2015 12:31:14 PM

Sample Name: 15006 #5

Instrument: HSGC#3


Operator: Elizabeth Wehner

Column: DB-ALC2


Location: Vial 28

Method: C:\HPCHEM\2\METHODS\SIMALC3.M


Sample Info:



#	Compound	Peak Area	RT (min)
1	Ethanol	948	1.019
2	n-Propanol	1607	1.746



Ethanol 0.128 g/100mL



n-Propanol 0.012 g/100mL

EW

Inj. Date: 1/13/2015 12:34:27 PM

Sample Name: POS CTRL (0.10)

Instrument: HSGC#3

Operator: Elizabeth Wehner


Column: DB-ALC2

Location: Vial 29


Method: C:\HPCHEM\2\METHODS\SIMALC3.M

Sample Info: POS CTRL: 0.10 g/100mL
 15006


->



#	Compound	Peak Area	RT (min)
1	Ethanol	740	1.019
2	n-Propanol	1553	1.745



Ethanol 0.104 g/100mL



n-Propanol 0.012 g/100mL

Handwritten signature

EW

Inj. Date: 1/13/2015 12:37:40 PM

Sample Name: NEG CTRL

Instrument: HSGC#3


Operator: Elizabeth Wehner

Column: DB-ALC2


Location: Vial 30

Method: C:\HPCHEM\2\METHODS\SIMALC3.M


Sample Info: 15006



#	Compound	Peak Area	RT (min)
1	Ethanol	0	0.000
2	n-Propanol	1525	1.746



Ethanol 0.000 g/100mL



n-Propanol 0.012 g/100mL

EW

EW

Sequence Parameters:

Operator: David Nguyen
 Data File Naming: Prefix/Counter
 Signal 1 Prefix: SIG1
 Counter: 0001
 Signal 2 Prefix: SIG2
 Counter: 0001
 Data Directory: C:\HPCHEM\2\DATA\
 Data Subdirectory: 150114DN
 Part of Methods to run: According to Runtime Checklist
 Barcode Reader: not used
 Shutdown Cmd/Macro: none

Sequence Comment:

CAL 1: 0.079 g/100mL - Lot#: E0814-01 Exp. 02/19/2015
 CAL 2: 0.158 g/100mL - Lot#: E0814-02 Exp. 02/19/2015
 CAL 3: 0.316 g/100mL - Lot#: E0814-03 Exp. 02/19/2015

 CTRL 1: 0.04 g/100mL - Lot#: FN05011301 Exp. 05/2018
 CTRL 2: 0.10 g/100mL - Lot#: FN08051301 Exp. 10/2018
 CTRL 3: 0.20 g/100mL - Lot#: FN03211401 Exp. 06/2019

 n-Propanol ISTD - Lot#: P1114 Exp. 02/20/2015

 Calibration vials 1-9 filed with 15004.

Sequence Table (Front Injector):

Method and Injection Info Part:

Line	Location	SampleName	Method	Inj	SampleType	InjVolume	DataFile
1	Vial 1	BLANK	SIMALC3	1	Sample		
2	Vial 2	CAL 1 (0.079)	SIMALC3	1	Calib		
3	Vial 3	CAL 2 (0.158)	SIMALC3	1	Calib		
4	Vial 4	CAL 3 (0.316)	SIMALC3	1	Calib		
5	Vial 5	NEG CTRL	SIMALC3	1	Ctrl Samp		
6	Vial 6	CTRL 1 (0.04)	SIMALC3	1	Ctrl Samp		
7	Vial 7	CTRL 2 (0.10)	SIMALC3	1	Ctrl Samp		
8	Vial 8	CTRL 3 (0.20)	SIMALC3	1	Ctrl Samp		
9	Vial 9	NEG CTRL	SIMALC3	1	Ctrl Samp		
10	Vial 10	15004 #1	SIMALC3	1	Sample		
11	Vial 11	15004 #2	SIMALC3	1	Sample		
12	Vial 12	15004 #3	SIMALC3	1	Sample		
13	Vial 13	15004 #4	SIMALC3	1	Sample		
14	Vial 14	15004 #5	SIMALC3	1	Sample		
15	Vial 15	POS CTRL (0.10)	SIMALC3	1	Ctrl Samp		
16	Vial 16	NEG CTRL	SIMALC3	1	Ctrl Samp		
17	Vial 17	15005 #1	SIMALC3	1	Sample		
18	Vial 18	15005 #2	SIMALC3	1	Sample		
19	Vial 19	15005 #3	SIMALC3	1	Sample		
20	Vial 20	15005 #4	SIMALC3	1	Sample		
21	Vial 21	15005 #5	SIMALC3	1	Sample		
22	Vial 22	POS CTRL (0.10)	SIMALC3	1	Ctrl Samp		
23	Vial 23	NEG CTRL	SIMALC3	1	Ctrl Samp		
24	Vial 24	15006 #1	SIMALC3	1	Sample		

15008
dn/2/15

dn

DN

Line	Location	SampleName	Method	Inj	SampleType	InjVolume	DataFile
25	Vial 25	15006 #2	SIMALC3	1	Sample		
26	Vial 26	15006 #3	SIMALC3	1	Sample		
27	Vial 27	15006 #4	SIMALC3	1	Sample		
28	Vial 28	15006 #5	SIMALC3	1	Sample		
29	Vial 29	POS CTRL (0.10)	SIMALC3	1	Ctrl Samp		
30	Vial 30	NEG CTRL	SIMALC3	1	Ctrl Samp		
31	Vial 31	15007 #1	SIMALC3	1	Sample		
32	Vial 32	15007 #2	SIMALC3	1	Sample		
33	Vial 33	15007 #3	SIMALC3	1	Sample		
34	Vial 34	15007 #4	SIMALC3	1	Sample		
35	Vial 35	15007 #5	SIMALC3	1	Sample		
36	Vial 36	POS CTRL (0.10)	SIMALC3	1	Ctrl Samp		
37	Vial 37	NEG CTRL	SIMALC3	1	Ctrl Samp		
38	Vial 38	15008 #1	SIMALC3	1	Sample		
39	Vial 39	15008 #2	SIMALC3	1	Sample		
40	Vial 40	15008 #3	SIMALC3	1	Sample		
41	Vial 41	15008 #4	SIMALC3	1	Sample		
42	Vial 42	15008 #5	SIMALC3	1	Sample		
43	Vial 43	POS CTRL (0.10)	SIMALC3	1	Ctrl Samp		
44	Vial 44	NEG CTRL	SIMALC3	1	Ctrl Samp		

Calibration Part:

Line	Location	SampleName	Method	CalLev	Update	RF	Update	RT	Interval
2	Vial 2	CAL 1 (0.079)	SIMALC3	1	Replace		Replace		
3	Vial 3	CAL 2 (0.158)	SIMALC3	2	Replace		Replace		
4	Vial 4	CAL 3 (0.316)	SIMALC3	3	Replace		Replace		

Sequence Table (Back Injector):

No entries - empty table!

15006
Jan/27/15

fn

DN

Washington State Patrol Toxicology Laboratory
 2203 Airport Way S Seattle, WA 98134

Inj. Date: 1/14/2015 10:45:47 AM

Sample Name: 15006 #1

Instrument: HSGC#3


Operator: David Nguyen

Column: DB-ALC2


Location: Vial 24

Method: C:\HPCHEM\2\METHODS\SIMALC3.M


Sample Info:



#	Compound	Peak Area	RT (min)
1	Ethanol	892	1.018
2	n-Propanol	1545	1.745



Ethanol . 0.127 g/100mL




n-Propanol 0.012 g/100mL


DN

Washington State Patrol Toxicology Laboratory
 2203 Airport Way S Seattle, WA 98134


Inj. Date: 1/14/2015 10:49:00 AM Sample Name: 15006 #2
 Instrument: HSGC#3 Operator: David Nguyen
 Column: DB-ALC2 Location: Vial 25
 Method: C:\HPCHEM\2\METHODS\SIMALC3.M
 Sample Info:



#	Compound	Peak Area	RT (min)
1	Ethanol	933	1.018
2	n-Propanol	1609	1.744



Ethanol 0.128 g/100mL



n-Propanol 0.012 g/100mL


fn

DZ


Washington State Patrol Toxicology Laboratory
 2203 Airport Way S Seattle, WA 98134

Inj. Date: 1/14/2015 10:52:13 AM Sample Name: 15006 #3
 Instrument: HSGC#3 Operator: David Nguyen
 Column: DB-ALC2 Location: Vial 26
 Method: C:\HPCHEM\2\METHODS\SIMALC3.M


Sample Info:



#	Compound	Peak Area	RT (min)
1	Ethanol	903	1.018
2	n-Propanol	1558	1.745



Ethanol 0.128 g/100mL




n-Propanol 0.012 g/100mL

DN


Washington State Patrol Toxicology Laboratory
 2203 Airport Way S Seattle, WA 98134

Inj. Date: 1/14/2015 10:55:27 AM Sample Name: 15006 #4
 Instrument: HSGC#3 Operator: David Nguyen
 Column: DB-ALC2 Location: Vial 27
 Method: C:\HPCHEM\2\METHODS\SIMALC3.M


Sample Info:



#	Compound	Peak Area	RT (min)
1	Ethanol	903	1.019
2	n-Propanol	1560	1.746



Ethanol 0.128 g/100mL




n-Propanol 0.012 g/100mL

for
DN


Washington State Patrol Toxicology Laboratory
 2203 Airport Way S Seattle, WA 98134

Inj. Date: 1/14/2015 10:58:40 AM Sample Name: 15006 #5
 Instrument: HSGC#3 Operator: David Nguyen
 Column: DB-ALC2 Location: Vial 28
 Method: C:\HPCHEM\2\METHODS\SIMALC3.M


Sample Info:



#	Compound	Peak Area	RT (min)
1	Ethanol	901	1.018
2	n-Propanol	1553	1.745



Ethanol 0.128 g/100mL




n-Propanol 0.012 g/100mL


for
DN

Washington State Patrol Toxicology Laboratory
 2203 Airport Way S Seattle, WA 98134


Inj. Date: 1/14/2015 11:01:53 AM Sample Name: POS CTRL (0.10)
 Instrument: HSGC#3 Operator: David Nguyen
 Column: DB-ALC2 Location: Vial 29
 Method: C:\HPCHEM\2\METHODS\SIMALC3.M
 Sample Info: POS CTRL: 0.10 g/100mL
 15006



#	Compound	Peak Area	RT (min)
1	Ethanol	757	1.018
2	n-Propanol	1628	1.744



Ethanol 0.102 g/100mL




n-Propanol 0.012 g/100mL

DN


DN

Washington State Patrol Toxicology Laboratory
2203 Airport Way S Seattle, WA 98134


Inj. Date: 1/14/2015 11:05:06 AM Sample Name: NEG CTRL
Instrument: HSGC#3 Operator: David Nguyen
Column: DB-ALC2 Location: Vial 30
Method: C:\HPCHEM\2\METHODS\SIMALC3.M
Sample Info: 15006



#	Compound	Peak Area	RT (min)
1	Ethanol	0	0.000
2	n-Propanol	1576	1.745



Ethanol 0.000 g/100mL



n-Propanol 0.012 g/100mL

fr
DN

Sequence Parameters:

Operator: Naziha Nuwayhid, PhD
 Data File Naming: Prefix/Counter
 Signal 1 Prefix: SIG1
 Counter: 0001
 Signal 2 Prefix: SIG2
 Counter: 0001
 Data Directory: C:\HPCHEM\2\DATA\
 Data Subdirectory: 150114NN
 Part of Methods to run: According to Runtime Checklist
 Barcode Reader: not used
 Shutdown Cmd/Macro: none

Sequence Comment:

CAL 1: 0.079 g/100mL - Lot#: E0814-01 Exp. 02/19/2015
 CAL 2: 0.158 g/100mL - Lot#: E0814-02 Exp. 02/19/2015
 CAL 3: 0.316 g/100mL - Lot#: E0814-03 Exp. 02/19/2015

 CTRL 1: 0.04 g/100mL - Lot#: FN05011301 Exp. 05/2018
 CTRL 2: 0.10 g/100mL - Lot#: FN08051301 Exp. 10/2018
 CTRL 3: 0.20 g/100mL - Lot#: FN03211401 Exp. 06/2019

 n-Propanol ISTD - Lot#: P1114 Exp. 02/20/2015

 Calibration vials 1-9 filed with 15004.

Sequence Table (Front Injector):

Method and Injection Info Part:

Line	Location	SampleName	Method	Inj	SampleType	InjVolume	DataFile
1	Vial 1	BLANK	SIMALC3	1	Sample		
2	Vial 2	CAL 1 (0.079)	SIMALC3	1	Calib		
3	Vial 3	CAL 2 (0.158)	SIMALC3	1	Calib		
4	Vial 4	CAL 3 (0.316)	SIMALC3	1	Calib		
5	Vial 5	NEG CTRL	SIMALC3	1	Ctrl Samp		
6	Vial 6	CTRL 1 (0.04)	SIMALC3	1	Ctrl Samp		
7	Vial 7	CTRL 2 (0.10)	SIMALC3	1	Ctrl Samp		
8	Vial 8	CTRL 3 (0.20)	SIMALC3	1	Ctrl Samp		
9	Vial 9	NEG CTRL	SIMALC3	1	Ctrl Samp		
10	Vial 10	15004 #1	SIMALC3	1	Sample		
11	Vial 11	15004 #2	SIMALC3	1	Sample		
12	Vial 12	15004 #3	SIMALC3	1	Sample		
13	Vial 13	15004 #4	SIMALC3	1	Sample		
14	Vial 14	15004 #5	SIMALC3	1	Sample		
15	Vial 15	POS CTRL (0.10)	SIMALC3	1	Ctrl Samp		
16	Vial 16	NEG CTRL	SIMALC3	1	Ctrl Samp		
17	Vial 17	15005 #1	SIMALC3	1	Sample		
18	Vial 18	15005 #2	SIMALC3	1	Sample		
19	Vial 19	15005 #3	SIMALC3	1	Sample		
20	Vial 20	15005 #4	SIMALC3	1	Sample		
21	Vial 21	15005 #5	SIMALC3	1	Sample		
22	Vial 22	POS CTRL (0.10)	SIMALC3	1	Ctrl Samp		
23	Vial 23	NEG CTRL	SIMALC3	1	Ctrl Samp		
24	Vial 24	15006 #1	SIMALC3	1	Sample		

~~1500~~
 Jan/27/15
~~1500~~
 15006 Jan/27/15
 Jn

Line	Location	SampleName	Method	Inj	SampleType	InjVolume	DataFile
25	Vial 25	15006 #2	SIMALC3	1	Sample		
26	Vial 26	15006 #3	SIMALC3	1	Sample		
27	Vial 27	15006 #4	SIMALC3	1	Sample		
28	Vial 28	15006 #5	SIMALC3	1	Sample		
29	Vial 29	POS CTRL (0.10)	SIMALC3	1	Ctrl Samp		
30	Vial 30	NEG CTRL	SIMALC3	1	Ctrl Samp		
31	Vial 31	15007 #1	SIMALC3	1	Sample		
32	Vial 32	15007 #2	SIMALC3	1	Sample		
33	Vial 33	15007 #3	SIMALC3	1	Sample		
34	Vial 34	15007 #4	SIMALC3	1	Sample		
35	Vial 35	15007 #5	SIMALC3	1	Sample		
36	Vial 36	POS CTRL (0.10)	SIMALC3	1	Ctrl Samp		
37	Vial 37	NEG CTRL	SIMALC3	1	Ctrl Samp		
38	Vial 38	15008 #1	SIMALC3	1	Sample		
39	Vial 39	15008 #2	SIMALC3	1	Sample		
40	Vial 40	15008 #3	SIMALC3	1	Sample		
41	Vial 41	15008 #4	SIMALC3	1	Sample		
42	Vial 42	15008 #5	SIMALC3	1	Sample		
43	Vial 43	POS CTRL (0.10)	SIMALC3	1	Ctrl Samp		
44	Vial 44	NEG CTRL	SIMALC3	1	Ctrl Samp		

Calibration Part:

Line	Location	SampleName	Method	CalLev	Update	RF	Update	RT	Interval
2	Vial 2	CAL 1 (0.079)	SIMALC3	1	Replace		Replace		
3	Vial 3	CAL 2 (0.158)	SIMALC3	2	Replace		Replace		
4	Vial 4	CAL 3 (0.316)	SIMALC3	3	Replace		Replace		

Sequence Table (Back Injector):

No entries - empty table!

15006
fn 1/27/15

fn

Washington State Patrol Toxicology Laboratory
2203 Airport Way S Seattle, WA 98134

Inj. Date: 1/14/2015 6:01:14 PM

Sample Name: 15006 #1

Instrument: HSGC#3


Operator: Naziha Nuwayhid, PhD

Column: DB-ALC2


Location: Vial 24

Method: C:\HPCHEM\2\METHODS\SIMALC3.M


Sample Info:



#	Compound	Peak Area	RT (min)
1	Ethanol	910	1.018
2	n-Propanol	1573	1.744



Ethanol 0.127 g/100mL



n-Propanol 0.012 g/100mL

Handwritten signature

Handwritten mark

Inj. Date: 1/14/2015 6:04:28 PM

Sample Name: 15006 #2

Instrument: HSGC#3


Operator: Naziha Nuwayhid, PhD

Column: DB-ALC2


Location: Vial 25

Method: C:\HPCHEM\2\METHODS\SIMALC3.M


Sample Info:



#	Compound	Peak Area	RT (min)
1	Ethanol	889	1.018
2	n-Propanol	1542	1.744



Ethanol 0.126 g/100mL



n-Propanol 0.012 g/100mL

Handwritten signature

Handwritten mark

Washington State Patrol Toxicology Laboratory
2203 Airport Way S Seattle, WA 98134

Inj. Date: 1/14/2015 6:07:41 PM

Sample Name: 15006 #3

Instrument: HSGC#3


Operator: Naziha Nuwayhid, PhD

Column: DB-ALC2


Location: Vial 26

Method: C:\HPCHEM\2\METHODS\SIMALC3.M


Sample Info:



#	Compound	Peak Area	RT (min)
1	Ethanol	913	1.018
2	n-Propanol	1577	1.744



Ethanol 0.127 g/100mL



n-Propanol 0.012 g/100mL

dr

M

Inj. Date: 1/14/2015 6:10:54 PM

Sample Name: 15006 #4

Instrument: HSGC#3


Operator: Naziha Nuwayhid, PhD

Column: DB-ALC2


Location: Vial 27

Method: C:\HPCHEM\2\METHODS\SIMALC3.M


Sample Info:



#	Compound	Peak Area	RT (min)
1	Ethanol	895	1.019
2	n-Propanol	1551	1.744



Ethanol 0.127 g/100mL



n-Propanol 0.012 g/100mL

Handwritten signature

Handwritten mark

Inj. Date: 1/14/2015 6:14:08 PM

Sample Name: 15006 #5

Instrument: HSGC#3


Operator: Naziha Nuwayhid, PhD

Column: DB-ALC2


Location: Vial 28

Method: C:\HPCHEM\2\METHODS\SIMALC3.M


Sample Info:



#	Compound	Peak Area	RT (min)
1	Ethanol	899	1.019
2	n-Propanol	1550	1.745



Ethanol 0.127 g/100mL



n-Propanol 0.012 g/100mL


Handwritten signature

Handwritten initials


Washington State Patrol Toxicology Laboratory
 2203 Airport Way S Seattle, WA 98134

Inj. Date: 1/14/2015 6:17:21 PM Sample Name: POS CTRL (0.10)
 Instrument: HSGC#3 Operator: Naziha Nuwayhid, PhD
 Column: DB-ALC2 Location: Vial 29
 Method: C:\HPCHEM\2\METHODS\SIMALC3.M


Sample Info: POS CTRL: 0.10 g/100mL
~~15004, 15005, 15006, 15007, 15008~~



#	Compound	Peak Area	RT (min)
1	Ethanol	722	1.018
2	n-Propanol	1551	1.743



Ethanol 0.102 g/100mL



n-Propanol 0.012 g/100mL

fn

MW

Washington State Patrol Toxicology Laboratory
 2203 Airport Way S Seattle, WA 98134

Inj. Date: 1/14/2015 6:20:34 PM

Sample Name: NEG CTRL

Instrument: HSGC#3

Operator: Naziha Nuwayhid, PhD

Column: DB-ALC2


Location: Vial 30

Method: C:\HPCHEM\2\METHODS\SIMALC3.M


Sample Info: ~~15004, 15005, 15006, 15007, 15008~~

1.15.15 MW


1.15.15 NW



#	Compound	Peak Area	RT (min)
1	Ethanol	0	0.000
2	n-Propanol	1507	1.742



Ethanol 0.000 g/100mL



n-Propanol 0.012 g/100mL

Sh

MW