



QUALITY ASSURANCE PROCEDURE SOLUTION TEST REPORT

BATCH REPORT: 15038

CUSTOMER INFORMATION

Washington State Patrol – Breath Test Program
811 East Roanoke SEATTLE, WA 98102

TESTING PROCEDURE USED: TLD Technical Manual, Chapter 4.0 Certification of Simulator Solutions; Headspace-Gas Chromatography.

TESTING ITEM INFORMATION

TARGET VAPOR CONCENTRATION: 0.20 g/210L
DATE PREPARED: 09/10/2015
BATCH UNITS: g/100mL

IDENTITY: QAP Solution
PREPARED BY: Naziha Nuwayhid

	NN	RF	AC
1	0.253	0.256	0.246
2	0.255	0.257	0.247
3	0.258	0.254	0.252
4	0.256	0.256	0.249
5	0.258	0.256	0.250
C	0.104	0.103	0.103

ETHANOL CONTROL INFORMATION

LOT NUMBER: FN08051301 EXPIRATION: 10/2018 CONCENTRATION: 0.10 g/100mL

RESULTS OF TESTING

AVERAGE SOLUTION CONCENTRATION: 0.2535 g/100mL PRECISION CV (%): 1.54
STANDARD DEVIATION: 0.00391 NUMBER OF TESTS: 15

EQUIVALENT VAPOR CONCENTRATION: **0.2061 g/210L**
EXPANDED UNCERTAINTY: ± 0.0048 (k=2, 95.45% confidence interval)




WASHINGTON STATE PATROL – TOXICOLOGY LABORATORY DIVISION



Lisa Noble Forensic Scientist Supervisor

9/25/15

DATE REPORT ISSUED

ANALYST	NAME	THIS TESTING WAS PERFORMED BY:		DATE TESTED
		SIGNATURE		
NN	Naziha Nuwayhid			09/10/2015
RF	Rebecca Flaherty			09/11/2015
AC	Amanda Chandler			09/17/2015

Washington State Patrol - Toxicology Laboratory Division
 QAP Test Report Calculation Record

QAP Solution Batch #: 15038

Date Prepared: 9/10/2015

Analyst:	NN	RF	AC
Date Tested:	9/10/2015	9/11/2015	9/17/2015
Instrument:	HSGC #3	HSGC #3	HSGC #3
1	0.253	0.256	0.246
2	0.255	0.257	0.247
3	0.258	0.254	0.252
4	0.256	0.256	0.249
5	0.258	0.256	0.250
C	0.104	0.103	0.103

CV^2_{COA}	$CV^2_{QAP\ Solution}$	$CV^2_{Control}$	$CV^2_{Part\ Coef}$
0.0000084100	0.0000158337	0.0000104058	0.0001016326

Ethanol Control Lot #: FN08051301
 Control Uncertainty (%): 0.29

Average Solution Concentration: 0.2535 g/100mL
 Standard Deviation: 0.00391 g/100mL
 Precision CV (%): 1.54
 Equivalent Vapor Concentration: 0.2061 g/210L
 Combined Standard Uncertainty (\pm): 0.0024 g/210L
 Expanded Uncertainty (\pm): 0.0048 coverage factor (k) =2 (95.45% level of confidence)

Calculations performed by: Lisa Noble [Signature] 9/21/15
 Name Signature Date

Calculations verified by: Amanda M. Black [Signature] 9-24-15 Method: Hand calculation
 Name Signature Date

Tech. review performed by: Lisa Noble [Signature] 9/21/15
 Name Signature Date

[Signature]

SIMULATOR SOLUTION DATA ENTRY REVIEW

Reviewer/s: Amanda M. Black Date: 9-24-15

Location: WSP-TLD(RSB) Seattle, WA Solution Batch Number: 15038

	YES	NO	N/A
Analysis dates do not precede preparation date:	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Declarations signed and properly dated:	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Data entry corresponds to all chromatograms:	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
All signatures present on Test Report:	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Average solution concentration correct:	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Standard deviation correct:	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
CV (%) correct:	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Equivalent vapor concentration correct:	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
All chromatograms and sequences included in file:	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Ethanol control information present: (lot # present & used within expiration)	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Complies with accuracy and precision requirements established by the State Toxicologist:	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

Comments:

Reviewer Signature:  Date: 9-24-15



SOLUTION CERTIFICATE REVIEW

Please check that the data on your chromatograms is the data entered into the Test Report, that the date to the right of your name is the date that you tested the solution, and then sign the Test Report.

Please initial and date below to affirm that you have:

- 1) Checked your data
- 2) Checked the date to the right of your name on the Test Report
- 3) Signed the Test Report

	Initials	Date
Amanda Chandler	AC	9/21/15
Andrew Gingras		
Asa Louis		
Brittany Thomas		
Christie Mitchell-Mata		
Christopher Johnston		
David Nguyen		
Dawn Sklerov		
Elizabeth Wehner		
Justin Knoy		
Katie Harris		
Lyndsey Lowe		
Naziha Nuwayhid	NW	9.22.15
Rebecca Flaherty	RF	9/21/15

Batch # 15038 9/21/15

JAY INSLEE
Governor



JOHN R. BATISTE
Chief

STATE OF WASHINGTON
WASHINGTON STATE PATROL
WASHINGTON STATE TOXICOLOGY LABORATORY

2203 Airport Way South, Suite 360 • Seattle, Washington 98134-2927 • (206) 262-6100 • FAX (206) 262-6145

**0.20 g/210 L QUALITY ASSURANCE PROCEDURE SOLUTION
CERTIFICATION FOR LOT 15038**


I, Rebecca Flaherty, do certify under penalty of perjury that:

I am employed by the Washington State Toxicology Laboratory, and my responsibilities include the preparation and certification of alcohol solutions for use with evidential breath test instruments.

I possess the following qualifications: BS degrees in Biochemistry and Psychobiology and MS degree in Forensic Science.

The quality assurance procedure (QAP) solution, Lot Number 15038, was prepared in the Washington State Toxicology Laboratory on 9/10/2015. I tested this solution and it was found to conform to those standards established by the State Toxicologist for the certification of simulator solution. It shall not be used to perform a quality assurance procedure after 9/10/2016.

Seattle, WA

 9/21/15

Rebecca Flaherty

Date

Forensic Scientist

JAY INSLEE
Governor



JOHN R. BATISTE
Chief

STATE OF WASHINGTON
WASHINGTON STATE PATROL
WASHINGTON STATE TOXICOLOGY LABORATORY

2203 Airport Way South, Suite 360 • Seattle, Washington 98134-2927 • (206) 262-6100 • FAX (206) 262-6145

**0.20 g/210 L QUALITY ASSURANCE PROCEDURE SOLUTION
CERTIFICATION FOR LOT 15038**


I, Amanda Chandler, do certify under penalty of perjury that:

I am employed by the Washington State Toxicology Laboratory, and my responsibilities include the preparation and certification of alcohol solutions for use with evidential breath test instruments.

I possess the following qualifications: MS degree in Forensic Toxicology.

The quality assurance procedure (QAP) solution, Lot Number 15038, was prepared in the Washington State Toxicology Laboratory on 9/10/2015. I tested this solution and it was found to conform to those standards established by the State Toxicologist for the certification of simulator solution. It shall not be used to perform a quality assurance procedure after 9/10/2016.

Seattle, WA

 9/21/15

Amanda Chandler

Date

Forensic Scientist

FILE A COPY IN THE BATCH FILE FOR EACH SOLUTION LISTED ON THE WORKSHEET

Preparation Date: 9.10.15 Expiration Date: 9.10.16 Initials of Preparer: NN

Lot # of 200-proof Ethanol used in preparation: 2DC0208

Date the 200-proof Ethanol bottle was opened: 6.16.15 & 9.10.15

After opening, each bottle of 200-proof Ethanol is approved for use for 6 months unless an extension is approved by the State Toxicologist. This timeframe applies to the 200-proof Ethanol only, not to simulator solutions which have a 1 year expiration.

Environmental conditions verified as acceptable:

Simulator Solution	Volume of Ethanol (mL)	Volume of Deionized Water (L)		Batch Number
QAP 0.04	11.2	18	<input checked="" type="checkbox"/>	<u>15034</u>
QAP 0.08	22.4	18	<input checked="" type="checkbox"/>	<u>15035</u>
QAP 0.10	28.1	18	<input checked="" type="checkbox"/>	<u>15036</u>
QAP 0.15	42.1	18	<input checked="" type="checkbox"/>	<u>15037</u>
QAP 0.20	56.1	18	<input checked="" type="checkbox"/>	<u>15038</u>
ESS	66.5	52	<input checked="" type="checkbox"/>	
			<input checked="" type="checkbox"/>	
		Stir bar is rotating	<input checked="" type="checkbox"/>	
		Stirred for minimum 30 minutes; 2 hours for ESS	<input checked="" type="checkbox"/>	
		Spigot purged	<input checked="" type="checkbox"/>	
		Aliquot taken	<input checked="" type="checkbox"/>	
		Batch labeled, packaged and sealed	<input checked="" type="checkbox"/>	<u>9.10.15</u> Date

NO ESS prepared 9.24.15 (NN)

If different ethanol lot numbers are used in the preparation of solutions, record them and the corresponding solution batch numbers in the comments section.

Comments:

QAP 0.04, 0.08, & 0.10 were prepared from bottle opened on 6.16.15
QAP 0.15, & 0.20 were prepared from ethanol bottle opened 9.10.15

Michelle Myler
 Analyst Signature

9.10.15
 Date

fr

Sequence Parameters:

Operator: Naziha Nuwayhid, PhD
 Data File Naming: Prefix/Counter
 Signal 1 Prefix: SIG1
 Counter: 0001
 Signal 2 Prefix: SIG2
 Counter: 0001
 Data Directory: C:\HPCHEM\2\DATA\
 Data Subdirectory: 150910NN
 Part of Methods to run: According to Runtime Checklist
 Barcode Reader: not used
 Shutdown Cmd/Macro: none

Sequence Comment:

Ethanol Calibrator 1, E0615-01 - Exp. 12/2/2015
 Ethanol Calibrator 2, E0615-02 - Exp. 12/2/2015
 Ethanol Calibrator 3, E0615-03 - Exp. 12/2/2015
 CTRL1 (0.04g/100mL), Lot # FN05011301 - Exp. 05/2018
 CTRL2 (0.10g/100mL), Lot # FN08051301 - Exp. 10/2018
 CTRL3 (0.20g/100mL), Lot # FN03211401 - Exp. 06/2019
 Internal Standard Lot#P0715 - Exp. 10/27/15

Calibration vials 1-9 filed with 15034.

Sequence Table (Front Injector):

Method and Injection Info Part:

Line	Location	SampleName	Method	Inj	SampleType	InjVolume	DataFile
1	Vial 1	BLANK	SIMALC3	1	Sample		
2	Vial 2	0.079 CAL 1	SIMALC3	1	Calib		
3	Vial 3	0.158 CAL 2	SIMALC3	1	Calib		
4	Vial 4	0.316 CAL 3	SIMALC3	1	Calib		
5	Vial 5	NEG CTRL	SIMALC3	1	Ctrl Samp		
6	Vial 6	0.04 CTRL	SIMALC3	1	Ctrl Samp		
7	Vial 7	0.10 CTRL	SIMALC3	1	Ctrl Samp		
8	Vial 8	0.20 CTRL	SIMALC3	1	Ctrl Samp		
9	Vial 9	NEG CTRL	SIMALC3	1	Ctrl Samp		
10	Vial 10	15034 #1	SIMALC3	1	Sample		
11	Vial 11	15034 #2	SIMALC3	1	Sample		
12	Vial 12	15034 #3	SIMALC3	1	Sample		
13	Vial 13	15034 #4	SIMALC3	1	Sample		
14	Vial 14	15034 #5	SIMALC3	1	Sample		
15	Vial 15	0.10 CTRL	SIMALC3	1	Ctrl Samp		
16	Vial 16	NEG CTRL	SIMALC3	1	Ctrl Samp		
17	Vial 17	15035 #1	SIMALC3	1	Sample		
18	Vial 18	15035 #2	SIMALC3	1	Sample		
19	Vial 19	15035 #3	SIMALC3	1	Sample		
20	Vial 20	15035 #4	SIMALC3	1	Sample		
21	Vial 21	15035 #5	SIMALC3	1	Sample		
22	Vial 22	0.10 CTRL	SIMALC3	1	Ctrl Samp		
23	Vial 23	NEG CTRL	SIMALC3	1	Ctrl Samp		
24	Vial 24	15036 #1	SIMALC3	1	Sample		
25	Vial 25	15036 #2	SIMALC3	1	Sample		
26	Vial 26	15036 #3	SIMALC3	1	Sample		

15038
Inj file

nn

Line	Location	SampleName	Method	Inj	SampleType	InjVolume	DataFile
27	Vial 27	15036 #4	SIMALC3	1	Sample		
28	Vial 28	15036 #5	SIMALC3	1	Sample		
29	Vial 29	0.10 CTRL	SIMALC3	1	Ctrl Samp		
30	Vial 30	NEG CTRL	SIMALC3	1	Ctrl Samp		
31	Vial 31	15037 #1	SIMALC3	1	Sample		
32	Vial 32	15037 #2	SIMALC3	1	Sample		
33	Vial 33	15037 #3	SIMALC3	1	Sample		
34	Vial 34	15037 #4	SIMALC3	1	Sample		
35	Vial 35	15037 #5	SIMALC3	1	Sample		
36	Vial 36	0.10 CTRL	SIMALC3	1	Ctrl Samp		
37	Vial 37	NEG CTRL	SIMALC3	1	Ctrl Samp		
38	Vial 38	15038 #1	SIMALC3	1	Sample		
39	Vial 39	15038 #2	SIMALC3	1	Sample		
40	Vial 40	15038 #3	SIMALC3	1	Sample		
41	Vial 41	15038 #4	SIMALC3	1	Sample		
42	Vial 42	15038 #5	SIMALC3	1	Sample		
43	Vial 43	0.10 CTRL	SIMALC3	1	Ctrl Samp		
44	Vial 44	NEG CTRL	SIMALC3	1	Ctrl Samp		

Calibration Part:

Line	Location	SampleName	Method	CalLev	Update	RF	Update	RT	Interval
2	Vial 2	0.079 CAL 1	SIMALC3	1	Replace		Replace		
3	Vial 3	0.158 CAL 2	SIMALC3	2	Replace		Replace		
4	Vial 4	0.316 CAL 3	SIMALC3	3	Replace		Replace		

Sequence Table (Back Injector):

No entries - empty table!

15038
for 9/21/15

Washington State Patrol Toxicology Laboratory
 2203 Airport Way S Seattle, WA 98134

Inj. Date: 9/10/2015 1:50:05 PM

Sample Name: 15038 #1

Instrument: HSGC#3


Operator: Naziha Nuwayhid, PhD

Column: DB-ALC2


Location: Vial 38

Method: C:\HPCHEM\2\METHODS\SIMALC3.M


Sample Info:



#	Compound	Peak Area	RT (min)
1	Ethanol	1767	1.022
2	n-Propanol	1593	1.748



Ethanol 0.253 g/100mL



n-Propanol 0.012 g/100mL

Handwritten initials

Handwritten signature

Washington State Patrol Toxicology Laboratory
 2203 Airport Way S Seattle, WA 98134

Inj. Date: 9/10/2015 1:53:18 PM

Sample Name: 15038 #2

Instrument: HSGC#3


Operator: Naziha Nuwayhid, PhD

Column: DB-ALC2


Location: Vial 39

Method: C:\HPCHEM\2\METHODS\SIMALC3.M


Sample Info:



#	Compound	Peak Area	RT (min)
1	Ethanol	1806	1.023
2	n-Propanol	1614	1.749



Ethanol 0.255 g/100mL



n-Propanol 0.012 g/100mL

Handwritten mark

Handwritten mark

Washington State Patrol Toxicology Laboratory
 2203 Airport Way S Seattle, WA 98134

Inj. Date: 9/10/2015 1:56:32 PM

Sample Name: 15038 #3

Instrument: HSGC#3


Operator: Naziha Nuwayhid, PhD

Column: DB-ALC2


Location: Vial 40

Method: C:\HPCHEM\2\METHODS\SIMALC3.M


Sample Info:



#	Compound	Peak Area	RT (min)
1	Ethanol	1896	1.023
2	n-Propanol	1680	1.748



Ethanol 0.258 g/100mL



n-Propanol 0.012 g/100mL

Handwritten signature

Handwritten signature

Washington State Patrol Toxicology Laboratory
 2203 Airport Way S Seattle, WA 98134

Inj. Date: 9/10/2015 1:59:45 PM

Sample Name: 15038 #4

Instrument: HSGC#3


Operator: Naziha Nuwayhid, PhD

Column: DB-ALC2


Location: Vial 41

Method: C:\HPCHEM\2\METHODS\SIMALC3.M


Sample Info:



#	Compound	Peak Area	RT (min)
1	Ethanol	1868	1.022
2	n-Propanol	1663	1.748



Ethanol 0.256 g/100mL



n-Propanol 0.012 g/100mL

Handwritten mark

Handwritten signature

Washington State Patrol Toxicology Laboratory
 2203 Airport Way S Seattle, WA 98134

Inj. Date: 9/10/2015 2:02:59 PM

Sample Name: 15038 #5

Instrument: HSGC#3


Operator: Naziha Nuwayhid, PhD

Column: DB-ALC2


Location: Vial 42

Method: C:\HPCHEM\2\METHODS\SIMALC3.M


Sample Info:



#	Compound	Peak Area	RT (min)
1	Ethanol	1967	1.023
2	n-Propanol	1742	1.749



Ethanol 0.258 g/100mL



n-Propanol 0.012 g/100mL

Handwritten signature

Handwritten initials

Washington State Patrol Toxicology Laboratory
 2203 Airport Way S Seattle, WA 98134

Inj. Date: 9/10/2015 2:06:12 PM

Sample Name: 0.10 CTRL

Instrument: HSGC#3


Operator: Naziha Nuwayhid, PhD

Column: DB-ALC2


Location: Vial 43

Method: C:\HPCHEM\2\METHODS\SIMALC3.M


Sample Info: 15038



#	Compound	Peak Area	RT (min)
1	Ethanol	798	1.022
2	n-Propanol	1747	1.749



Ethanol 0.104 g/100mL



n-Propanol 0.012 g/100mL

Handwritten initials

Handwritten initials

Washington State Patrol Toxicology Laboratory
2203 Airport Way S Seattle, WA 98134

Inj. Date: 9/10/2015 2:09:28 PM

Sample Name: NEG CTRL

Instrument: HSGC#3


Operator: Naziha Nuwayhid, PhD

Column: DB-ALC2


Location: Vial 44

Method: C:\HPCHEM\2\METHODS\SIMALC3.M


Sample Info: 15038



#	Compound	Peak Area	RT (min)
1	Ethanol	0	0.000
2	n-Propanol	1706	1.748



Ethanol 0.000 g/100mL



n-Propanol 0.012 g/100mL

Handwritten signature

Handwritten initials

Sequence Parameters:

Operator: Rebecca Flaherty
 Data File Naming: Prefix/Counter
 Signal 1 Prefix: SIG1
 Counter: 0001
 Signal 2 Prefix: SIG2
 Counter: 0001
 Data Directory: C:\HPCHEM\2\DATA\
 Data Subdirectory: 150911RF
 Part of Methods to run: According to Runtime Checklist
 Barcode Reader: not used
 Shutdown Cmd/Macro: none

Sequence Comment:

Ethanol Calibrator 1, E0615-01 - Exp. 12/2/2015
 Ethanol Calibrator 2, E0615-02 - Exp. 12/2/2015
 Ethanol Calibrator 3, E0615-03 - Exp. 12/2/2015
 CTRL1 (0.04g/100mL), Lot # FN05011301 - Exp. 05/2018
 CTRL2 (0.10g/100mL), Lot # FN08051301 - Exp. 10/2018
 CTRL3 (0.20g/100mL), Lot # FN03211401 - Exp. 06/2019
 Internal Standard Lot#P0715 - Exp. 10/27/15

Calibration vials 1-9 filed with 15034.

*Run inadvertently stopped
 after vial 34 due to wrong
 vial count downloaded to
 auto sampler. changed to
 changed vial #s in autosampler
 control & sequence automatically
 resumed RF 9.11.15*

Sequence Table (Front Injector):

Method and Injection Info Part:

Line	Location	SampleName	Method	Inj	SampleType	InjVolume	DataFile
1	Vial 1	BLANK	SIMALC3	1	Sample		
2	Vial 2	0.079 CAL 1	SIMALC3	1	Calib		
3	Vial 3	0.158 CAL 2	SIMALC3	1	Calib		
4	Vial 4	0.316 CAL 3	SIMALC3	1	Calib		
5	Vial 5	NEG CTRL	SIMALC3	1	Ctrl Samp		
6	Vial 6	0.04 CTRL	SIMALC3	1	Ctrl Samp		
7	Vial 7	0.10 CTRL	SIMALC3	1	Ctrl Samp		
8	Vial 8	0.20 CTRL	SIMALC3	1	Ctrl Samp		
9	Vial 9	NEG CTRL	SIMALC3	1	Ctrl Samp		
10	Vial 10	15034 #1	SIMALC3	1	Sample		
11	Vial 11	15034 #2	SIMALC3	1	Sample		
12	Vial 12	15034 #3	SIMALC3	1	Sample		
13	Vial 13	15034 #4	SIMALC3	1	Sample		
14	Vial 14	15034 #5	SIMALC3	1	Sample		
15	Vial 15	0.10 CTRL	SIMALC3	1	Ctrl Samp		
16	Vial 16	NEG CTRL	SIMALC3	1	Ctrl Samp		
17	Vial 17	15035 #1	SIMALC3	1	Sample		
18	Vial 18	15035 #2	SIMALC3	1	Sample		
19	Vial 19	15035 #3	SIMALC3	1	Sample		
20	Vial 20	15035 #4	SIMALC3	1	Sample		
21	Vial 21	15035 #5	SIMALC3	1	Sample		
22	Vial 22	0.10 CTRL	SIMALC3	1	Ctrl Samp		
23	Vial 23	NEG CTRL	SIMALC3	1	Ctrl Samp		
24	Vial 24	15036 #1	SIMALC3	1	Sample		
25	Vial 25	15036 #2	SIMALC3	1	Sample		
26	Vial 26	15036 #3	SIMALC3	1	Sample		

15038

9/11/15

RF

Line	Location	SampleName	Method	Inj	SampleType	InjVolume	DataFile
27	Vial 27	15036 #4	SIMALC3	1	Sample		
28	Vial 28	15036 #5	SIMALC3	1	Sample		
29	Vial 29	0.10 CTRL	SIMALC3	1	Ctrl Samp		
30	Vial 30	NEG CTRL	SIMALC3	1	Ctrl Samp		
31	Vial 31	15037 #1	SIMALC3	1	Sample		
32	Vial 32	15037 #2	SIMALC3	1	Sample		
33	Vial 33	15037 #3	SIMALC3	1	Sample		
34	Vial 34	15037 #4	SIMALC3	1	Sample		
35	Vial 35	15037 #5	SIMALC3	1	Sample		
36	Vial 36	0.10 CTRL	SIMALC3	1	Ctrl Samp		
37	Vial 37	NEG CTRL	SIMALC3	1	Ctrl Samp		
38	Vial 38	15038 #1	SIMALC3	1	Sample		
39	Vial 39	15038 #2	SIMALC3	1	Sample		
40	Vial 40	15038 #3	SIMALC3	1	Sample		
41	Vial 41	15038 #4	SIMALC3	1	Sample		
42	Vial 42	15038 #5	SIMALC3	1	Sample		
43	Vial 43	0.10 CTRL	SIMALC3	1	Ctrl Samp		
44	Vial 44	NEG CTRL	SIMALC3	1	Ctrl Samp		

Calibration Part:

Line	Location	SampleName	Method	CalLev	Update	RF	Update	RT	Interval
2	Vial 2	0.079 CAL 1	SIMALC3	1	Replace		Replace		
3	Vial 3	0.158 CAL 2	SIMALC3	2	Replace		Replace		
4	Vial 4	0.316 CAL 3	SIMALC3	3	Replace		Replace		

Sequence Table (Back Injector):

No entries - empty table!

15038
Inj/2/15

RF

Washington State Patrol Toxicology Laboratory
2203 Airport Way S Seattle, WA 98134

Inj. Date: 9/11/2015 2:56:25 PM

Sample Name: 15038 #1

Instrument: HSGC#3


Operator: Rebecca Flaherty

Column: DB-ALC2


Location: Vial 38

Method: C:\HPCHEM\2\METHODS\SIMALC3.M


Sample Info:



#	Compound	Peak Area	RT (min)
1	Ethanol	1799	1.021
2	n-Propanol	1606	1.748



Ethanol 0.256 g/100mL



n-Propanol 0.012 g/100mL

RF

RF

Washington State Patrol Toxicology Laboratory
 2203 Airport Way S Seattle, WA 98134

Inj. Date: 9/11/2015 2:59:38 PM

Sample Name: 15038 #2

Instrument: HSGC#3


Operator: Rebecca Flaherty

Column: DB-ALC2


Location: Vial 39

Method: C:\HPCHEM\2\METHODS\SIMALC3.M


Sample Info:



#	Compound	Peak Area	RT (min)
1	Ethanol	1833	1.022
2	n-Propanol	1630	1.749



Ethanol 0.257 g/100mL



n-Propanol 0.012 g/100mL

RF

RF

Washington State Patrol Toxicology Laboratory
 2203 Airport Way S Seattle, WA 98134

Inj. Date: 9/11/2015 3:02:51 PM

Sample Name: 15038 #3

Instrument: HSGC#3


Operator: Rebecca Flaherty

Column: DB-ALC2


Location: Vial 40

Method: C:\HPCHEM\2\METHODS\SIMALC3.M


Sample Info:



#	Compound	Peak Area	RT (min)
1	Ethanol	1766	1.021
2	n-Propanol	1585	1.748



Ethanol 0.254 g/100mL



n-Propanol 0.012 g/100mL

RF

RF

Washington State Patrol Toxicology Laboratory
 2203 Airport Way S Seattle, WA 98134

Inj. Date: 9/11/2015 3:06:05 PM

Sample Name: 15038 #4

Instrument: HSGC#3


Operator: Rebecca Flaherty

Column: DB-ALC2


Location: Vial 41

Method: C:\HPCHEM\2\METHODS\SIMALC3.M


Sample Info:



#	Compound	Peak Area	RT (min)
1	Ethanol	1830	1.022
2	n-Propanol	1630	1.749



Ethanol 0.256 g/100mL



n-Propanol 0.012 g/100mL

RF

RF

Washington State Patrol Toxicology Laboratory
 2203 Airport Way S Seattle, WA 98134

Inj. Date: 9/11/2015 3:09:19 PM

Sample Name: 15038 #5

Instrument: HSGC#3


Operator: Rebecca Flaherty

Column: DB-ALC2


Location: Vial 42

Method: C:\HPCHEM\2\METHODS\SIMALC3.M


Sample Info:



#	Compound	Peak Area	RT (min)
1	Ethanol	1806	1.022
2	n-Propanol	1610	1.749



Ethanol 0.256 g/100mL



n-Propanol 0.012 g/100mL

RF

RF

Washington State Patrol Toxicology Laboratory
 2203 Airport Way S Seattle, WA 98134

Inj. Date: 9/11/2015 3:12:33 PM

Sample Name: 0.10 CTRL

Instrument: HSGC#3


Operator: Rebecca Flaherty

Column: DB-ALC2


Location: Vial 43

Method: C:\HPCHEM\2\METHODS\SIMALC3.M


Sample Info: 15038



#	Compound	Peak Area	RT (min)
1	Ethanol	732	1.022
2	n-Propanol	1626	1.750



Ethanol 0.103 g/100mL



n-Propanol 0.012 g/100mL

RF

RF

Inj. Date: 9/11/2015 3:15:45 PM

Sample Name: NEG CTRL

Instrument: HSGC#3


Operator: Rebecca Flaherty

Column: DB-ALC2


Location: Vial 44

Method: C:\HPCHEM\2\METHODS\SIMALC3.M


Sample Info: 15038



#	Compound	Peak Area	RT (min)
1	Ethanol	0	0.000
2	n-Propanol	1656	1.754



Ethanol 0.000 g/100mL



n-Propanol 0.012 g/100mL

RF

RF

Sequence Parameters:

Operator: Amanda Chandler
 Data File Naming: Prefix/Counter
 Signal 1 Prefix: SIG1
 Counter: 0001
 Signal 2 Prefix: SIG2
 Counter: 0001
 Data Directory: C:\HPCHEM\2\DATA\
 Data Subdirectory: 150917AC
 Part of Methods to run: According to Runtime Checklist
 Barcode Reader: not used
 Shutdown Cmd/Macro: none

Sequence Comment:

Ethanol Calibrator 1, E0615-01 - Exp. 12/2/2015
 Ethanol Calibrator 2, E0615-02 - Exp. 12/2/2015
 Ethanol Calibrator 3, E0615-03 - Exp. 12/2/2015
 CTRL1 (0.04g/100mL), Lot # FN05011301 - Exp. 05/2018
 CTRL2 (0.10g/100mL), Lot # FN08051301 - Exp. 10/2018
 CTRL3 (0.20g/100mL), Lot # FN03211401 - Exp. 06/2019
 Internal Standard Lot#P0715 - Exp. 10/27/15

Calibration vials 1-9 filed with 15034.

Sequence Table (Front Injector):

Method and Injection Info Part:

Line	Location	SampleName	Method	Inj	SampleType	InjVolume	DataFile
1	Vial 1	BLANK	SIMALC3	1	Sample		
2	Vial 2	0.079 CAL 1	SIMALC3	1	Calib		
3	Vial 3	0.158 CAL 2	SIMALC3	1	Calib		
4	Vial 4	0.316 CAL 3	SIMALC3	1	Calib		
5	Vial 5	NEG CTRL	SIMALC3	1	Ctrl Samp		
6	Vial 6	0.04 CTRL	SIMALC3	1	Ctrl Samp		
7	Vial 7	0.10 CTRL	SIMALC3	1	Ctrl Samp		
8	Vial 8	0.20 CTRL	SIMALC3	1	Ctrl Samp		
9	Vial 9	NEG CTRL	SIMALC3	1	Ctrl Samp		
10	Vial 10	15034 #1	SIMALC3	1	Sample		
11	Vial 11	15034 #2	SIMALC3	1	Sample		
12	Vial 12	15034 #3	SIMALC3	1	Sample		
13	Vial 13	15034 #4	SIMALC3	1	Sample		
14	Vial 14	15034 #5	SIMALC3	1	Sample		
15	Vial 15	0.10 CTRL	SIMALC3	1	Ctrl Samp		
16	Vial 16	NEG CTRL	SIMALC3	1	Ctrl Samp		
17	Vial 17	15035 #1	SIMALC3	1	Sample		
18	Vial 18	15035 #2	SIMALC3	1	Sample		
19	Vial 19	15035 #3	SIMALC3	1	Sample		
20	Vial 20	15035 #4	SIMALC3	1	Sample		
21	Vial 21	15035 #5	SIMALC3	1	Sample		
22	Vial 22	0.10 CTRL	SIMALC3	1	Ctrl Samp		
23	Vial 23	NEG CTRL	SIMALC3	1	Ctrl Samp		
24	Vial 24	15036 #1	SIMALC3	1	Sample		
25	Vial 25	15036 #2	SIMALC3	1	Sample		
26	Vial 26	15036 #3	SIMALC3	1	Sample		

15038
 12/2/15

AL

Line	Location	SampleName	Method	Inj	SampleType	InjVolume	DataFile
27	Vial 27	15036 #4	SIMALC3	1	Sample		
28	Vial 28	15036 #5	SIMALC3	1	Sample		
29	Vial 29	0.10 CTRL	SIMALC3	1	Ctrl Samp		
30	Vial 30	NEG CTRL	SIMALC3	1	Ctrl Samp		
31	Vial 31	15038 #1	SIMALC3	1	Sample		
32	Vial 32	15038 #2	SIMALC3	1	Sample		
33	Vial 33	15038 #3	SIMALC3	1	Sample		
34	Vial 34	15038 #4	SIMALC3	1	Sample		
35	Vial 35	15038 #5	SIMALC3	1	Sample		
36	Vial 36	0.10 CTRL	SIMALC3	1	Ctrl Samp		
37	Vial 37	NEG CTRL	SIMALC3	1	Ctrl Samp		
38	Vial 38	15039 #1	SIMALC3	1	Sample		
39	Vial 39	15039 #2	SIMALC3	1	Sample		
40	Vial 40	15039 #3	SIMALC3	1	Sample		
41	Vial 41	15039 #4	SIMALC3	1	Sample		
42	Vial 42	15039 #5	SIMALC3	1	Sample		
43	Vial 43	0.10 CTRL	SIMALC3	1	Ctrl Samp		
44	Vial 44	NEG CTRL	SIMALC3	1	Ctrl Samp		

Calibration Part:

Line	Location	SampleName	Method	CalLev	Update	RF	Update	RT	Interval
2	Vial 2	0.079 CAL 1	SIMALC3	1	Replace		Replace		
3	Vial 3	0.158 CAL 2	SIMALC3	2	Replace		Replace		
4	Vial 4	0.316 CAL 3	SIMALC3	3	Replace		Replace		

Sequence Table (Back Injector):

No entries - empty table!

15038
Jan/21/15

AR

Washington State Patrol Toxicology Laboratory
 2203 Airport Way S Seattle, WA 98134

Inj. Date: 9/17/2015 12:54:12 PM

Sample Name: 15038 #1

Instrument: HSGC#3


Operator: Amanda Chandler

Column: DB-ALC2


Location: Vial 31

Method: C:\HPCHEM\2\METHODS\SIMALC3.M


Sample Info:



#	Compound	Peak Area	RT (min)
1	Ethanol	1879	1.023
2	n-Propanol	1746	1.751



Ethanol 0.246 g/100mL



n-Propanol 0.012 g/100mL

Handwritten signature

Handwritten 'R'

Washington State Patrol Toxicology Laboratory
 2203 Airport Way S Seattle, WA 98134

Inj. Date: 9/17/2015 12:57:25 PM

Sample Name: 15038 #2

Instrument: HSGC#3


Operator: Amanda Chandler

Column: DB-ALC2


Location: Vial 32

Method: C:\HPCHEM\2\METHODS\SIMALC3.M


Sample Info:



#	Compound	Peak Area	RT (min)
1	Ethanol	1893	1.023
2	n-Propanol	1748	1.750



Ethanol 0.247 g/100mL



n-Propanol 0.012 g/100mL

fr

AR

Washington State Patrol Toxicology Laboratory
 2203 Airport Way S Seattle, WA 98134

Inj. Date: 9/17/2015 1:00:39 PM

Sample Name: 15038 #3

Instrument: HSGC#3


Operator: Amanda Chandler

Column: DB-ALC2


Location: Vial 33

Method: C:\HPCHEM\2\METHODS\SIMALC3.M


Sample Info:



#	Compound	Peak Area	RT (min)
1	Ethanol	1962	1.023
2	n-Propanol	1776	1.751



Ethanol 0.252 g/100mL



n-Propanol 0.012 g/100mL

Handwritten mark

Handwritten mark

Washington State Patrol Toxicology Laboratory
2203 Airport Way S Seattle, WA 98134

Inj. Date: 9/17/2015 1:03:52 PM

Sample Name: 15038 #4

Instrument: HSGC#3


Operator: Amanda Chandler

Column: DB-ALC2


Location: Vial 34

Method: C:\HPCHEM\2\METHODS\SIMALC3.M


Sample Info:



#	Compound	Peak Area	RT (min)
1	Ethanol	1924	1.023
2	n-Propanol	1762	1.751



Ethanol 0.249 g/100mL



n-Propanol 0.012 g/100mL

Handwritten mark

Handwritten mark

Washington State Patrol Toxicology Laboratory
 2203 Airport Way S Seattle, WA 98134

Inj. Date: 9/17/2015 1:07:06 PM

Sample Name: 15038 #5

Instrument: HSGC#3


Operator: Amanda Chandler

Column: DB-ALC2


Location: Vial 35

Method: C:\HPCHEM\2\METHODS\SIMALC3.M


Sample Info:



#	Compound	Peak Area	RT (min)
1	Ethanol	1944	1.023
2	n-Propanol	1774	1.750



Ethanol 0.250 g/100mL



n-Propanol 0.012 g/100mL

30

AR

Washington State Patrol Toxicology Laboratory
 2203 Airport Way S Seattle, WA 98134

Inj. Date: 9/17/2015 1:10:19 PM

Sample Name: 0.10 CTRL

Instrument: HSGC#3


Operator: Amanda Chandler

Column: DB-ALC2


Location: Vial 36

Method: C:\HPCHEM\2\METHODS\SIMALC3.M


Sample Info: 15038



#	Compound	Peak Area	RT (min)
1	Ethanol	808	1.023
2	n-Propanol	1809	1.752



Ethanol 0.103 g/100mL



n-Propanol 0.012 g/100mL

Handwritten mark

Handwritten mark

Washington State Patrol Toxicology Laboratory
2203 Airport Way S Seattle, WA 98134

Inj. Date: 9/17/2015 1:13:32 PM

Sample Name: NEG CTRL

Instrument: HSGC#3


Operator: Amanda Chandler

Column: DB-ALC2


Location: Vial 37

Method: C:\HPCHEM\2\METHODS\SIMALC3.M


Sample Info: 15038



#	Compound	Peak Area	RT (min)
1	Ethanol	0	0.000
2	n-Propanol	1841	1.752



Ethanol 0.000 g/100mL



n-Propanol 0.012 g/100mL

Handwritten signature

Handwritten signature