



QUALITY ASSURANCE PROCEDURE SOLUTION TEST REPORT

BATCH REPORT: 15052

CUSTOMER INFORMATION

Washington State Patrol – Breath Test Program
811 East Roanoke SEATTLE, WA 98102

TESTING PROCEDURE USED: TLD Technical Manual, Chapter 4.0 Certification of Simulator Solutions; Headspace-Gas Chromatography.

TESTING ITEM INFORMATION

TARGET VAPOR CONCENTRATION: 0.04 g/210L
DATE PREPARED: 11/13/2015
BATCH UNITS: g/100mL

IDENTITY: QAP Solution
PREPARED BY: Justin L. Knoy

	JLK	AG	CM
1	0.051	0.050	0.049
2	0.051	0.050	0.050
3	0.051	0.050	0.050
4	0.051	0.050	0.049
5	0.051	0.051	0.049
C	0.101	0.102	0.101

ETHANOL CONTROL INFORMATION

LOT NUMBER: FN08051301 EXPIRATION: 10/2018 CONCENTRATION: 0.10 g/100mL

RESULTS OF TESTING

AVERAGE SOLUTION CONCENTRATION: 0.0502 g/100mL PRECISION CV (%): 1.54
STANDARD DEVIATION: 0.00077 NUMBER OF TESTS: 15

EQUIVALENT VAPOR CONCENTRATION: **0.0408 g/210L**
EXPANDED UNCERTAINTY: ± 0.0010 (k=2, 95.45% confidence interval)




WASHINGTON STATE PATROL – TOXICOLOGY LABORATORY DIVISION



Lisa Noble Forensic Scientist Supervisor

1/13/16
DATE REPORT ISSUED

THIS TESTING WAS PERFORMED BY:

ANALYST	NAME	SIGNATURE	DATE TESTED
JLK	Justin L. Knoy		11/13/2015
AG	Andrew Gingras		11/13/2015
CM	Christie Mitchell-Mata		11/18/2015

Washington State Patrol - Toxicology Laboratory Division
QAP Test Report Calculation Record

QAP Solution Batch #: 15052

Date Prepared: 11/13/2015

Analyst:	JLK	AG	CM
Date Tested:	11/13/2015	11/13/2015	11/18/2015
Instrument:	HSGC #1	HSGC #1	HSGC #1
1	0.051	0.050	0.049
2	0.051	0.050	0.050
3	0.051	0.050	0.050
4	0.051	0.050	0.049
5	0.051	0.051	0.049
C	0.101	0.102	0.101

CV^2_{COA}	$CV^2_{QAP\ Solution}$	$CV^2_{Control}$	$CV^2_{Part\ Coef}$
0.0000084100	0.0000158728	0.0000108206	0.0001016326

Ethanol Control Lot #: FN08051301
Control Uncertainty (%): 0.29

Average Solution Concentration: 0.0502 g/100mL
Standard Deviation: 0.00077 g/100mL
Precision CV (%): 1.54
Equivalent Vapor Concentration: 0.0408 g/210L
Combined Standard Uncertainty (\pm): 0.0005 g/210L
Expanded Uncertainty (\pm): 0.0010 coverage factor (k) =2 (95.45% level of confidence)

Calculations performed by: Lisa Noble [Signature] 12/10/15
Name Signature Date

Calculations verified by: Amanda M. Black [Signature] 1-12-16
Name Signature Date

Method: Hand Calculation

Tech. review performed by: Lisa Noble [Signature] 12/10/15
Name Signature Date

[Signature]

SIMULATOR SOLUTION DATA ENTRY REVIEW


Reviewer/s: Amanda M. Black Date: 1-12-16

Location: WSP-FLSB Seattle, WA Solution Batch Number: 15052

	YES	NO	N/A
Analysis dates do not precede preparation date:	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Declarations signed and properly dated:	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Data entry corresponds to all chromatograms:	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
All signatures present on Test Report:	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Average solution concentration correct:	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Standard deviation correct:	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
CV (%) correct:	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Equivalent vapor concentration correct:	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
All chromatograms and sequences included in file:	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Ethanol control information present: (lot # present & used within expiration)	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Complies with accuracy and precision requirements established by the State Toxicologist:	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

Comments:

Reviewer Signature: _____



Date: 1-12-16



SOLUTION CERTIFICATE REVIEW

Please check that the data on your chromatograms is the data entered into the Test Report, that the date to the right of your name is the date that you tested the solution, and then sign the Test Report.

Please initial and date below to affirm that you have:

- 1) Checked your data
- 2) Checked the date to the right of your name on the Test Report
- 3) Signed the Test Report

	Initials	Date
Amanda Chandler		
Andrew Gingras	<i>AG</i>	12/14/15
Asa Louis		
Brittany Thomas		
Christie Mitchell-Mata	<i>CM</i>	12/11/15
Christopher Johnston		
David Nguyen		
Dawn Sklerov		
Elizabeth Wehner		
Justin Knoy	<i>JK</i>	12-11-15
Katie Harris		
Lyndsey Lowe		
Naziha Nuwayhid		
Rebecca Flaherty		

Batch # 15052 *for 12/10/15*

[Handwritten signature]

JAY INSLEE
Governor



JOHN R. BATISTE
Chief

STATE OF WASHINGTON
WASHINGTON STATE PATROL
WASHINGTON STATE TOXICOLOGY LABORATORY

2203 Airport Way South, Suite 360 • Seattle, Washington 98134-2027 • (206) 262-6100 • FAX (206) 262-6145

**0.04 g/210 L QUALITY ASSURANCE PROCEDURE SOLUTION
CERTIFICATION FOR LOT 15052**


I, Justin L. Knoy, do certify under penalty of perjury that:

I am employed by the Washington State Toxicology Laboratory, and my responsibilities include the preparation and certification of alcohol solutions for use with evidential breath test instruments.

I possess the following qualifications: BS degree in Biology, and MS degree in Forensic Science.

The quality assurance procedure (QAP) solution, Lot Number 15052, was prepared in the Washington State Toxicology Laboratory on 11/13/2015. I tested this solution and it was found to conform to those standards established by the State Toxicologist for the certification of simulator solution. It shall not be used to perform a quality assurance procedure after 11/13/2016.

Seattle, WA

 12-11-15
Justin L. Knoy Date
Forensic Scientist

JAY INSLEE
Governor



JOHN R. BATISTE
Chief

STATE OF WASHINGTON
WASHINGTON STATE PATROL
WASHINGTON STATE TOXICOLOGY LABORATORY

2203 Airport Way South, Suite 360 • Seattle, Washington 98134-2027 • (206) 262-6100 • FAX (206) 262-6145

**0.04 g/210 L QUALITY ASSURANCE PROCEDURE SOLUTION
CERTIFICATION FOR LOT 15052**


I, Andrew Gingras, do certify under penalty of perjury that:

I am employed by the Washington State Toxicology Laboratory, and my responsibilities include the preparation and certification of alcohol solutions for use with evidential breath test instruments.

I possess the following qualifications: BS degree in Cell and Molecular Biology and MS degree in Forensic Science.

The quality assurance procedure (QAP) solution, Lot Number 15052, was prepared in the Washington State Toxicology Laboratory on 11/13/2015. I tested this solution and it was found to conform to those standards established by the State Toxicologist for the certification of simulator solution. It shall not be used to perform a quality assurance procedure after 11/13/2016.

Seattle, WA

 12/14/15

Andrew Gingras Date
Forensic Scientist

JAY INSLEE
Governor



JOHN R. BATISTE
Chief

STATE OF WASHINGTON
WASHINGTON STATE PATROL
WASHINGTON STATE TOXICOLOGY LABORATORY

2203 Airport Way South, Suite 360 • Seattle, Washington 98134-2027 • (206) 262-6100 • FAX (206) 262-6145

**0.04 g/210 L QUALITY ASSURANCE PROCEDURE SOLUTION
CERTIFICATION FOR LOT 15052**


I, Christie Mitchell-Mata, do certify under penalty of perjury that:

I am employed by the Washington State Toxicology Laboratory, and my responsibilities include the preparation and certification of alcohol solutions for use with evidential breath test instruments.

I possess the following qualifications: BA degree in Chemistry, MFS degree in Forensic Science, and over four years experience in forensic toxicology.

The quality assurance procedure (QAP) solution, Lot Number 15052, was prepared in the Washington State Toxicology Laboratory on 11/13/2015. I tested this solution and it was found to conform to those standards established by the State Toxicologist for the certification of simulator solution. It shall not be used to perform a quality assurance procedure after 11/13/2016.

Seattle, WA

 12/11/2015

Christie Mitchell-Mata

Date

Forensic Scientist

FILE A COPY IN THE BATCH FILE FOR EACH SOLUTION LISTED ON THE WORKSHEET

Preparation Date: 11-13-15 Expiration Date: 11-13-16 Initials of Preparer: JKLot # of 200-proof Ethanol used in preparation: 2EA0437Date the 200-proof Ethanol bottle was opened: 11-13-15

After opening, each bottle of 200-proof Ethanol is approved for use for 6 months unless an extension is approved by the State Toxicologist. This timeframe applies to the 200-proof Ethanol only, not to simulator solutions which have a 1 year expiration.

Environmental conditions verified as acceptable:

Simulator Solution	Volume of Ethanol (mL)	Volume of Deionized Water (L)		Batch Number
QAP 0.04	11.2	18	<input checked="" type="checkbox"/>	<u>15052</u>
QAP 0.08	22.4	18	<input checked="" type="checkbox"/>	<u>15053</u>
QAP 0.10	28.1	18	<input checked="" type="checkbox"/>	<u>15054</u>
QAP 0.15	42.1	18	<input checked="" type="checkbox"/>	<u>15055</u>
QAP 0.20	56.1	18	<input checked="" type="checkbox"/>	<u>15056</u>
ESS	66.5	52	<input type="checkbox"/>	<u> </u>

Stir bar is rotating Stirred for minimum 30 minutes; 2 hours for ESS Spigot purged Aliquot taken Batch labeled, packaged and sealed 11-13-15
Date

If different ethanol lot numbers are used in the preparation of solutions, record them and the corresponding solution batch numbers in the comments section.

Comments:

Justin J. Z...
Analyst Signature

11-13-15
Date

JK

Sequence Parameters:

Operator: Justin Knoy
 Data File Naming: Prefix/Counter
 Signal 1 Prefix: SIG1
 Counter: 0001
 Signal 2 Prefix: SIG2
 Counter: 0001
 Data Directory: C:\HPCHEM\1\DATA\
 Data Subdirectory: 151113JK
 Part of Methods to run: According to Runtime Checklist
 Barcode Reader: not used
 Shutdown Cmd/Macro: none

Sequence Comment:

Ethanol Calibrator 1, E0615-01 - Exp. 12/02/2015
 Ethanol Calibrator 2, E0615-02 - Exp. 12/02/2015
 Ethanol Calibrator 3, E0615-03 - Exp. 12/02/2015
 CTRL1 (0.04g/100mL), Lot # FN05011301 - Exp. 05/2018
 CTRL2 (0.10g/100mL), Lot # FN08051301 - Exp. 10/2018
 CTRL3 (0.20g/100mL), Lot # FN03211401 - Exp. 06/2019
 Internal Standard Lot#P0915 - Exp. 12/18/2015

Calibration vials 1-9 filed with 15052.

Sequence Table (Front Injector):

Method and Injection Info Part:

Line	Location	SampleName	Method	Inj	SampleType	InjVolume	DataFile
1	Vial 1	BLANK	SIMALC1	1	Sample		
2	Vial 2	0.079 CAL 1	SIMALC1	1	Calib		
3	Vial 3	0.158 CAL 2	SIMALC1	1	Calib		
4	Vial 4	0.316 CAL 3	SIMALC1	1	Calib		
5	Vial 5	NEG CTRL	SIMALC1	1	Ctrl Samp		
6	Vial 6	0.04 CTRL	SIMALC1	1	Ctrl Samp		
7	Vial 7	0.10 CTRL	SIMALC1	1	Ctrl Samp		
8	Vial 8	0.20 CTRL	SIMALC1	1	Ctrl Samp		
9	Vial 9	NEG CTRL	SIMALC1	1	Ctrl Samp		
10	Vial 10	15052-1	SIMALC1	1	Sample		
11	Vial 11	15052-2	SIMALC1	1	Sample		
12	Vial 12	15052-3	SIMALC1	1	Sample		
13	Vial 13	15052-4	SIMALC1	1	Sample		
14	Vial 14	15052-5	SIMALC1	1	Sample		
15	Vial 15	0.10 CTRL	SIMALC1	1	Ctrl Samp		
16	Vial 16	NEG CTRL	SIMALC1	1	Ctrl Samp		
17	Vial 17	15053-1	SIMALC1	1	Sample		
18	Vial 18	15053-2	SIMALC1	1	Sample		
19	Vial 19	15053-3	SIMALC1	1	Sample		
20	Vial 20	15053-4	SIMALC1	1	Sample		
21	Vial 21	15053-5	SIMALC1	1	Sample		
22	Vial 22	0.10 CTRL	SIMALC1	1	Ctrl Samp		
23	Vial 23	NEG CTRL	SIMALC1	1	Ctrl Samp		
24	Vial 24	15054-1	SIMALC1	1	Sample		
25	Vial 25	15054-2	SIMALC1	1	Sample		
26	Vial 26	15054-3	SIMALC1	1	Sample		

15052

Justin Knoy

JK

Line	Location	SampleName	Method	Inj	SampleType	InjVolume	DataFile
27	Vial 27	15054-4	SIMALC1	1	Sample		
28	Vial 28	15054-5	SIMALC1	1	Sample		
29	Vial 29	0.10 CTRL	SIMALC1	1	Ctrl Samp		
30	Vial 30	NEG CTRL	SIMALC1	1	Ctrl Samp		
31	Vial 31	15055-1	SIMALC1	1	Sample		
32	Vial 32	15055-2	SIMALC1	1	Sample		
33	Vial 33	15055-3	SIMALC1	1	Sample		
34	Vial 34	15055-4	SIMALC1	1	Sample		
35	Vial 35	15055-5	SIMALC1	1	Sample		
36	Vial 36	0.10 CTRL	SIMALC1	1	Ctrl Samp		
37	Vial 37	NEG CTRL	SIMALC1	1	Ctrl Samp		
38	Vial 38	15056-1	SIMALC1	1	Sample		
39	Vial 39	15056-2	SIMALC1	1	Sample		
40	Vial 40	15056-3	SIMALC1	1	Sample		
41	Vial 41	15056-4	SIMALC1	1	Sample		
42	Vial 42	15056-5	SIMALC1	1	Sample		
43	Vial 43	0.10 CTRL	SIMALC1	1	Ctrl Samp		
44	Vial 44	NEG CTRL	SIMALC1	1	Ctrl Samp		

Calibration Part:

Line	Location	SampleName	Method	CalLev	Update	RF	Update	RT	Interval
2	Vial 2	0.079 CAL 1	SIMALC1	1	Replace		Replace		
3	Vial 3	0.158 CAL 2	SIMALC1	2	Replace		Replace		
4	Vial 4	0.316 CAL 3	SIMALC1	3	Replace		Replace		

Sequence Table (Back Injector):

No entries - empty table!

15052
Justin

JK

=====
Calibration Table
=====

Calib. Data Modified : Friday, November 13, 2015 11:50:48 AM

Calculate : Internal Standard
Based on : Peak Area

Rel. Reference Window : 5.000 %
Abs. Reference Window : 0.050 min
Rel. Non-ref. Window : 5.000 %
Abs. Non-ref. Window : 0.050 min
Multiplier : 1.0000
Dilution : 1.0000
Sample Amount : 0.00000

Use Multiplier & Dilution Factor with ISTDs

Uncalibrated Peaks : not reported
Partial Calibration : No recalibration if peaks missing

Curve Type : Linear
Origin : Included
Weight : Equal

Recalibration Settings:
Average Response : No Update
Average Retention Time: No Update

Calibration Report Options :
Printout of recalibrations within a sequence:
Normal Report after Recalibration

Sample ISTD Information:

ISTD #	ISTD Amount [g/100mL]	Name
1	1.20000e-2	n-Propanol

Signal 1: FID1 A,

RetTime [min]	Lvl Sig	Amount [g/100mL]	Area	Amt/Area	Ref Grp Name
1.090	1 1	7.97800e-2	918.98016	8.68136e-5	1 Ethanol
	2	1.60980e-1	1798.47375	8.95092e-5	
	3	3.18440e-1	3554.04883	8.95992e-5	
1.754	1 1	1.20000e-2	2614.82617	4.58922e-6	I1 n-Propanol
	2	1.20000e-2	2574.15698	4.66172e-6	
	3	1.20000e-2	2550.39453	4.70515e-6	


15052
Jmz1015

=====
Peak Sum Table
=====


No Entries in table
=====

JR

=====
Calibration Curves
=====



Ethanol at exp. RT: 1.090
FID1 A,
Correlation: 0.99998
Residual Std. Dev.: 0.00417
Formula: $y = mx + b$
m: 5.24566e-2
b: -2.08714e-4
x: Amount Ratio
y: Area Ratio



n-Propanol at exp. RT: 1.754
FID1 A,
Correlation: 1.00000
Residual Std. Dev.: 0.00000
Formula: $y = mx + b$
m: 1.00000
b: 0.00000
x: Amount Ratio
y: Area Ratio

=====
15052
Inzholak

Washington State Patrol Toxicology Laboratory
2203 Airport Way S Seattle, WA 98134

Inj. Date: 11/13/2015 11:23:29 AM

Sample Name: BLANK

Instrument: HSGC#1


Operator: Justin Knoy

Column: DB-ALC1


Location: Vial 1

Method: C:\HPCHEM\1\METHODS\SIMALC1.M

Sample Info: 15052




#	Compound	Peak Area	RT (min)
1	Ethanol	0	0.000
2	n-Propanol	0	0.000



Ethanol 0.000 g/100mL

JK



n-Propanol 0.000 g/100mL

JK

Washington State Patrol Toxicology Laboratory
2203 Airport Way S Seattle, WA 98134

Inj. Date: 11/13/2015 11:26:49 AM

Sample Name: 0.079 CAL 1

Instrument: HSGC#1


Operator: Justin Knoy

Column: DB-ALC1


Location: Vial 2

Method: C:\HPCHEM\1\METHODS\SIMALC1.M

Sample Info: 15052




#	Compound	Peak Area	RT (min)
1	Ethanol	919	1.090
2	n-Propanol	2615	1.754



Ethanol 0.080 g/100mL

JK



n-Propanol 0.012 g/100mL

JK

Washington State Patrol Toxicology Laboratory
 2203 Airport Way S Seattle, WA 98134

Inj. Date: 11/13/2015 11:30:04 AM

Sample Name: 0.158 CAL 2

Instrument: HSGC#1


Operator: Justin Knoy

Column: DB-ALC1


Location: Vial 3

Method: C:\HPCHEM\1\METHODS\SIMALC1.M

Sample Info: 15052




#	Compound	Peak Area	RT (min)
1	Ethanol	1798	1.091
2	n-Propanol	2574	1.755



Ethanol 0.160 g/100mL

Handwritten signature




n-Propanol 0.012 g/100mL

Handwritten initials


Washington State Patrol Toxicology Laboratory
2203 Airport Way S Seattle, WA 98134

Inj. Date: 11/13/2015 11:33:21 AM
Instrument: HSGC#1
Column: DB-ALC1
Method: C:\HPCHEM\1\METHODS\SIMALC1.M
Sample Info: 15052


Sample Name: 0.316 CAL 3
Operator: Justin Knoy
Location: Vial 4



#	Compound	Peak Area	RT (min)
1	Ethanol	3554	1.090
2	n-Propanol	2550	1.754



Ethanol 0.319 g/100mL



n-Propanol 0.012 g/100mL


Washington State Patrol Toxicology Laboratory
2203 Airport Way S Seattle, WA 98134

Inj. Date: 11/13/2015 11:36:35 AM
Instrument: HSGC#1


Sample Name: NEG CTRL
Operator: Justin Knoy
Location: Vial 5

Column: DB-ALC1
Method: C:\HPCHEM\1\METHODS\SIMALC1.M


Sample Info: 15052



#	Compound	Peak Area	RT (min)
1	Ethanol	0	0.000
2	n-Propanol	2576	1.753



Ethanol 0.000 g/100mL



n-Propanol 0.012 g/100mL

JK

JK

Washington State Patrol Toxicology Laboratory
 2203 Airport Way S Seattle, WA 98134

Inj. Date: 11/13/2015 11:39:48 AM

Sample Name: 0.04 CTRL

Instrument: HSGC#1


Operator: Justin Knoy

Column: DB-ALC1


Location: Vial 6

Method: C:\HPCHEM\1\METHODS\SIMALC1.M


Sample Info: 15052



#	Compound	Peak Area	RT (min)
1	Ethanol	465	1.090
2	n-Propanol	2599	1.754



Ethanol 0.041 g/100mL




n-Propanol 0.012 g/100mL

JK


JK

Washington State Patrol Toxicology Laboratory
2203 Airport Way S Seattle, WA 98134

Inj. Date: 11/13/2015 11:43:01 AM Sample Name: 0.10 CTRL
Instrument: HSGC#1 Operator: Justin Knoy
Column: DB-ALC1 Location: Vial 7
Method: C:\HPCHEM\1\METHODS\SIMALC1.M
Sample Info: 15052




#	Compound	Peak Area	RT (min)
1	Ethanol	1130	1.091
2	n-Propanol	2565	1.754



Ethanol 0.101 g/100mL

fr



n-Propanol 0.012 g/100mL

JV

Washington State Patrol Toxicology Laboratory
 2203 Airport Way S Seattle, WA 98134

Inj. Date: 11/13/2015 11:46:15 AM

Sample Name: 0.20 CTRL

Instrument: HSGC#1


Operator: Justin Knoy

Column: DB-ALC1


Location: Vial 8

Method: C:\HPCHEM\1\METHODS\SIMALC1.M

Sample Info: 15052




#	Compound	Peak Area	RT (min)
1	Ethanol	2300	1.091
2	n-Propanol	2591	1.755



Ethanol 0.203 g/100mL

Handwritten signature




n-Propanol 0.012 g/100mL


Handwritten initials JTK

Washington State Patrol Toxicology Laboratory
2203 Airport Way S Seattle, WA 98134

Inj. Date: 11/13/2015 11:49:27 AM Sample Name: NEG CTRL
Instrument: HSGC#1 Operator: Justin Knoy
Column: DB-ALC1 Location: Vial 9
Method: C:\HPCHEM\1\METHODS\SIMALC1.M
Sample Info: 15052




#	Compound	Peak Area	RT (min)
1	Ethanol	0	0.000
2	n-Propanol	2572	1.753



Ethanol 0.000 g/100mL

Handwritten signature



n-Propanol 0.012 g/100mL

Handwritten initials

Washington State Patrol Toxicology Laboratory
 2203 Airport Way S Seattle, WA 98134

Inj. Date: 11/13/2015 11:52:39 AM

Sample Name: 15052-1

Instrument: HSGC#1


Operator: Justin Knoy

Column: DB-ALC1


Location: Vial 10

Method: C:\HPCHEM\1\METHODS\SIMALC1.M

Sample Info:




#	Compound	Peak Area	RT (min)
1	Ethanol	578	1.091
2	n-Propanol	2614	1.755



Ethanol 0.051 g/100mL

JK



n-Propanol 0.012 g/100mL

JK

Washington State Patrol Toxicology Laboratory
 2203 Airport Way S Seattle, WA 98134

Inj. Date: 11/13/2015 11:55:51 AM

Sample Name: 15052-2

Instrument: HSGC#1


Operator: Justin Knoy

Column: DB-ALC1


Location: Vial 11

Method: C:\HPCHEM\1\METHODS\SIMALC1.M


Sample Info:



#	Compound	Peak Area	RT (min)
1	Ethanol	577	1.091
2	n-Propanol	2605	1.755



Ethanol 0.051 g/100mL



n-Propanol 0.012 g/100mL

JK

JK

Washington State Patrol Toxicology Laboratory
 2203 Airport Way S Seattle, WA 98134

Inj. Date: 11/13/2015 11:59:04 AM

Sample Name: 15052-3

Instrument: HSGC#1


Operator: Justin Knoy

Column: DB-ALC1


Location: Vial 12

Method: C:\HPCHEM\1\METHODS\SIMALC1.M


Sample Info:



#	Compound	Peak Area	RT (min)
1	Ethanol	578	1.092
2	n-Propanol	2605	1.756



Ethanol 0.051 g/100mL



n-Propanol 0.012 g/100mL

JK

JK

Washington State Patrol Toxicology Laboratory
2203 Airport Way S Seattle, WA 98134

Inj. Date: 11/13/2015 12:02:17 PM

Sample Name: 15052-4

Instrument: HSGC#1


Operator: Justin Knoy

Column: DB-ALC1


Location: Vial 13

Method: C:\HPCHEM\1\METHODS\SIMALC1.M


Sample Info:



#	Compound	Peak Area	RT (min)
1	Ethanol	584	1.091
2	n-Propanol	2630	1.755



Ethanol 0.051 g/100mL



n-Propanol 0.012 g/100mL

Handwritten signature


Handwritten initials

Washington State Patrol Toxicology Laboratory
 2203 Airport Way S Seattle, WA 98134


Inj. Date: 11/13/2015 12:05:31 PM
 Instrument: HSGC#1
 Column: DB-ALC1
 Method: C:\HPCHEM\1\METHODS\SIMALC1.M

Sample Name: 15052-5
 Operator: Justin Knoy
 Location: Vial 14


Sample Info:



#	Compound	Peak Area	RT (min)
1	Ethanol	582	1.092
2	n-Propanol	2625	1.756



Ethanol 0.051 g/100mL



n-Propanol 0.012 g/100mL


JK

JK


Washington State Patrol Toxicology Laboratory
 2203 Airport Way S Seattle, WA 98134

Inj. Date: 11/13/2015 12:08:44 PM
 Instrument: HSGC#1
 Column: DB-ALC1
 Method: C:\HPCHEM\1\METHODS\SIMALC1.M
 Sample Info: 15052


Sample Name: 0.10 CTRL
 Operator: Justin Knoy
 Location: Vial 15



#	Compound	Peak Area	RT (min)
1	Ethanol	1157	1.091
2	n-Propanol	2611	1.755



Ethanol 0.101 g/100mL




n-Propanol 0.012 g/100mL

Handwritten signature


Handwritten mark

Washington State Patrol Toxicology Laboratory
2203 Airport Way S Seattle, WA 98134


Inj. Date: 11/13/2015 12:11:57 PM Sample Name: NEG CTRL
Instrument: HSGC#1 Operator: Justin Knoy
Column: DB-ALC1 Location: Vial 16
Method: C:\HPCHEM\1\METHODS\SIMALC1.M
Sample Info: 15052



#	Compound	Peak Area	RT (min)
1	Ethanol	0	0.000
2	n-Propanol	2574	1.753



Ethanol 0.000 g/100mL



n-Propanol 0.012 g/100mL

JK

JK

Sequence Parameters:

Operator: Andrew Gingras
 Data File Naming: Prefix/Counter
 Signal 1 Prefix: SIG1
 Counter: 0001
 Signal 2 Prefix: SIG2
 Counter: 0001
 Data Directory: C:\HPCHEM\1\DATA\
 Data Subdirectory: 151113AG
 Part of Methods to run: According to Runtime Checklist
 Barcode Reader: not used
 Shutdown Cmd/Macro: none

Sequence Comment:

Ethanol Calibrator 1, E0615-01 - Exp. 12/02/2015
 Ethanol Calibrator 2, E0615-02 - Exp. 12/02/2015
 Ethanol Calibrator 3, E0615-03 - Exp. 12/02/2015
 CTRL1 (0.04g/100mL), Lot # FN05011301 - Exp. 05/2018
 CTRL2 (0.10g/100mL), Lot # FN08051301 - Exp. 10/2018
 CTRL3 (0.20g/100mL), Lot # FN03211401 - Exp. 06/2019
 Internal Standard Lot#P0915 - Exp. 12/18/2015

Calibration vials 1-9 filed with 15052.

Sequence Table (Front Injector):

Method and Injection Info Part:

Line	Location	SampleName	Method	Inj	SampleType	InjVolume	DataFile
1	Vial 1	BLANK	SIMALC1	1	Sample		
2	Vial 2	0.079 CAL 1	SIMALC1	1	Calib		
3	Vial 3	0.158 CAL 2	SIMALC1	1	Calib		
4	Vial 4	0.316 CAL 3	SIMALC1	1	Calib		
5	Vial 5	NEG CTRL	SIMALC1	1	Ctrl Samp		
6	Vial 6	0.04 CTRL	SIMALC1	1	Ctrl Samp		
7	Vial 7	0.10 CTRL	SIMALC1	1	Ctrl Samp		
8	Vial 8	0.20 CTRL	SIMALC1	1	Ctrl Samp		
9	Vial 9	NEG CTRL	SIMALC1	1	Ctrl Samp		
10	Vial 10	15052-1	SIMALC1	1	Sample		
11	Vial 11	15052-2	SIMALC1	1	Sample		
12	Vial 12	15052-3	SIMALC1	1	Sample		
13	Vial 13	15052-4	SIMALC1	1	Sample		
14	Vial 14	15052-5	SIMALC1	1	Sample		
15	Vial 15	0.10 CTRL	SIMALC1	1	Ctrl Samp		
16	Vial 16	NEG CTRL	SIMALC1	1	Ctrl Samp		
17	Vial 17	15053-1	SIMALC1	1	Sample		
18	Vial 18	15053-2	SIMALC1	1	Sample		
19	Vial 19	15053-3	SIMALC1	1	Sample		
20	Vial 20	15053-4	SIMALC1	1	Sample		
21	Vial 21	15053-5	SIMALC1	1	Sample		
22	Vial 22	0.10 CTRL	SIMALC1	1	Ctrl Samp		
23	Vial 23	NEG CTRL	SIMALC1	1	Ctrl Samp		
24	Vial 24	15054-1	SIMALC1	1	Sample		
25	Vial 25	15054-2	SIMALC1	1	Sample		
26	Vial 26	15054-3	SIMALC1	1	Sample		

15052

Injections

Line	Location	SampleName	Method	Inj	SampleType	InjVolume	DataFile
27	Vial 27	15054-4	SIMALC1	1	Sample		
28	Vial 28	15054-5	SIMALC1	1	Sample		
29	Vial 29	0.10 CTRL	SIMALC1	1	Ctrl Samp		
30	Vial 30	NEG CTRL	SIMALC1	1	Ctrl Samp		
31	Vial 31	15055-1	SIMALC1	1	Sample		
32	Vial 32	15055-2	SIMALC1	1	Sample		
33	Vial 33	15055-3	SIMALC1	1	Sample		
34	Vial 34	15055-4	SIMALC1	1	Sample		
35	Vial 35	15055-5	SIMALC1	1	Sample		
36	Vial 36	0.10 CTRL	SIMALC1	1	Ctrl Samp		
37	Vial 37	NEG CTRL	SIMALC1	1	Ctrl Samp		
38	Vial 38	15056-1	SIMALC1	1	Sample		
39	Vial 39	15056-2	SIMALC1	1	Sample		
40	Vial 40	15056-3	SIMALC1	1	Sample		
41	Vial 41	15056-4	SIMALC1	1	Sample		
42	Vial 42	15056-5	SIMALC1	1	Sample		
43	Vial 43	0.10 CTRL	SIMALC1	1	Ctrl Samp		
44	Vial 44	NEG CTRL	SIMALC1	1	Ctrl Samp		

Calibration Part:

Line	Location	SampleName	Method	CalLev	Update	RF	Update	RT	Interval
2	Vial 2	0.079 CAL 1	SIMALC1	1	Replace		Replace		
3	Vial 3	0.158 CAL 2	SIMALC1	2	Replace		Replace		
4	Vial 4	0.316 CAL 3	SIMALC1	3	Replace		Replace		

Sequence Table (Back Injector):

No entries - empty table!

15052
12/10/15

=====
 Calibration Table
 =====

Calib. Data Modified : Friday, November 13, 2015 2:13:54 PM

Calculate : Internal Standard
 Based on : Peak Area

Rel. Reference Window : 5.000 %
 Abs. Reference Window : 0.050 min
 Rel. Non-ref. Window : 5.000 %
 Abs. Non-ref. Window : 0.050 min
 Multiplier : 1.0000
 Dilution : 1.0000
 Sample Amount : 0.00000

Use Multiplier & Dilution Factor with ISTDs
 Uncalibrated Peaks : not reported
 Partial Calibration : No recalibration if peaks missing

Curve Type : Linear
 Origin : Included
 Weight : Equal

Recalibration Settings:
 Average Response : No Update
 Average Retention Time: No Update

Calibration Report Options :
 Printout of recalibrations within a sequence:
 Normal Report after Recalibration

Sample ISTD Information:

ISTD #	ISTD Amount [g/100mL]	Name
1	1.20000e-2	n-Propanol

Signal 1: FID1 A,


RetTime [min]	Lvl Sig	Amount [g/100mL]	Area	Amt/Area	Ref	Grp Name
1.092	1 1	7.97800e-2	931.52808	8.56442e-5	1	Ethanol
		1.60980e-1	1829.50879	8.79908e-5		
		3.18440e-1	3647.98975	8.72919e-5		
1.755	1 1	1.20000e-2	2640.18164	4.54514e-6	I1	n-Propanol
		1.20000e-2	2612.95264	4.59251e-6		
		1.20000e-2	2612.21509	4.59380e-6		

=====
 Peak Sum Table
 =====


No Entries in table
 =====

15052
 L
 Muzikis

=====
Calibration Curves
=====



Ethanol at exp. RT: 1.092
FID1 A,
Correlation: 0.99998
Residual Std. Dev.: 0.00440
Formula: $y = mx + b$
m: 5.25608e-2
b: 4.32853e-5
x: Amount Ratio
y: Area Ratio




n-Propanol at exp. RT: 1.755
FID1 A,
Correlation: 1.00000
Residual Std. Dev.: 0.00000
Formula: $y = mx + b$
m: 1.00000
b: 0.00000
x: Amount Ratio
y: Area Ratio

=====
15052


Included

Washington State Patrol Toxicology Laboratory
2203 Airport Way S Seattle, WA 98134

Inj. Date: 11/13/2015 2:01:48 PM Sample Name: BLANK
Instrument: HSGC#1 Operator: Andrew Gingras
Column: DB-ALC1 Location: Vial 1
Method: C:\HPCHEM\1\METHODS\SIMALC1.M
Sample Info: 15052




#	Compound	Peak Area	RT (min)
1	Ethanol	0	0.000
2	n-Propanol	0	0.000



Ethanol 0.000 g/100mL

fn



n-Propanol 0.000 g/100mL

6

Washington State Patrol Toxicology Laboratory
 2203 Airport Way S Seattle, WA 98134

Inj. Date: 11/13/2015 2:05:07 PM

Sample Name: 0.079 CAL 1

Instrument: HSGC#1


Operator: Andrew Gingras

Column: DB-ALC1


Location: Vial 2

Method: C:\HPCHEM\1\METHODS\SIMALC1.M


Sample Info: 15052



#	Compound	Peak Area	RT (min)
1	Ethanol	932	1.092
2	n-Propanol	2640	1.756



Ethanol 0.081 g/100mL




n-Propanol 0.012 g/100mL

Handwritten mark


Handwritten signature

Washington State Patrol Toxicology Laboratory
2203 Airport Way S Seattle, WA 98134


Inj. Date: 11/13/2015 2:08:24 PM Sample Name: 0.158 CAL 2
Instrument: HSGC#1 Operator: Andrew Gingras
Column: DB-ALC1 Location: Vial 3
Method: C:\HPCHEM\1\METHODS\SIMALC1.M
Sample Info: 15052



#	Compound	Peak Area	RT (min)
1	Ethanol	1830	1.092
2	n-Propanol	2613	1.756



Ethanol 0.160 g/100mL



n-Propanol 0.012 g/100mL

Handwritten signature

Handwritten signature


Washington State Patrol Toxicology Laboratory
 2203 Airport Way S Seattle, WA 98134

Inj. Date: 11/13/2015 2:11:41 PM
 Instrument: HSGC#1


Sample Name: 0.316 CAL 3
 Operator: Andrew Gingras
 Location: Vial 4

Column: DB-ALC1
 Method: C:\HPCHEM\1\METHODS\SIMALC1.M


Sample Info: 15052



#	Compound	Peak Area	RT (min)
1	Ethanol	3648	1.092
2	n-Propanol	2612	1.755



Ethanol 0.319 g/100mL



n-Propanol 0.012 g/100mL

Handwritten signature

Handwritten signature


Washington State Patrol Toxicology Laboratory
2203 Airport Way S Seattle, WA 98134

Inj. Date: 11/13/2015 2:14:55 PM
Instrument: HSGC#1


Sample Name: NEG CTRL
Operator: Andrew Gingras
Location: Vial 5

Column: DB-ALC1
Method: C:\HPCHEM\1\METHODS\SIMALC1.M


Sample Info: 15052



#	Compound	Peak Area	RT (min)
1	Ethanol	0	0.000
2	n-Propanol	2604	1.755



Ethanol 0.000 g/100mL



n-Propanol 0.012 g/100mL


Handwritten signature

Handwritten signature


Washington State Patrol Toxicology Laboratory
 2203 Airport Way S Seattle, WA 98134

Inj. Date: 11/13/2015 2:18:08 PM
 Instrument: HSGC#1
 Column: DB-ALC1
 Method: C:\HPCHEM\1\METHODS\SIMALC1.M
 Sample Info: 15052


Sample Name: 0.04 CTRL
 Operator: Andrew Gingras
 Location: Vial 6



#	Compound	Peak Area	RT (min)
1	Ethanol	471	1.092
2	n-Propanol	2622	1.755



Ethanol 0.041 g/100mL



n-Propanol 0.012 g/100mL

Handwritten initials

Handwritten signature


Washington State Patrol Toxicology Laboratory
 2203 Airport Way S Seattle, WA 98134

Inj. Date: 11/13/2015 2:21:21 PM
 Instrument: HSGC#1


Sample Name: 0.10 CTRL
 Operator: Andrew Gingras
 Location: Vial 7

Column: DB-ALC1
 Method: C:\HPCHEM\1\METHODS\SIMALC1.M


Sample Info: 15052



#	Compound	Peak Area	RT (min)
1	Ethanol	1139	1.091
2	n-Propanol	2564	1.754



Ethanol 0.101 g/100mL



n-Propanol 0.012 g/100mL

J

JC

Washington State Patrol Toxicology Laboratory
 2203 Airport Way S Seattle, WA 98134

Inj. Date: 11/13/2015 2:24:34 PM

Sample Name: 0.20 CTRL

Instrument: HSGC#1


Operator: Andrew Gingras

Column: DB-ALC1


Location: Vial 8

Method: C:\HPCHEM\1\METHODS\SIMALC1.M


Sample Info: 15052



#	Compound	Peak Area	RT (min)
1	Ethanol	2256	1.089
2	n-Propanol	2545	1.753



Ethanol 0.202 g/100mL



n-Propanol 0.012 g/100mL

Handwritten mark

Handwritten signature

Washington State Patrol Toxicology Laboratory
2203 Airport Way S Seattle, WA 98134

Inj. Date: 11/13/2015 2:27:48 PM

Sample Name: NEG CTRL

Instrument: HSGC#1


Operator: Andrew Gingras

Column: DB-ALC1


Location: Vial 9

Method: C:\HPCHEM\1\METHODS\SIMALC1.M

Sample Info: 15052




#	Compound	Peak Area	RT (min)
1	Ethanol	0	0.000
2	n-Propanol	2624	1.755



Ethanol 0.000 g/100mL

Handwritten signature



n-Propanol 0.012 g/100mL


Handwritten signature

Washington State Patrol Toxicology Laboratory
 2203 Airport Way S Seattle, WA 98134


Inj. Date: 11/13/2015 2:31:01 PM
 Instrument: HSGC#1
 Column: DB-ALC1
 Method: C:\HPCHEM\1\METHODS\SIMALC1.M

Sample Name: 15052-1
 Operator: Andrew Gingras
 Location: Vial 10

Sample Info:




#	Compound	Peak Area	RT (min)
1	Ethanol	580	1.091
2	n-Propanol	2624	1.755



Ethanol 0.050 g/100mL

Handwritten signature



n-Propanol 0.012 g/100mL

Handwritten signature


Washington State Patrol Toxicology Laboratory
2203 Airport Way S Seattle, WA 98134

Inj. Date: 11/13/2015 2:34:14 PM
Instrument: HSGC#1
Column: DB-ALC1


Sample Name: 15052-2
Operator: Andrew Gingras
Location: Vial 11

Method: C:\HPCHEM\1\METHODS\SIMALC1.M

Sample Info:




#	Compound	Peak Area	RT (min)
1	Ethanol	581	1.092
2	n-Propanol	2629	1.756



Ethanol 0.050 g/100mL

Handwritten signature



n-Propanol 0.012 g/100mL

Handwritten signature

Washington State Patrol Toxicology Laboratory
 2203 Airport Way S Seattle, WA 98134

Inj. Date: 11/13/2015 2:37:27 PM

Sample Name: 15052-3

Instrument: HSGC#1


Operator: Andrew Gingras

Column: DB-ALC1


Location: Vial 12

Method: C:\HPCHEM\1\METHODS\SIMALC1.M

Sample Info:




#	Compound	Peak Area	RT (min)
1	Ethanol	583	1.092
2	n-Propanol	2637	1.756



Ethanol 0.050 g/100mL

Handwritten signature



n-Propanol 0.012 g/100mL

Handwritten signature

Washington State Patrol Toxicology Laboratory
 2203 Airport Way S Seattle, WA 98134

Inj. Date: 11/13/2015 2:40:41 PM

Sample Name: 15052-4

Instrument: HSGC#1


Operator: Andrew Gingras

Column: DB-ALC1


Location: Vial 13

Method: C:\HPCHEM\1\METHODS\SIMALC1.M


Sample Info:



#	Compound	Peak Area	RT (min)
1	Ethanol	585	1.091
2	n-Propanol	2644	1.755



Ethanol 0.050 g/100mL



n-Propanol 0.012 g/100mL

Handwritten signature

Handwritten signature

Washington State Patrol Toxicology Laboratory
2203 Airport Way S Seattle, WA 98134

Inj. Date: 11/13/2015 2:43:56 PM

Sample Name: 15052-5

Instrument: HSGC#1


Operator: Andrew Gingras

Column: DB-ALC1


Location: Vial 14

Method: C:\HPCHEM\1\METHODS\SIMALC1.M


Sample Info:



#	Compound	Peak Area	RT (min)
1	Ethanol	586	1.092
2	n-Propanol	2642	1.756



Ethanol 0.051 g/100mL



n-Propanol 0.012 g/100mL

Handwritten signature

Handwritten signature

Washington State Patrol Toxicology Laboratory
 2203 Airport Way S Seattle, WA 98134

Inj. Date: 11/13/2015 2:47:07 PM

Sample Name: 0.10 CTRL

Instrument: HSGC#1


Operator: Andrew Gingras

Column: DB-ALC1


Location: Vial 15

Method: C:\HPCHEM\1\METHODS\SIMALC1.M


Sample Info: 15052



#	Compound	Peak Area	RT (min)
1	Ethanol	1178	1.092
2	n-Propanol	2630	1.756



Ethanol 0.102 g/100mL




n-Propanol 0.012 g/100mL

Handwritten signature


Handwritten signature

Washington State Patrol Toxicology Laboratory
2203 Airport Way S Seattle, WA 98134

Inj. Date: 11/13/2015 2:50:20 PM Sample Name: NEG CTRL
Instrument: HSGC#1 Operator: Andrew Gingras
Column: DB-ALC1 Location: Vial 16
Method: C:\HPCHEM\1\METHODS\SIMALC1.M
Sample Info: 15052




#	Compound	Peak Area	RT (min)
1	Ethanol	0	0.000
2	n-Propanol	2612	1.755



Ethanol 0.000 g/100mL

Handwritten signature



n-Propanol 0.012 g/100mL

Handwritten signature

Sequence Parameters:

Operator: Christie Mitchell-Mata

Data File Naming: Prefix/Counter

Signal 1 Prefix: SIG1

Counter: 0001

Signal 2 Prefix: SIG2

Counter: 0001

Data Directory: C:\HPCHEM\1\DATA\

Data Subdirectory: 151118CM

Part of Methods to run: According to Runtime Checklist

Barcode Reader: not used

Shutdown Cmd/Macro: none

Sequence Comment:

Cal 1 (0.079 g/100mL) - Lot#E1015-01 - Exp 04/29/16

Cal 2 (0.158 g/100mL) - Lot#E1015-02 - Exp 04/29/16

Cal 3 (0.316 g/100mL) - Lot#E1015-03 - Exp 04/29/16

CTRL 1 (0.04 g/100mL) - Lot#FN05011301 - Exp 05/2018

CTRL 2 (0.10 g/100mL) - Lot#FN08051301 - Exp 10/2018

CTRL 3 (0.20 g/100mL) - Lot#FN03211401 - Exp 06/2019

n-Propanol ISTD - Lot# P0915 - Exp 12/18/15

Calibration vials 1-9 filed with 15052.

Sequence Table (Front Injector):

Method and Injection Info Part:

Line	Location	SampleName	Method	Inj	SampleType	InjVolume	DataFile
1	Vial 1	BLANK	SIMALC1	1	Sample		
2	Vial 2	0.079 CAL 1	SIMALC1	1	Calib		
3	Vial 3	0.158 CAL 2	SIMALC1	1	Calib		
4	Vial 4	0.316 CAL 3	SIMALC1	1	Calib		
5	Vial 5	NEG CTRL	SIMALC1	1	Ctrl Samp		
6	Vial 6	0.04 CTRL	SIMALC1	1	Ctrl Samp		
7	Vial 7	0.10 CTRL	SIMALC1	1	Ctrl Samp		
8	Vial 8	0.20 CTRL	SIMALC1	1	Ctrl Samp		
9	Vial 9	NEG CTRL	SIMALC1	1	Ctrl Samp		
10	Vial 10	15052-1	SIMALC1	1	Sample		
11	Vial 11	15052-2	SIMALC1	1	Sample		
12	Vial 12	15052-3	SIMALC1	1	Sample		
13	Vial 13	15052-4	SIMALC1	1	Sample		
14	Vial 14	15052-5	SIMALC1	1	Sample		
15	Vial 15	0.10 CTRL	SIMALC1	1	Ctrl Samp		
16	Vial 16	NEG CTRL	SIMALC1	1	Ctrl Samp		
17	Vial 17	15053-1	SIMALC1	1	Sample		
18	Vial 18	15053-2	SIMALC1	1	Sample		
19	Vial 19	15053-3	SIMALC1	1	Sample		
20	Vial 20	15053-4	SIMALC1	1	Sample		
21	Vial 21	15053-5	SIMALC1	1	Sample		
22	Vial 22	0.10 CTRL	SIMALC1	1	Ctrl Samp		
23	Vial 23	NEG CTRL	SIMALC1	1	Ctrl Samp		
24	Vial 24	15054-1	SIMALC1	1	Sample		
25	Vial 25	15054-2	SIMALC1	1	Sample		
26	Vial 26	15054-3	SIMALC1	1	Sample		

15052

Injections

u

Line	Location	SampleName	Method	Inj	SampleType	InjVolume	DataFile
27	Vial 27	15054-4	SIMALC1	1	Sample		
28	Vial 28	15054-5	SIMALC1	1	Sample		
29	Vial 29	0.10 CTRL	SIMALC1	1	Ctrl Samp		
30	Vial 30	NEG CTRL	SIMALC1	1	Ctrl Samp		
31	Vial 31	15055-1	SIMALC1	1	Sample		
32	Vial 32	15055-2	SIMALC1	1	Sample		
33	Vial 33	15055-3	SIMALC1	1	Sample		
34	Vial 34	15055-4	SIMALC1	1	Sample		
35	Vial 35	15055-5	SIMALC1	1	Sample		
36	Vial 36	0.10 CTRL	SIMALC1	1	Ctrl Samp		
37	Vial 37	NEG CTRL	SIMALC1	1	Ctrl Samp		
38	Vial 38	15056-1	SIMALC1	1	Sample		
39	Vial 39	15056-2	SIMALC1	1	Sample		
40	Vial 40	15056-3	SIMALC1	1	Sample		
41	Vial 41	15056-4	SIMALC1	1	Sample		
42	Vial 42	15056-5	SIMALC1	1	Sample		
43	Vial 43	0.10 CTRL	SIMALC1	1	Ctrl Samp		
44	Vial 44	NEG CTRL	SIMALC1	1	Ctrl Samp		

Calibration Part:

Line	Location	SampleName	Method	CalLev	Update	RF	Update	RT	Interval
2	Vial 2	0.079 CAL 1	SIMALC1	1	Replace		Replace		
3	Vial 3	0.158 CAL 2	SIMALC1	2	Replace		Replace		
4	Vial 4	0.316 CAL 3	SIMALC1	3	Replace		Replace		

Sequence Table (Back Injector):

No entries - empty table!

15052

Smith

CM 11/18/15
15111807

u

=====
 Calibration Table
 =====

Calib. Data Modified : Wednesday, November 18, 2015 8:03:28 AM

Calculate : Internal Standard
 Based on : Peak Area

Rel. Reference Window : 5.000 %
 Abs. Reference Window : 0.050 min
 Rel. Non-ref. Window : 5.000 %
 Abs. Non-ref. Window : 0.050 min
 Multiplier : 1.0000
 Dilution : 1.0000
 Sample Amount : 0.00000

Use Multiplier & Dilution Factor with ISTDs

Uncalibrated Peaks : not reported
 Partial Calibration : No recalibration if peaks missing

Curve Type : Linear
 Origin : Included
 Weight : Equal

Recalibration Settings:
 Average Response : No Update
 Average Retention Time: No Update

Calibration Report Options :
 Printout of recalibrations within a sequence:
 Normal Report after Recalibration

Sample ISTD Information:

ISTD #	ISTD Amount [g/100mL]	Name
1	1.20000e-2	n-Propanol

Signal 1: FID1 A,

RetTime [min]	Lvl Sig	Amount [g/100mL]	Area	Amt/Area	Ref Grp Name
1.091	1 1	7.96600e-2	902.57788	8.82583e-5	1 Ethanol
		2 1.61370e-1	1938.90002	8.32276e-5	
		3 3.22930e-1	3567.88892	9.05101e-5	
1.754	1 1	1.20000e-2	2546.05176	4.71318e-6	I1 n-Propanol
		2 1.20000e-2	2628.15601	4.56594e-6	
		3 1.20000e-2	2512.61621	4.77590e-6	

15052

Christie Mitchell

=====
 Peak Sum Table
 =====

No Entries in table
 =====

CM


Washington State Patrol Toxicology Laboratory
2203 Airport Way S Seattle, WA 98134

Inj. Date: 11/18/2015 7:51:22 AM
Instrument: HSGC#1


Sample Name: BLANK
Operator: Christie Mitchell-Mata
Location: Vial 1

Column: DB-ALC1
Method: C:\HPCHEM\1\METHODS\SIMALC1.M


Sample Info: 15052



#	Compound	Peak Area	RT (min)
1	Ethanol	0	0.000
2	n-Propanol	0	0.000



Ethanol 0.000 g/100mL



n-Propanol 0.000 g/100mL

fm

am

Washington State Patrol Toxicology Laboratory
 2203 Airport Way S Seattle, WA 98134

Inj. Date: 11/18/2015 7:54:41 AM

Sample Name: 0.079 CAL 1

Instrument: HSGC#1


Operator: Christie Mitchell-Mata

Column: DB-ALC1


Location: Vial 2

Method: C:\HPCHEM\1\METHODS\SIMALC1.M

Sample Info: 15052




#	Compound	Peak Area	RT (min)
1	Ethanol	903	1.091
2	n-Propanol	2546	1.754



Ethanol 0.079 g/100mL

Handwritten mark



n-Propanol 0.012 g/100mL

Handwritten mark


Washington State Patrol Toxicology Laboratory
 2203 Airport Way S Seattle, WA 98134

Inj. Date: 11/18/2015 7:57:57 AM
 Instrument: HSGC#1


Sample Name: 0.158 CAL 2
 Operator: Christie Mitchell-Mata
 Location: Vial 3

Column: DB-ALC1
 Method: C:\HPCHEM\1\METHODS\SIMALC1.M

Sample Info: 15052




#	Compound	Peak Area	RT (min)
1	Ethanol	1939	1.093
2	n-Propanol	2628	1.756



Ethanol 0.166 g/100mL

Handwritten mark



n-Propanol 0.012 g/100mL

Handwritten mark

Washington State Patrol Toxicology Laboratory
 2203 Airport Way S Seattle, WA 98134

Inj. Date: 11/18/2015 8:01:14 AM

Sample Name: 0.316 CAL 3

Instrument: HSGC#1


Operator: Christie Mitchell-Mata

Column: DB-ALC1


Location: Vial 4

Method: C:\HPCHEM\1\METHODS\SIMALC1.M

Sample Info: 15052




#	Compound	Peak Area	RT (min)
1	Ethanol	3568	1.091
2	n-Propanol	2513	1.754



Ethanol 0.321 g/100mL

h



n-Propanol 0.012 g/100mL

u

Washington State Patrol Toxicology Laboratory
2203 Airport Way S Seattle, WA 98134

Inj. Date: 11/18/2015 8:04:28 AM

Sample Name: NEG CTRL

Instrument: HSGC#1


Operator: Christie Mitchell-Mata

Column: DB-ALC1


Location: Vial 5

Method: C:\HPCHEM\1\METHODS\SIMALC1.M

Sample Info: 15052




#	Compound	Peak Area	RT (min)
1	Ethanol	0	0.000
2	n-Propanol	2554	1.754



Ethanol 0.000 g/100mL

Handwritten signature



n-Propanol 0.012 g/100mL

Handwritten signature

Washington State Patrol Toxicology Laboratory
2203 Airport Way S Seattle, WA 98134

Inj. Date: 11/18/2015 8:07:41 AM

Sample Name: 0.04 CTRL

Instrument: HSGC#1


Operator: Christie Mitchell-Mata

Column: DB-ALC1


Location: Vial 6

Method: C:\HPCHEM\1\METHODS\SIMALC1.M


Sample Info: 15052



#	Compound	Peak Area	RT (min)
1	Ethanol	459	1.091
2	n-Propanol	2551	1.755



Ethanol 0.039 g/100mL



n-Propanol 0.012 g/100mL

Handwritten mark

Handwritten mark

Washington State Patrol Toxicology Laboratory
 2203 Airport Way S Seattle, WA 98134

Inj. Date: 11/18/2015 8:10:55 AM

Sample Name: 0.10 CTRL

Instrument: HSGC#1


Operator: Christie Mitchell-Mata

Column: DB-ALC1


Location: Vial 7

Method: C:\HPCHEM\1\METHODS\SIMALC1.M


Sample Info: 15052



#	Compound	Peak Area	RT (min)
1	Ethanol	1151	1.092
2	n-Propanol	2546	1.755



Ethanol 0.101 g/100mL



n-Propanol 0.012 g/100mL

Handwritten signature

Handwritten signature

Washington State Patrol Toxicology Laboratory
 2203 Airport Way S Seattle, WA 98134

Inj. Date: 11/18/2015 8:14:08 AM

Sample Name: 0.20 CTRL

Instrument: HSGC#1


Operator: Christie Mitchell-Mata

Column: DB-ALC1


Location: Vial 8

Method: C:\HPCHEM\1\METHODS\SIMALC1.M


Sample Info: 15052



#	Compound	Peak Area	RT (min)
1	Ethanol	2270	1.091
2	n-Propanol	2532	1.754



Ethanol 0.202 g/100mL



n-Propanol 0.012 g/100mL

Handwritten signature

Handwritten signature

Washington State Patrol Toxicology Laboratory
2203 Airport Way S Seattle, WA 98134

Inj. Date: 11/18/2015 8:17:21 AM

Sample Name: NEG CTRL

Instrument: HSGC#1


Operator: Christie Mitchell-Mata

Column: DB-ALC1


Location: Vial 9

Method: C:\HPCHEM\1\METHODS\SIMALC1.M


Sample Info: 15052



#	Compound	Peak Area	RT (min)
1	Ethanol	0	0.000
2	n-Propanol	2546	1.755



Ethanol 0.000 g/100mL



n-Propanol 0.012 g/100mL

Handwritten signature


Handwritten signature

Washington State Patrol Toxicology Laboratory
 2203 Airport Way S Seattle, WA 98134


Inj. Date: 11/18/2015 8:20:35 AM
 Instrument: HSGC#1
 Column: DB-ALC1
 Method: C:\HPCHEM\1\METHODS\SIMALC1.M

Sample Name: 15052-1
 Operator: Christie Mitchell-Mata
 Location: Vial 10


Sample Info:



#	Compound	Peak Area	RT (min)
1	Ethanol	587	1.093
2	n-Propanol	2616	1.756



Ethanol 0.049 g/100mL



n-Propanol 0.012 g/100mL

Handwritten mark

Handwritten mark

Washington State Patrol Toxicology Laboratory
 2203 Airport Way S Seattle, WA 98134

Inj. Date: 11/18/2015 8:23:48 AM

Sample Name: 15052-2

Instrument: HSGC#1


Operator: Christie Mitchell-Mata

Column: DB-ALC1


Location: Vial 11

Method: C:\HPCHEM\1\METHODS\SIMALC1.M

Sample Info:




#	Compound	Peak Area	RT (min)
1	Ethanol	608	1.093
2	n-Propanol	2696	1.756



Ethanol 0.050 g/100mL

Handwritten signature



n-Propanol 0.012 g/100mL

Handwritten mark

Washington State Patrol Toxicology Laboratory
2203 Airport Way S Seattle, WA 98134

Inj. Date: 11/18/2015 8:27:01 AM

Sample Name: 15052-3

Instrument: HSGC#1


Operator: Christie Mitchell-Mata

Column: DB-ALC1


Location: Vial 12

Method: C:\HPCHEM\1\METHODS\SIMALC1.M


Sample Info:



#	Compound	Peak Area	RT (min)
1	Ethanol	613	1.092
2	n-Propanol	2710	1.756



Ethanol 0.050 g/100mL



n-Propanol 0.012 g/100mL

Handwritten signature

Handwritten signature

Washington State Patrol Toxicology Laboratory
2203 Airport Way S Seattle, WA 98134

Inj. Date: 11/18/2015 8:30:14 AM

Sample Name: 15052-4

Instrument: HSGC#1


Operator: Christie Mitchell-Mata

Column: DB-ALC1


Location: Vial 13

Method: C:\HPCHEM\1\METHODS\SIMALC1.M


Sample Info:



#	Compound	Peak Area	RT (min)
1	Ethanol	571	1.091
2	n-Propanol	2557	1.755



Ethanol 0.049 g/100mL



n-Propanol 0.012 g/100mL

Washington State Patrol Toxicology Laboratory
 2203 Airport Way S Seattle, WA 98134

Inj. Date: 11/18/2015 8:33:28 AM

Sample Name: 15052-5

Instrument: HSGC#1


Operator: Christie Mitchell-Mata

Column: DB-ALC1


Location: Vial 14

Method: C:\HPCHEM\1\METHODS\SIMALC1.M


Sample Info:



#	Compound	Peak Area	RT (min)
1	Ethanol	573	1.091
2	n-Propanol	2564	1.754



Ethanol 0.049 g/100mL



n-Propanol 0.012 g/100mL

Handwritten signature

Handwritten signature


Washington State Patrol Toxicology Laboratory
 2203 Airport Way S Seattle, WA 98134

Inj. Date: 11/18/2015 8:36:41 AM
 Instrument: HSGC#1


Sample Name: 0.10 CTRL
 Operator: Christie Mitchell-Mata
 Location: Vial 15

Column: DB-ALC1
 Method: C:\HPCHEM\1\METHODS\SIMALC1.M


Sample Info: 15052



#	Compound	Peak Area	RT (min)
1	Ethanol	1162	1.093
2	n-Propanol	2571	1.756



Ethanol 0.101 g/100mL



n-Propanol 0.012 g/100mL

h

ur

Washington State Patrol Toxicology Laboratory
2203 Airport Way S Seattle, WA 98134

Inj. Date: 11/18/2015 8:39:54 AM

Sample Name: NEG CTRL

Instrument: HSGC#1


Operator: Christie Mitchell-Mata

Column: DB-ALC1


Location: Vial 16

Method: C:\HPCHEM\1\METHODS\SIMALC1.M

Sample Info: 15052




#	Compound	Peak Area	RT (min)
1	Ethanol	0	0.000
2	n-Propanol	2548	1.755



Ethanol 0.000 g/100mL

Handwritten signature



n-Propanol 0.012 g/100mL

Handwritten signature