



QUALITY ASSURANCE PROCEDURE SOLUTION TEST REPORT

BATCH REPORT: 16005

CUSTOMER INFORMATION

Washington State Patrol – Breath Test Program
811 East Roanoke SEATTLE, WA 98102

TESTING PROCEDURE USED: TLD Technical Manual, Chapter 4.0 Certification of Simulator Solutions;
Headspace-Gas Chromatography.

TESTING ITEM INFORMATION

TARGET VAPOR CONCENTRATION: 0.15 g/210L
DATE PREPARED: 02/10/2016
BATCH UNITS: g/100mL

IDENTITY: QAP Solution
PREPARED BY: Elizabeth Wehner

	EW	JLK	AG
1	0.190	0.191	0.190
2	0.190	0.191	0.191
3	0.191	0.192	0.191
4	0.191	0.192	0.192
5	0.192	0.191	0.190
C	0.104	0.103	0.102

ETHANOL CONTROL INFORMATION

LOT NUMBER: FN08051301 EXPIRATION: 10/2018 CONCENTRATION: 0.10 g/100mL

RESULTS OF TESTING

AVERAGE SOLUTION CONCENTRATION: 0.1910 g/100mL PRECISION CV (%): 0.40
STANDARD DEVIATION: 0.00076 NUMBER OF TESTS: 15

EQUIVALENT VAPOR CONCENTRATION: **0.1553 g/210L**
EXPANDED UNCERTAINTY: ± 0.0038 (k=2, 95.45% confidence interval)

WASHINGTON STATE PATROL – TOXICOLOGY LABORATORY DIVISION






Lisa Noble Forensic Scientist Supervisor

3/18/16

DATE REPORT ISSUED

THIS TESTING WAS PERFORMED BY:

ANALYST	NAME	SIGNATURE	DATE TESTED
EW	Elizabeth Wehner		02/10/2016
JLK	Justin L. Knoy		02/11/2016
AG	Andrew Gingras		02/12/2016

This report applies only to the item being tested and shall not be reproduced except in full, without the written approval of the WSP Toxicology Laboratory Division. Page 1 of 1

Washington State Patrol - Toxicology Laboratory Division
QAP Test Report Calculation Record

QAP Solution Batch #: 16005

Date Prepared: 2/10/2016

Analyst:	EW	JLK	AG
Date Tested:	2/10/2016	2/11/2016	2/12/2016
Instrument:	HSGC #3	HSGC #3	HSGC #3
1	0.190	0.191	0.190
2	0.190	0.191	0.191
3	0.191	0.192	0.191
4	0.191	0.192	0.192
5	0.192	0.191	0.190
C	0.104	0.103	0.102

CV^2_{COA}	$CV^2_{QAP\ Solution}$	$CV^2_{Control}$	$CV^2_{Part\ Coef}$
0.0000084100	0.0000010442	0.0000314199	0.0001016326

Ethanol Control Lot #: FN08051301
Control Uncertainty (%): 0.29

Average Solution Concentration: 0.1910 g/100mL
Standard Deviation: 0.00076 g/100mL
Precision CV (%): 0.40
Equivalent Vapor Concentration: 0.1553 g/210L
Combined Standard Uncertainty (\pm): 0.0019 g/210L
Expanded Uncertainty (\pm): 0.0038 coverage factor (k) =2 (95.45% level of confidence)

Calculations performed by: Lisa Noble [Signature] 3/7/16
Name Signature Date

Calculations verified by: Amanda M. Black [Signature] 3-17-16
Name Signature Date

Method: Hand Calculation

Tech. review performed by: Lisa Noble [Signature] 3/7/16
Name Signature Date

[Handwritten mark]

SIMULATOR SOLUTION DATA ENTRY REVIEW

Reviewer/s: Amanda M. Black Date: 3-17-16

Location: WSP-FLSB Seattle, WA Solution Batch Number: 16005

	YES	NO	N/A
Analysis dates do not precede preparation date:	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Declarations signed and properly dated:	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Data entry corresponds to all chromatograms:	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
All signatures present on Test Report:	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Average solution concentration correct:	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Standard deviation correct:	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
CV (%) correct:	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Equivalent vapor concentration correct:	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
All chromatograms and sequences included in file:	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Ethanol control information present: (lot # present & used within expiration)	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Complies with accuracy and precision requirements established by the State Toxicologist:	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

Comments:

Reviewer Signature: _____



Date: _____

3-17-16



SOLUTION CERTIFICATE REVIEW

Please check that the data on your chromatograms is the data entered into the Test Report, that the date to the right of your name is the date that you tested the solution, and then sign the Test Report.

Please initial and date below to affirm that you have:

- 1) Checked your data
- 2) Checked the date to the right of your name on the Test Report
- 3) Signed the Test Report

	Initials	Date
Amanda Chandler		
Andrew Gingras	<i>AG</i>	3/10/16
Asa Louis		
Brittany Thomas		
Christie Mitchell-Mata		
Christopher Johnston		
David Nguyen		
Dawn Sklerov		
Elizabeth Wehner	<i>EW</i>	03/09/16
Justin Knoy	<i>JK</i>	3.10.16
Katie Harris		
Lyndsey Lowe		
Naziha Nuwayhid		
Rebecca Flaherty		

Batch # 16005 *Ln 3/17/16*

JAY INSLEE
Governor



JOHN R. BATISTE
Chief

STATE OF WASHINGTON
WASHINGTON STATE PATROL
WASHINGTON STATE TOXICOLOGY LABORATORY

2203 Airport Way South, Suite 360 • Seattle, Washington 98134-2027 • (206) 262-6100 • FAX (206) 262-6145

**0.15 g/210 L QUALITY ASSURANCE PROCEDURE SOLUTION
CERTIFICATION FOR LOT 16005**

I, Elizabeth Wehner, do certify under penalty of perjury that:

I am employed by the Washington State Toxicology Laboratory, and my responsibilities include the preparation and certification of alcohol solutions for use with evidential breath test instruments.

I possess the following qualifications: BS degree in Biochemistry.

The quality assurance procedure (QAP) solution, Lot Number 16005, was prepared in the Washington State Toxicology Laboratory on 2/10/2016. I tested this solution and it was found to conform to those standards established by the State Toxicologist for the certification of simulator solution. It shall not be used to perform a quality assurance procedure after 2/10/2017.

Seattle, WA

Elizabeth Wehner 03/09/14

Elizabeth Wehner

Date

Forensic Scientist

JAY INSLEE
Governor



JOHN R. BATISTE
Chief

STATE OF WASHINGTON
WASHINGTON STATE PATROL
WASHINGTON STATE TOXICOLOGY LABORATORY

2203 Airport Way South, Suite 360 • Seattle, Washington 98134-2027 • (206) 262-6100 • FAX (206) 262-6145

**0.15 g/210 L QUALITY ASSURANCE PROCEDURE SOLUTION
CERTIFICATION FOR LOT 16005**


I, Justin L. Knoy, do certify under penalty of perjury that:

I am employed by the Washington State Toxicology Laboratory, and my responsibilities include the preparation and certification of alcohol solutions for use with evidential breath test instruments.

I possess the following qualifications: BS degree in Biology, and MS degree in Forensic Science.

The quality assurance procedure (QAP) solution, Lot Number 16005, was prepared in the Washington State Toxicology Laboratory on 2/10/2016. I tested this solution and it was found to conform to those standards established by the State Toxicologist for the certification of simulator solution. It shall not be used to perform a quality assurance procedure after 2/10/2017.

Seattle, WA

 3.10.16
Justin L. Knoy Date
Forensic Scientist



JAY INSLEE
Governor



JOHN R. BATISTE
Chief

STATE OF WASHINGTON
WASHINGTON STATE PATROL
WASHINGTON STATE TOXICOLOGY LABORATORY

2203 Airport Way South, Suite 360 • Seattle, Washington 98134-2027 • (206) 262-6100 • FAX (206) 262-6145

**0.15 g/210 L QUALITY ASSURANCE PROCEDURE SOLUTION
CERTIFICATION FOR LOT 16005**

I, Andrew Gingras, do certify under penalty of perjury that:

I am employed by the Washington State Toxicology Laboratory, and my responsibilities include the preparation and certification of alcohol solutions for use with evidential breath test instruments.

I possess the following qualifications: BS degree in Cell and Molecular Biology and MS degree in Forensic Science.

The quality assurance procedure (QAP) solution, Lot Number 16005, was prepared in the Washington State Toxicology Laboratory on 2/10/2016. I tested this solution and it was found to conform to those standards established by the State Toxicologist for the certification of simulator solution. It shall not be used to perform a quality assurance procedure after 2/10/2017.

Seattle, WA

 3/10/2016

Andrew Gingras

Date

Forensic Scientist

FILE A COPY IN THE BATCH FILE FOR EACH SOLUTION LISTED ON THE WORKSHEET

Preparation Date: 02/10/16 Expiration Date: 02/10/17 Initials of Preparer: EW

Lot # of 200-proof Ethanol used in preparation: 2EA0437

Date the 200-proof Ethanol bottle was opened: 11/13/15

After opening, each bottle of 200-proof Ethanol is approved for use for 6 months unless an extension is approved by the State Toxicologist. This timeframe applies to the 200-proof Ethanol only, not to simulator solutions which have a 1 year expiration.

Environmental conditions verified as acceptable:

Simulator Solution	Volume of Ethanol (mL)	Volume of Deionized Water (L)		Batch Number
QAP 0.04	11.2	18	<input checked="" type="checkbox"/>	<u>16002</u>
QAP 0.08	22.4	18	<input checked="" type="checkbox"/>	<u>16003</u>
QAP 0.10	28.1	18	<input checked="" type="checkbox"/>	<u>16004</u>
QAP 0.15	42.1	18	<input checked="" type="checkbox"/>	<u>16005</u>
QAP 0.20	56.1	18	<input type="checkbox"/>	
ESS	66.5	52	<input type="checkbox"/>	

Stir bar is rotating

Stirred for minimum 30 minutes; 2 hours for ESS

Spigot purged

Aliquot taken

Batch labeled, packaged and sealed

02/10/16
Date

If different ethanol lot numbers are used in the preparation of solutions, record them and the corresponding solution batch numbers in the comments section.

Comments:

Elizabeth Wehner
Analyst Signature

02/10/16
Date

EW

Sequence Parameters:

Operator: Elizabeth Wehner
 Data File Naming: Prefix/Counter
 Signal 1 Prefix: SIG1
 Counter: 0001
 Signal 2 Prefix: SIG2
 Counter: 0001
 Data Directory: C:\HPCHEM\2\DATA\
 Data Subdirectory: 160210EW
 Part of Methods to run: According to Runtime Checklist
 Barcode Reader: not used
 Shutdown Cmd/Macro: none

Sequence Comment:

CAL 1: 0.079 g/100mL - Lot#: E1015-01 Exp. 04/29/2016
 CAL 2: 0.158 g/100mL - Lot#: E1015-02 Exp. 04/29/2016
 CAL 3: 0.316 g/100mL - Lot#: E1015-03 Exp. 04/29/2016

 CTRL 1: 0.04 g/100mL - Lot#: FN05011301 Exp. 05/2018
 CTRL 2: 0.10 g/100mL - Lot#: FN08051301 Exp. 10/2018
 CTRL 3: 0.20 g/100mL - Lot#: FN03211401 Exp. 06/2019

 n-Propanol ISTD - Lot#: P0216 Exp: 05/02/2016

 Calibration vials 1-9 are filed with Batch 16002.

Sequence Table (Front Injector):

Method and Injection Info Part:

Line	Location	SampleName	Method	Inj	SampleType	InjVolume	DataFile
1	Vial 1	BLANK	SIMALC3	1	Sample		
2	Vial 2	CAL 1 (0.079)	SIMALC3	1	Calib		
3	Vial 3	CAL 2 (0.158)	SIMALC3	1	Calib		
4	Vial 4	CAL 3 (0.316)	SIMALC3	1	Calib		
5	Vial 5	NEG CTRL	SIMALC3	1	Ctrl Samp		
6	Vial 6	CTRL 1 (0.04)	SIMALC3	1	Ctrl Samp		
7	Vial 7	CTRL 2 (0.10)	SIMALC3	1	Ctrl Samp		
8	Vial 8	CTRL 3 (0.20)	SIMALC3	1	Ctrl Samp		
9	Vial 9	NEG CTRL	SIMALC3	1	Ctrl Samp		
10	Vial 10	16002 #1	SIMALC3	1	Sample		
11	Vial 11	16002 #2	SIMALC3	1	Sample		
12	Vial 12	16002 #3	SIMALC3	1	Sample		
13	Vial 13	16002 #4	SIMALC3	1	Sample		
14	Vial 14	16002 #5	SIMALC3	1	Sample		
15	Vial 15	POS CTRL (0.10)	SIMALC3	1	Ctrl Samp		
16	Vial 16	NEG CTRL	SIMALC3	1	Ctrl Samp		
17	Vial 17	16003 #1	SIMALC3	1	Sample		
18	Vial 18	16003 #2	SIMALC3	1	Sample		
19	Vial 19	16003 #3	SIMALC3	1	Sample		
20	Vial 20	16003 #4	SIMALC3	1	Sample		
21	Vial 21	16003 #5	SIMALC3	1	Sample		
22	Vial 22	POS CTRL (0.10)	SIMALC3	1	Ctrl Samp		
23	Vial 23	NEG CTRL	SIMALC3	1	Ctrl Samp		
24	Vial 24	16004 #1	SIMALC3	1	Sample		

16005
 In 2/10/16

EW

Line	Location	SampleName	Method	Inj	SampleType	InjVolume	DataFile
25	Vial 25	16004 #2	SIMALC3	1	Sample		
26	Vial 26	16004 #3	SIMALC3	1	Sample		
27	Vial 27	16004 #4	SIMALC3	1	Sample		
28	Vial 28	16004 #5	SIMALC3	1	Sample		
29	Vial 29	POS CTRL (0.10)	SIMALC3	1	Ctrl Samp		
30	Vial 30	NEG CTRL	SIMALC3	1	Ctrl Samp		
31	Vial 31	16005 #1	SIMALC3	1	Sample		
32	Vial 32	16005 #2	SIMALC3	1	Sample		
33	Vial 33	16005 #3	SIMALC3	1	Sample		
34	Vial 34	16005 #4	SIMALC3	1	Sample		
35	Vial 35	16005 #5	SIMALC3	1	Sample		
36	Vial 36	POS CTRL (0.10)	SIMALC3	1	Ctrl Samp		
37	Vial 37	NEG CTRL	SIMALC3	1	Ctrl Samp		

Calibration Part:

Line	Location	SampleName	Method	CalLev	Update	RF	Update	RT	Interval
2	Vial 2	CAL 1 (0.079)	SIMALC3	1	Replace		Replace		
3	Vial 3	CAL 2 (0.158)	SIMALC3	2	Replace		Replace		
4	Vial 4	CAL 3 (0.316)	SIMALC3	3	Replace		Replace		

Sequence Table (Back Injector):

No entries - empty table!

16005

EW
3/7/16

EW

Washington State Patrol Toxicology Laboratory
2203 Airport Way S Seattle, WA 98134

Inj. Date: 2/10/2016 11:45:15 AM

Sample Name: 16005 #1

Instrument: HSGC#3


Operator: Elizabeth Wehner

Column: DB-ALC2


Location: Vial 31

Method: C:\HPCHEM\2\METHODS\SIMALC3.M


Sample Info:



#	Compound	Peak Area	RT (min)
1	Ethanol	1364	1.023
2	n-Propanol	1619	1.751



Ethanol 0.190 g/100mL



n-Propanol 0.012 g/100mL

Ym

EW

Inj. Date: 2/10/2016 11:48:28 AM

Sample Name: 16005 #2

Instrument: HSGC#3


Operator: Elizabeth Wehner

Column: DB-ALC2


Location: Vial 32

Method: C:\HPCHEM\2\METHODS\SIMALC3.M


Sample Info:



#	Compound	Peak Area	RT (min)
1	Ethanol	1352	1.023
2	n-Propanol	1600	1.751



Ethanol 0.190 g/100mL



n-Propanol 0.012 g/100mL

EW

EW

Inj. Date: 2/10/2016 11:51:42 AM

Sample Name: 16005 #3

Instrument: HSGC#3


Operator: Elizabeth Wehner

Column: DB-ALC2


Location: Vial 33

Method: C:\HPCHEM\2\METHODS\SIMALC3.M


Sample Info:



#	Compound	Peak Area	RT (min)
1	Ethanol	1376	1.023
2	n-Propanol	1622	1.752



Ethanol 0.191 g/100mL



n-Propanol 0.012 g/100mL

EW

EW

Inj. Date: 2/10/2016 11:54:54 AM

Sample Name: 16005 #4

Instrument: HSGC#3


Operator: Elizabeth Wehner

Column: DB-ALC2


Location: Vial 34

Method: C:\HPCHEM\2\METHODS\SIMALC3.M


Sample Info:



#	Compound	Peak Area	RT (min)
1	Ethanol	1371	1.023
2	n-Propanol	1619	1.751



Ethanol 0.191 g/100mL



n-Propanol 0.012 g/100mL

Handwritten signature

EW

Washington State Patrol Toxicology Laboratory
 2203 Airport Way S Seattle, WA 98134

Inj. Date: 2/10/2016 11:58:08 AM

Sample Name: 16005 #5

Instrument: HSGC#3


Operator: Elizabeth Wehner

Column: DB-ALC2


Location: Vial 35

Method: C:\HPCHEM\2\METHODS\SIMALC3.M


Sample Info:



#	Compound	Peak Area	RT (min)
1	Ethanol	1386	1.023
2	n-Propanol	1626	1.751



Ethanol 0.192 g/100mL



n-Propanol 0.012 g/100mL

EW

EW

Washington State Patrol Toxicology Laboratory
 2203 Airport Way S Seattle, WA 98134

Inj. Date: 2/10/2016 12:01:22 PM

Sample Name: POS CTRL (0.10)

Instrument: HSGC#3

Operator: Elizabeth Wehner


Column: DB-ALC2

Location: Vial 36


Method: C:\HPCHEM\2\METHODS\SIMALC3.M

Sample Info: POS CTRL: 0.10 g/100mL
 16005


->



#	Compound	Peak Area	RT (min)
1	Ethanol	780	1.023
2	n-Propanol	1690	1.751



Ethanol 0.104 g/100mL



n-Propanol 0.012 g/100mL

Handwritten signature

EW

Inj. Date: 2/10/2016 12:04:36 PM

Sample Name: NEG CTRL

Instrument: HSGC#3


Operator: Elizabeth Wehner

Column: DB-ALC2


Location: Vial 37

Method: C:\HPCHEM\2\METHODS\SIMALC3.M


Sample Info: 16005



#	Compound	Peak Area	RT (min)
1	Ethanol	0	0.000
2	n-Propanol	1703	1.752



Ethanol 0.000 g/100mL



n-Propanol 0.012 g/100mL

Handwritten initials

EW

Sequence Parameters:

Operator: Justin Knoy
 Data File Naming: Prefix/Counter
 Signal 1 Prefix: SIG1
 Counter: 0001
 Signal 2 Prefix: SIG2
 Counter: 0001
 Data Directory: C:\HPCHEM\2\DATA\
 Data Subdirectory: 160211JK
 Part of Methods to run: According to Runtime Checklist
 Barcode Reader: not used
 Shutdown Cmd/Macro: none

Sequence Comment:

Ethanol Calibrator 1, E1015-01 - Exp. 04/29/2016
 Ethanol Calibrator 2, E1015-02 - Exp. 04/29/2016
 Ethanol Calibrator 3, E1015-03 - Exp. 04/29/2016
 CTRL1 (0.04g/100mL), Lot # FN05011301 - Exp. 05/2018
 CTRL2 (0.10g/100mL), Lot # FN08051301 - Exp. 10/2018
 CTRL3 (0.20g/100mL), Lot # FN03211401 - Exp. 06/2019

 Internal Standard Lot#P0216 - Exp. 05/02/2016

 Calibration vials 1-9 filed with 16002.

Sequence Table (Front Injector):

Method and Injection Info Part:

Line	Location	SampleName	Method	Inj	SampleType	InjVolume	DataFile
1	Vial 1	BLANK	SIMALC3	1	Sample		
2	Vial 2	CAL1 0.079	SIMALC3	1	Calib		
3	Vial 3	CAL2 0.158	SIMALC3	1	Calib		
4	Vial 4	CAL3 0.316	SIMALC3	1	Calib		
5	Vial 5	NEG CTRL	SIMALC3	1	Ctrl Samp		
6	Vial 6	CTRL1 (0.04)	SIMALC3	1	Ctrl Samp		
7	Vial 7	CTRL2 (0.10)	SIMALC3	1	Ctrl Samp		
8	Vial 8	CTRL3 (0.20)	SIMALC3	1	Ctrl Samp		
9	Vial 9	NEG CTRL	SIMALC3	1	Ctrl Samp		
10	Vial 10	16002-1	SIMALC3	1	Sample		
11	Vial 11	16002-2	SIMALC3	1	Sample		
12	Vial 12	16002-3	SIMALC3	1	Sample		
13	Vial 13	16002-4	SIMALC3	1	Sample		
14	Vial 14	16002-5	SIMALC3	1	Sample		
15	Vial 15	CTRL2 (0.10)	SIMALC3	1	Ctrl Samp		
16	Vial 16	NEG CTRL	SIMALC3	1	Ctrl Samp		
17	Vial 17	16003-1	SIMALC3	1	Sample		
18	Vial 18	16003-2	SIMALC3	1	Sample		
19	Vial 19	16003-3	SIMALC3	1	Sample		
20	Vial 20	16003-4	SIMALC3	1	Sample		
21	Vial 21	16003-5	SIMALC3	1	Sample		
22	Vial 22	CTRL2 (0.10)	SIMALC3	1	Ctrl Samp		
23	Vial 23	NEG CTRL	SIMALC3	1	Ctrl Samp		
24	Vial 24	16004-1	SIMALC3	1	Sample		
25	Vial 25	16004-2	SIMALC3	1	Sample		

16005
Jn 3/7/16

JY

Line	Location	SampleName	Method	Inj	SampleType	InjVolume	DataFile
26	Vial 26	16004-3	SIMALC3	1	Sample		
27	Vial 27	16004-4	SIMALC3	1	Sample		
28	Vial 28	16004-5	SIMALC3	1	Sample		
29	Vial 29	CTRL2 (0.10)	SIMALC3	1	Ctrl Samp		
30	Vial 30	NEG CTRL	SIMALC3	1	Ctrl Samp		
31	Vial 31	16005-1	SIMALC3	1	Sample		
32	Vial 32	16005-2	SIMALC3	1	Sample		
33	Vial 33	16005-3	SIMALC3	1	Sample		
34	Vial 34	16005-4	SIMALC3	1	Sample		
35	Vial 35	16005-5	SIMALC3	1	Sample		
36	Vial 36	CTRL2 (0.10)	SIMALC3	1	Ctrl Samp		
37	Vial 37	NEG CTRL	SIMALC3	1	Ctrl Samp		

Calibration Part:

Line	Location	SampleName	Method	CalLev	Update	RF	Update	RT	Interval
2	Vial 2	CAL1 0.079	SIMALC3	1	Replace		Replace		
3	Vial 3	CAL2 0.158	SIMALC3	2	Replace		Replace		
4	Vial 4	CAL3 0.316	SIMALC3	3	Replace		Replace		

Sequence Table (Back Injector):

No entries - empty table!

16005

JK

JK

Inj. Date: 2/11/2016 5:43:48 PM

Sample Name: 16005-1

Instrument: HSGC#3


Operator: Justin Knoy

Column: DB-ALC2


Location: Vial 31

Method: C:\HPCHEM\2\METHODS\SIMALC3.M


Sample Info:



#	Compound	Peak Area	RT (min)
1	Ethanol	1365	1.022
2	n-Propanol	1603	1.749



Ethanol 0.191 g/100mL



n-Propanol 0.012 g/100mL

JK

JK

Inj. Date: 2/11/2016 5:47:01 PM

Sample Name: 16005-2

Instrument: HSGC#3


Operator: Justin Knoy

Column: DB-ALC2


Location: Vial 32

Method: C:\HPCHEM\2\METHODS\SIMALC3.M


Sample Info:



#	Compound	Peak Area	RT (min)
1	Ethanol	1379	1.022
2	n-Propanol	1611	1.751



Ethanol 0.191 g/100mL



n-Propanol 0.012 g/100mL

fn

JX

Washington State Patrol Toxicology Laboratory
 2203 Airport Way S Seattle, WA 98134

Inj. Date: 2/11/2016 5:50:14 PM

Sample Name: 16005-3

Instrument: HSGC#3


Operator: Justin Knoy

Column: DB-ALC2


Location: Vial 33

Method: C:\HPCHEM\2\METHODS\SIMALC3.M


Sample Info:



#	Compound	Peak Area	RT (min)
1	Ethanol	1379	1.023
2	n-Propanol	1606	1.751



Ethanol 0.192 g/100mL



n-Propanol 0.012 g/100mL

JK

JK

Washington State Patrol Toxicology Laboratory
 2203 Airport Way S Seattle, WA 98134

Inj. Date: 2/11/2016 5:53:27 PM

Sample Name: 16005-4

Instrument: HSGC#3


Operator: Justin Knoy

Column: DB-ALC2


Location: Vial 34

Method: C:\HPCHEM\2\METHODS\SIMALC3.M


Sample Info:



#	Compound	Peak Area	RT (min)
1	Ethanol	1385	1.023
2	n-Propanol	1613	1.752



Ethanol 0.192 g/100mL



n-Propanol 0.012 g/100mL

JK

JK

Inj. Date: 2/11/2016 5:56:41 PM

Sample Name: 16005-5

Instrument: HSGC#3


Operator: Justin Knoy

Column: DB-ALC2


Location: Vial 35

Method: C:\HPCHEM\2\METHODS\SIMALC3.M


Sample Info:



#	Compound	Peak Area	RT (min)
1	Ethanol	1368	1.023
2	n-Propanol	1601	1.751



Ethanol 0.191 g/100mL




n-Propanol 0.012 g/100mL


JN

JZ


Inj. Date: 2/11/2016 5:59:55 PM Sample Name: CTRL2 (0.10)
Instrument: HSGC#3 Operator: Justin Knoy
Column: DB-ALC2 Location: Vial 36
Method: C:\HPCHEM\2\METHODS\SIMALC3.M
Sample Info: 0.10g/100mL ; 16005



#	Compound	Peak Area	RT (min)
1	Ethanol	754	1.023
2	n-Propanol	1636	1.752



Ethanol 0.103 g/100mL




n-Propanol 0.012 g/100mL

JK


JK

Washington State Patrol Toxicology Laboratory
2203 Airport Way S Seattle, WA 98134


Inj. Date: 2/11/2016 6:03:09 PM Sample Name: NEG CTRL
Instrument: HSGC#3 Operator: Justin Knoy
Column: DB-ALC2 Location: Vial 37
Method: C:\HPCHEM\2\METHODS\SIMALC3.M
Sample Info: 16005



#	Compound	Peak Area	RT (min)
1	Ethanol	0	0.000
2	n-Propanol	1667	1.752



Ethanol 0.000 g/100mL



n-Propanol 0.012 g/100mL

Handwritten signature

Handwritten initials

Sequence Parameters:

Sequence run by:

Operator: Andrew Gingras
 Data File Naming: Prefix/Counter
 Signal 1 Prefix: SIG1
 Counter: 0001
 Signal 2 Prefix: SIG2
 Counter: 0001
 Data Directory: C:\HPCHEM\2\DATA\
 Data Subdirectory: 160212AG
 Part of Methods to run: According to Runtime Checklist
 Barcode Reader: not used
 Shutdown Cmd/Macro: none

AG
3/10/16

Sequence Comment:

Ethanol Calibrator 1, E1015-01 - Exp. 04/29/2016
 Ethanol Calibrator 2, E1015-02 - Exp. 04/29/2016
 Ethanol Calibrator 3, E1015-03 - Exp. 04/29/2016
 CTRL1 (0.04g/100mL), Lot # FN05011301 - Exp. 05/2018
 CTRL2 (0.10g/100mL), Lot # FN08051301 - Exp. 10/2018
 CTRL3 (0.20g/100mL), Lot # FN03211401 - Exp. 06/2019

 Internal Standard Lot#P0216 - Exp. 05/02/2016

 Calibration vials 1-9 filed with 16002.

Sequence Table (Front Injector):

Method and Injection Info Part:

Line	Location	SampleName	Method	Inj	SampleType	InjVolume	DataFile
1	Vial 1	BLANK	SIMALC3	1	Sample		
2	Vial 2	CAL1 0.079	SIMALC3	1	Calib		
3	Vial 3	CAL2 0.158	SIMALC3	1	Calib		
4	Vial 4	CAL3 0.316	SIMALC3	1	Calib		
5	Vial 5	NEG CTRL	SIMALC3	1	Ctrl Samp		
6	Vial 6	CTRL1 (0.04)	SIMALC3	1	Ctrl Samp		
7	Vial 7	CTRL2 (0.10)	SIMALC3	1	Ctrl Samp		
8	Vial 8	CTRL3 (0.20)	SIMALC3	1	Ctrl Samp		
9	Vial 9	NEG CTRL	SIMALC3	1	Ctrl Samp		
10	Vial 10	16002-1	SIMALC3	1	Sample		
11	Vial 11	16002-2	SIMALC3	1	Sample		
12	Vial 12	16002-3	SIMALC3	1	Sample		
13	Vial 13	16002-4	SIMALC3	1	Sample		
14	Vial 14	16002-5	SIMALC3	1	Sample		
15	Vial 15	CTRL2 (0.10)	SIMALC3	1	Ctrl Samp		
16	Vial 16	NEG CTRL	SIMALC3	1	Ctrl Samp		
17	Vial 17	16003-1	SIMALC3	1	Sample		
18	Vial 18	16003-2	SIMALC3	1	Sample		
19	Vial 19	16003-3	SIMALC3	1	Sample		
20	Vial 20	16003-4	SIMALC3	1	Sample		
21	Vial 21	16003-5	SIMALC3	1	Sample		
22	Vial 22	CTRL2 (0.10)	SIMALC3	1	Ctrl Samp		
23	Vial 23	NEG CTRL	SIMALC3	1	Ctrl Samp		
24	Vial 24	16004-1	SIMALC3	1	Sample		
25	Vial 25	16004-2	SIMALC3	1	Sample		

16005
fn 3/17/16

AG

Line	Location	SampleName	Method	Inj	SampleType	InjVolume	DataFile
26	Vial 26	16004-3	SIMALC3	1	Sample		
27	Vial 27	16004-4	SIMALC3	1	Sample		
28	Vial 28	16004-5	SIMALC3	1	Sample		
29	Vial 29	CTRL2 (0.10)	SIMALC3	1	Ctrl Samp		
30	Vial 30	NEG CTRL	SIMALC3	1	Ctrl Samp		
31	Vial 31	16005-1	SIMALC3	1	Sample		
32	Vial 32	16005-2	SIMALC3	1	Sample		
33	Vial 33	16005-3	SIMALC3	1	Sample		
34	Vial 34	16005-4	SIMALC3	1	Sample		
35	Vial 35	16005-5	SIMALC3	1	Sample		
36	Vial 36	CTRL2 (0.10)	SIMALC3	1	Ctrl Samp		
37	Vial 37	NEG CTRL	SIMALC3	1	Ctrl Samp		

Calibration Part:

Line	Location	SampleName	Method	CalLev	Update	RF	Update	RT	Interval
2	Vial 2	CAL1 0.079	SIMALC3	1	Replace		Replace		
3	Vial 3	CAL2 0.158	SIMALC3	2	Replace		Replace		
4	Vial 4	CAL3 0.316	SIMALC3	3	Replace		Replace		

Sequence Table (Back Injector):

No entries - empty table!

16005

Handwritten signature

Handwritten signature

Washington State Patrol Toxicology Laboratory
2203 Airport Way S Seattle, WA 98134

Inj. Date: 2/12/2016 9:53:33 AM

Sample Name: 16005-1

Instrument: HSGC#3


Operator: Andrew Gingras

Column: DB-ALC2


Location: Vial 31

Method: C:\HPCHEM\2\METHODS\SIMALC3.M


Sample Info:



#	Compound	Peak Area	RT (min)
1	Ethanol	1379	1.022
2	n-Propanol	1625	1.749



Ethanol 0.190 g/100mL



n-Propanol 0.012 g/100mL

h

AC

Inj. Date: 2/12/2016 9:56:45 AM

Sample Name: 16005-2

Instrument: HSGC#3


Operator: Andrew Gingras

Column: DB-ALC2


Location: Vial 32

Method: C:\HPCHEM\2\METHODS\SIMALC3.M


Sample Info:



#	Compound	Peak Area	RT (min)
1	Ethanol	1402	1.022
2	n-Propanol	1643	1.750



Ethanol 0.191 g/100mL



n-Propanol 0.012 g/100mL

Handwritten mark

Handwritten signature

Inj. Date: 2/12/2016 9:59:59 AM

Sample Name: 16005-3

Instrument: HSGC#3


Operator: Andrew Gingras

Column: DB-ALC2


Location: Vial 33

Method: C:\HPCHEM\2\METHODS\SIMALC3.M


Sample Info:



#	Compound	Peak Area	RT (min)
1	Ethanol	1388	1.022
2	n-Propanol	1633	1.749



Ethanol 0.191 g/100mL



n-Propanol 0.012 g/100mL

Handwritten initials

Handwritten signature

Washington State Patrol Toxicology Laboratory
 2203 Airport Way S Seattle, WA 98134

Inj. Date: 2/12/2016 10:03:12 AM

Sample Name: 16005-4

Instrument: HSGC#3


Operator: Andrew Gingras

Column: DB-ALC2


Location: Vial 34

Method: C:\HPCHEM\2\METHODS\SIMALC3.M


Sample Info:



#	Compound	Peak Area	RT (min)
1	Ethanol	1396	1.022
2	n-Propanol	1634	1.750



Ethanol 0.192 g/100mL



n-Propanol 0.012 g/100mL

Handwritten signature

Handwritten signature

Washington State Patrol Toxicology Laboratory
 2203 Airport Way S Seattle, WA 98134

Inj. Date: 2/12/2016 10:06:26 AM

Sample Name: 16005-5

Instrument: HSGC#3


Operator: Andrew Gingras

Column: DB-ALC2


Location: Vial 35

Method: C:\HPCHEM\2\METHODS\SIMALC3.M


Sample Info:



#	Compound	Peak Area	RT (min)
1	Ethanol	1447	1.022
2	n-Propanol	1710	1.750



Ethanol 0.190 g/100mL



n-Propanol 0.012 g/100mL

hr

hr

Washington State Patrol Toxicology Laboratory
 2203 Airport Way S Seattle, WA 98134

Inj. Date: 2/12/2016 10:09:39 AM

Sample Name: CTRL2 (0.10)

Instrument: HSGC#3


Operator: Andrew Gingras

Column: DB-ALC2


Location: Vial 36

Method: C:\HPCHEM\2\METHODS\SIMALC3.M


Sample Info: 0.10g/100mL ; 16005



#	Compound	Peak Area	RT (min)
1	Ethanol	753	1.022
2	n-Propanol	1658	1.748



Ethanol 0.102 g/100mL



n-Propanol 0.012 g/100mL

Handwritten signature

Handwritten signature

Washington State Patrol Toxicology Laboratory
 2203 Airport Way S Seattle, WA 98134

Inj. Date: 2/12/2016 10:12:54 AM

Sample Name: NEG CTRL

Instrument: HSGC#3


Operator: Andrew Gingras

Column: DB-ALC2


Location: Vial 37

Method: C:\HPCHEM\2\METHODS\SIMALC3.M


Sample Info: 16005



#	Compound	Peak Area	RT (min)
1	Ethanol	0	0.000
2	n-Propanol	1665	1.749



Ethanol 0.000 g/100mL



n-Propanol 0.012 g/100mL

ln

SC