



QUALITY ASSURANCE PROCEDURE SOLUTION TEST REPORT

BATCH REPORT: 17036

CUSTOMER INFORMATION

Washington State Patrol – Breath Test Program
811 East Roanoke SEATTLE, WA 98102

TESTING PROCEDURE USED: TLD Technical Manual, Chapter 4.0 Certification of Simulator Solutions;
Headspace-Gas Chromatography.

TESTING ITEM INFORMATION

TARGET VAPOR CONCENTRATION: 0.10 g/210L
DATE PREPARED: 04/12/2017
BATCH UNITS: g/100mL

IDENTITY: QAP Solution
PREPARED BY: Justin L. Knoy

	JLK	AG	CM
1	0.125	0.125	0.127
2	0.125	0.125	0.123
3	0.126	0.125	0.122
4	0.125	0.125	0.123
5	0.125	0.125	0.126
C	0.101	0.102	0.099

ETHANOL CONTROL INFORMATION

LOT NUMBER: FN08051301 EXPIRATION: 10/2018 CONCENTRATION: 0.10 g/100mL

RESULTS OF TESTING

AVERAGE SOLUTION CONCENTRATION: 0.1248 g/100mL PRECISION CV (%): 1.01
STANDARD DEVIATION: 0.00126 NUMBER OF TESTS: 15

EQUIVALENT VAPOR CONCENTRATION: **0.1015 g/210L**
EXPANDED UNCERTAINTY: ± 0.0028 (k=2, 95.45% confidence interval)

WASHINGTON STATE PATROL – TOXICOLOGY LABORATORY DIVISION

Brittany Thomas
Brittany Thomas Forensic Scientist Supervisor

4/20/17
DATE REPORT ISSUED

THIS TESTING WAS PERFORMED BY:

ANALYST	NAME	SIGNATURE	DATE TESTED
JLK	Justin L. Knoy	<i>Justin L. Knoy</i>	04/12/2017
AG	Andrew Gingras	<i>Andrew Gingras</i>	04/13/2017
CM	Christi Mitchell-Mata	<i>Christi Mitchell-Mata</i>	04/17/2017

SIMULATOR SOLUTION DATA ENTRY REVIEW

Reviewer/s: Amanda M. Black Date: 5-2-17

Location: WSP-FLSB Seattle, WA Solution Batch Number: 17036

	YES	NO	N/A
Analysis dates do not precede preparation date:	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Declarations signed and properly dated:	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Data entry corresponds to all chromatograms:	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
All signatures present on Test Report:	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Average solution concentration correct:	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Standard deviation correct:	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
CV (%) correct:	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Equivalent vapor concentration correct:	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
All chromatograms and sequences included in file:	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Ethanol control information present: (lot # present & used within expiration)	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Complies with accuracy and precision requirements established by the State Toxicologist:	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

Comments:

Reviewer Signature: _____



Date: 5-2-17

Washington State Patrol - Toxicology Laboratory Division
QAP Test Report Calculation Record

QAP Solution Batch #: 17036

Date Prepared: 4/12/2017

Analyst:	JLK	AG	CM
Date Tested:	4/12/2017	4/13/2017	4/17/2017
Instrument:	HSGC #1	HSGC #1	HSGC #1
1	0.125	0.125	0.127
2	0.125	0.125	0.123
3	0.126	0.125	0.122
4	0.125	0.125	0.123
5	0.125	0.125	0.126
C	0.101	0.102	0.099

CV^2_{COA}	$CV^2_{QAP\ Solution}$	$CV^2_{Control}$	$CV^2_{Part\ Coef}$
0.0000084100	0.0000068486	0.0000767510	0.0001016326

Ethanol Control Lot #: FN08051301
Control Uncertainty (%): 0.29

Average Solution Concentration: 0.1248 g/100mL
Standard Deviation: 0.00126 g/100mL
Precision CV (%): 1.01
Equivalent Vapor Concentration: 0.1015 g/210L
Combined Standard Uncertainty (\pm): 0.0014 g/210L
Expanded Uncertainty (\pm): 0.0028 coverage factor (k) =2 (95.45% level of confidence)

Calculations performed by: Brittany Thomas Brittany Thomas 4/12/17
Name Signature Date

Calculations verified by: Amanda M. Black [Signature] 5-2-17 Method: Hand calculation
Name Signature Date

Tech. review performed by: Brittany Thomas Brittany Thomas 4/19/17
Name Signature Date

SOLUTION CERTIFICATE REVIEW

Please check that the data on your chromatograms is the data entered into the Test Report, that the date to the right of your name is the date that you tested the solution, and then sign the Test Report.

Please initial and date below to affirm that you have:

- 1) Checked your data
- 2) Checked the date to the right of your name on the Test Report
- 3) Signed the Test Report

	Initials	Date
Amanda Chandler		
Andrew Gingras	<i>AG</i>	<i>4/20/17</i>
Asa Louis		
Brittany Thomas		
Christie Mitchell-Mata	<i>CM</i>	<i>4/20/17</i>
Christopher Johnston		
David Nguyen		
Dawn Sklerov		
Elizabeth Wehner		
Justin Knoy	<i>JK</i>	<i>4.20.17</i>
Katie Harris		
Lyndsey Knoy		
Naziha Nuwayhid		
Rebecca Flaherty		

Batch # 17036

BT 4/20/17

JAY INSLEE
Governor



JOHN R. BATISTE
Chief

STATE OF WASHINGTON
WASHINGTON STATE PATROL
WASHINGTON STATE TOXICOLOGY LABORATORY

2203 Airport Way South, Suite 360 • Seattle, Washington 98134-2027 • (206) 262-6100 • FAX (206) 262-6145

**0.10 g/210 L QUALITY ASSURANCE PROCEDURE SOLUTION
CERTIFICATION FOR LOT 17036**


I, Justin L. Knoy, do certify under penalty of perjury that:

I am employed by the Washington State Toxicology Laboratory, and my responsibilities include the preparation and certification of alcohol solutions for use with evidential breath test instruments.

I possess the following qualifications: BS degree in Biology, MS degree in Forensic Science, and am certified as a Diplomate in Forensic Toxicology by the American Board of Forensic Toxicology.

The quality assurance procedure (QAP) solution, Lot Number 17036, was prepared in the Washington State Toxicology Laboratory on 4/12/2017. I tested this solution and it was found to conform to those standards established by the State Toxicologist for the certification of simulator solution. It shall not be used to perform a quality assurance procedure after 4/12/2018.

Seattle, WA

 4.20.17
Justin L. Knoy Date
Forensic Toxicologist

JAY INSLEE
Governor



JOHN R. BATISTE
Chief

STATE OF WASHINGTON
WASHINGTON STATE PATROL
WASHINGTON STATE TOXICOLOGY LABORATORY

2203 Airport Way South, Suite 360 • Seattle, Washington 98134-2027 • (206) 262-6100 • FAX (206) 262-6145

**0.10 g/210 L QUALITY ASSURANCE PROCEDURE SOLUTION
CERTIFICATION FOR LOT 17036**

I, Christie Mitchell-Mata, do certify under penalty of perjury that:

I am employed by the Washington State Toxicology Laboratory, and my responsibilities include the preparation and certification of alcohol solutions for use with evidential breath test instruments.

I possess the following qualifications: BA degree in Chemistry, MFS degree in Forensic Science.

The quality assurance procedure (QAP) solution, Lot Number 17036, was prepared in the Washington State Toxicology Laboratory on 4/12/2017. I tested this solution and it was found to conform to those standards established by the State Toxicologist for the certification of simulator solution. It shall not be used to perform a quality assurance procedure after 4/12/2018.

Seattle, WA

Handwritten signature of Christie Mitchell-Mata, dated 4/20/2017.

Christie Mitchell-Mata

Date

Forensic Toxicologist



JAY INSLEE
Governor



JOHN R. BATISTE
Chief

STATE OF WASHINGTON
WASHINGTON STATE PATROL
WASHINGTON STATE TOXICOLOGY LABORATORY

2203 Airport Way South, Suite 360 • Seattle, Washington 98134-2027 • (206) 262-6100 • FAX (206) 262-6145

**0.10 g/210 L QUALITY ASSURANCE PROCEDURE SOLUTION
CERTIFICATION FOR LOT 17036**



I, Andrew Gingras, do certify under penalty of perjury that:

I am employed by the Washington State Toxicology Laboratory, and my responsibilities include the preparation and certification of alcohol solutions for use with evidential breath test instruments.

I possess the following qualifications: BS degree in Cell and Molecular Biology and MS degree in Forensic Science.

The quality assurance procedure (QAP) solution, Lot Number 17036, was prepared in the Washington State Toxicology Laboratory on 4/12/2017. I tested this solution and it was found to conform to those standards established by the State Toxicologist for the certification of simulator solution. It shall not be used to perform a quality assurance procedure after 4/12/2018.

Seattle, WA

Andrew Gingras

Date

Forensic Scientist



FILE A COPY IN THE BATCH FILE FOR EACH SOLUTION LISTED ON THE WORKSHEET

Preparation Date: 4-12-17 Expiration Date: 4-12-18 Initials of Preparer: JK

Lot # of 200-proof Ethanol used in preparation: 2FE0139

Date the 200-proof Ethanol bottle was opened: 4-12-17

After opening, each bottle of 200-proof Ethanol is approved for use for 6 months unless an extension is approved by the State Toxicologist. This timeframe applies to the 200-proof Ethanol only, not to simulator solutions which have a 1 year expiration.

Environmental conditions verified as acceptable:

Simulator Solution	Volume of Ethanol (mL)	Volume of Deionized Water (L)		Batch Number
QAP 0.04	11.2	18	<input checked="" type="checkbox"/>	<u>17034</u>
QAP 0.08	22.4	18	<input checked="" type="checkbox"/>	<u>17035</u>
QAP 0.10	28.1	18	<input checked="" type="checkbox"/>	<u>17036</u>
QAP 0.15	42.1	18	<input checked="" type="checkbox"/>	<u>17037</u>
QAP 0.20 ^{0.08}	56.1 ^{22.4 mL}	18	<input checked="" type="checkbox"/>	<u>17038</u>
ESS	66.5	52	<input type="checkbox"/>	<u>N/A</u>

Stir bar is rotating

Stirred for minimum 30 minutes; 2 hours for ESS

Spigot purged

Aliquot taken

Batch labeled, packaged and sealed Date 4-12-17

If different ethanol lot numbers are used in the preparation of solutions, record them and the corresponding solution batch numbers in the comments section.

Comments:

Justin S. King
Analyst Signature

4-12-17
Date

^{PT 4/20/17}
~~17035~~
17036
^{PT}
4/20/17

Sequence Parameters:

Operator: Justin Knoy
 Data File Naming: Prefix/Counter
 Signal 1 Prefix: SIG1
 Counter: 0001
 Signal 2 Prefix: SIG2
 Counter: 0001
 Data Directory: C:\HPCHEM\1\DATA\
 Data Subdirectory: 170412JK
 Part of Methods to run: According to Runtime Checklist
 Barcode Reader: not used
 Shutdown Cmd/Macro: none

Sequence Comment:

Ethanol Calibrator 1, E0217-01 - Exp. 08/21/2017
 Ethanol Calibrator 2, E0217-02 - Exp. 08/21/2017
 Ethanol Calibrator 3, E0217-03 - Exp. 08/21/2017
 CTRL1 (0.04g/100mL), Lot # FN12181501 - Exp. 12/2020
 CTRL2 (0.10g/100mL), Lot # FN08051301 - Exp. 10/2018
 CTRL3 (0.20g/100mL), Lot # FN08101505 - Exp. 02/2021
 Internal Standard Lot#P0317 - Exp. 06/13/2017

Calibration vials 1-9 filed with 17034.

Diluter #1 JK 4.20.17 ✓ Reviewed Pat 4/20/17

Sequence Table (Front Injector):

Method and Injection Info Part:

Line	Location	SampleName	Method	Inj	SampleType	InjVolume	DataFile
1	Vial 1	BLANK	SIMALC1	1	Sample		
2	Vial 2	0.079 CAL 1	SIMALC1	1	Calib		
3	Vial 3	0.158 CAL 2	SIMALC1	1	Calib		
4	Vial 4	0.316 CAL 3	SIMALC1	1	Calib		
5	Vial 5	NEG CTRL	SIMALC1	1	Ctrl Samp		
6	Vial 6	0.04 CTRL	SIMALC1	1	Ctrl Samp		
7	Vial 7	0.10 CTRL	SIMALC1	1	Ctrl Samp		
8	Vial 8	0.20 CTRL	SIMALC1	1	Ctrl Samp		
9	Vial 9	NEG CTRL	SIMALC1	1	Ctrl Samp		
10	Vial 10	17034-1	SIMALC1	1	Sample		
11	Vial 11	17034-2	SIMALC1	1	Sample		
12	Vial 12	17034-3	SIMALC1	1	Sample		
13	Vial 13	17034-4	SIMALC1	1	Sample		
14	Vial 14	17034-5	SIMALC1	1	Sample		
15	Vial 15	0.10 CTRL	SIMALC1	1	Ctrl Samp		
16	Vial 16	NEG CTRL	SIMALC1	1	Ctrl Samp		
17	Vial 17	17035-1	SIMALC1	1	Sample		
18	Vial 18	17035-2	SIMALC1	1	Sample		
19	Vial 19	17035-3	SIMALC1	1	Sample		
20	Vial 20	17035-4	SIMALC1	1	Sample		
21	Vial 21	17035-5	SIMALC1	1	Sample		
22	Vial 22	0.10 CTRL	SIMALC1	1	Ctrl Samp		
23	Vial 23	NEG CTRL	SIMALC1	1	Ctrl Samp		
24	Vial 24	17036-1	SIMALC1	1	Sample		
25	Vial 25	17036-2	SIMALC1	1	Sample		
26	Vial 26	17036-3	SIMALC1	1	Sample		

*17036
Pat
4/19/17*

JK

Line	Location	SampleName	Method	Inj	SampleType	InjVolume	DataFile
27	Vial 27	17036-4	SIMALC1	1	Sample		
28	Vial 28	17036-5	SIMALC1	1	Sample		
29	Vial 29	0.10 CTRL	SIMALC1	1	Ctrl Samp		
30	Vial 30	NEG CTRL	SIMALC1	1	Ctrl Samp		
31	Vial 31	17037-1	SIMALC1	1	Sample		
32	Vial 32	17037-2	SIMALC1	1	Sample		
33	Vial 33	17037-3	SIMALC1	1	Sample		
34	Vial 34	17037-4	SIMALC1	1	Sample		
35	Vial 35	17037-5	SIMALC1	1	Sample		
36	Vial 36	0.10 CTRL	SIMALC1	1	Ctrl Samp		
37	Vial 37	NEG CTRL	SIMALC1	1	Ctrl Samp		
38	Vial 38	17038-1	SIMALC1	1	Sample		
39	Vial 39	17038-2	SIMALC1	1	Sample		
40	Vial 40	17038-3	SIMALC1	1	Sample		
41	Vial 41	17038-4	SIMALC1	1	Sample		
42	Vial 42	17038-5	SIMALC1	1	Sample		
43	Vial 43	0.10 CTRL	SIMALC1	1	Ctrl Samp		
44	Vial 44	NEG CTRL	SIMALC1	1	Ctrl Samp		

Calibration Part:

Line	Location	SampleName	Method	CalLev	Update	RF	Update	RT	Interval
2	Vial 2	0.079 CAL 1	SIMALC1	1	Replace		Replace		
3	Vial 3	0.158 CAL 2	SIMALC1	2	Replace		Replace		
4	Vial 4	0.316 CAL 3	SIMALC1	3	Replace		Replace		

Sequence Table (Back Injector):

No entries - empty table!

17036
RT
4/19/17

JK

JK

Washington State Patrol Toxicology Laboratory
 2203 Airport Way S Seattle, WA 98134

Inj. Date: 4/12/2017 12:19:21 PM

Sample Name: 17036-1

Instrument: HSGC#1


Operator: Justin Knoy

Column: DB-ALC1


Location: Vial 24

Method: C:\HPCHEM\1\METHODS\SIMALC1.M

Sample Info:




#	Compound	Peak Area	RT (min)
1	Ethanol	1491	1.084
2	n-Propanol	2636	1.764



Ethanol 0.125 g/100mL

JK



n-Propanol 0.012 g/100mL

JK

Washington State Patrol Toxicology Laboratory
 2203 Airport Way S Seattle, WA 98134

Inj. Date: 4/12/2017 12:22:34 PM

Sample Name: 17036-2

Instrument: HSGC#1


Operator: Justin Knoy

Column: DB-ALC1


Location: Vial 25

Method: C:\HPCHEM\1\METHODS\SIMALC1.M

Sample Info:




#	Compound	Peak Area	RT (min)
1	Ethanol	1513	1.083
2	n-Propanol	2671	1.763



Ethanol 0.125 g/100mL

BT



n-Propanol 0.012 g/100mL

JT

Washington State Patrol Toxicology Laboratory
2203 Airport Way S Seattle, WA 98134

Inj. Date: 4/12/2017 12:25:47 PM

Sample Name: 17036-3

Instrument: HSGC#1


Operator: Justin Knoy

Column: DB-ALC1


Location: Vial 26

Method: C:\HPCHEM\1\METHODS\SIMALC1.M

Sample Info:




#	Compound	Peak Area	RT (min)
1	Ethanol	1517	1.085
2	n-Propanol	2661	1.765



Ethanol 0.126 g/100mL

JK



n-Propanol 0.012 g/100mL

JK

Washington State Patrol Toxicology Laboratory
 2203 Airport Way S Seattle, WA 98134

Inj. Date: 4/12/2017 12:29:00 PM

Sample Name: 17036-4

Instrument: HSGC#1


Operator: Justin Knoy

Column: DB-ALC1


Location: Vial 27

Method: C:\HPCHEM\1\METHODS\SIMALC1.M

Sample Info:




#	Compound	Peak Area	RT (min)
1	Ethanol	1517	1.084
2	n-Propanol	2677	1.763



Ethanol 0.125 g/100mL

JK



n-Propanol 0.012 g/100mL

JK

Washington State Patrol Toxicology Laboratory
2203 Airport Way S Seattle, WA 98134

Inj. Date: 4/12/2017 12:32:13 PM

Sample Name: 17036-5

Instrument: HSGC#1


Operator: Justin Knoy

Column: DB-ALC1


Location: Vial 28

Method: C:\HPCHEM\1\METHODS\SIMALC1.M

Sample Info:




#	Compound	Peak Area	RT (min)
1	Ethanol	1492	1.084
2	n-Propanol	2635	1.763



Ethanol 0.125 g/100mL

JK



n-Propanol 0.012 g/100mL

JK

Washington State Patrol Toxicology Laboratory
 2203 Airport Way S Seattle, WA 98134

Inj. Date: 4/12/2017 12:35:26 PM

Sample Name: 0.10 CTRL

Instrument: HSGC#1


Operator: Justin Knoy

Column: DB-ALC1


Location: Vial 29

Method: C:\HPCHEM\1\METHODS\SIMALC1.M

Sample Info: 17036




#	Compound	Peak Area	RT (min)
1	Ethanol	1242	1.085
2	n-Propanol	2709	1.764



Ethanol 0.101 g/100mL

JK



n-Propanol 0.012 g/100mL

JK

Washington State Patrol Toxicology Laboratory
 2203 Airport Way S Seattle, WA 98134

Inj. Date: 4/12/2017 12:38:40 PM

Sample Name: NEG CTRL

Instrument: HSGC#1


Operator: Justin Knoy

Column: DB-ALC1


Location: Vial 30

Method: C:\HPCHEM\1\METHODS\SIMALC1.M

Sample Info: 17036




#	Compound	Peak Area	RT (min)
1	Ethanol	0	0.000
2	n-Propanol	2673	1.762



Ethanol 0.000 g/100mL

POT



n-Propanol 0.012 g/100mL

JK

Sequence Parameters:

Operator: Andrew Gingras
 Data File Naming: Prefix/Counter
 Signal 1 Prefix: SIG1
 Counter: 0001
 Signal 2 Prefix: SIG2
 Counter: 0001
 Data Directory: C:\HPCHEM\1\DATA\
 Data Subdirectory: 170413AG
 Part of Methods to run: According to Runtime Checklist
 Barcode Reader: not used
 Shutdown Cmd/Macro: none

Sequence Comment:

CAL 1 (0.079g/100mL) - LOT# E0217-01 - EXP 8/21/2017
 CAL 2 (0.158g/100mL) - LOT# E0217-02 - EXP 8/21/2017
 CAL 3 (0.316g/100mL) - LOT# E0217-03 - EXP 8/21/2017

n-Propanol ISTD - LOT# P0117 - 4/20/2017
 CTRL 1 (0.04g/100mL) - LOT# FN12181501 - EXP 12/2020
 CTRL 2 (0.10g/100mL) - LOT# FN08051301 - EXP 10/2018
 CTRL 3 (0.20g/100mL) - LOT# FN08101505 - EXP 2/2021

Calibrators and controls filed with 17034
 Dilutor #3

AG
 4/13/17

Sequence Table (Front Injector):

Method and Injection Info Part:

Line	Location	SampleName	Method	Inj	SampleType	InjVolume	DataFile
1	Vial 1	BLANK	SIMALC1	1	Sample		
2	Vial 2	0.079 CAL 1	SIMALC1	1	Calib		
3	Vial 3	0.158 CAL 2	SIMALC1	1	Calib		
4	Vial 4	0.316 CAL 3	SIMALC1	1	Calib		
5	Vial 5	NEG CTRL	SIMALC1	1	Ctrl Samp		
6	Vial 6	0.04 CTRL	SIMALC1	1	Ctrl Samp		
7	Vial 7	0.10 CTRL	SIMALC1	1	Ctrl Samp		
8	Vial 8	0.20 CTRL	SIMALC1	1	Ctrl Samp		
9	Vial 9	NEG CTRL	SIMALC1	1	Ctrl Samp		
10	Vial 10	QAP 17034 #1	SIMALC1	1	Sample		
11	Vial 11	QAP 17034 #2	SIMALC1	1	Sample		
12	Vial 12	QAP 17034 #3	SIMALC1	1	Sample		
13	Vial 13	QAP 17034 #4	SIMALC1	1	Sample		
14	Vial 14	QAP 17034 #5	SIMALC1	1	Sample		
15	Vial 15	0.10 CTRL	SIMALC1	1	Ctrl Samp		
16	Vial 16	NEG CTRL	SIMALC1	1	Ctrl Samp		
17	Vial 17	QAP 17035 #1	SIMALC1	1	Sample		
18	Vial 18	QAP 17035 #2	SIMALC1	1	Sample		
19	Vial 19	QAP 17035 #3	SIMALC1	1	Sample		
20	Vial 20	QAP 17035 #4	SIMALC1	1	Sample		
21	Vial 21	QAP 17035 #5	SIMALC1	1	Sample		
22	Vial 22	0.10 CTRL	SIMALC1	1	Ctrl Samp		
23	Vial 23	NEG CTRL	SIMALC1	1	Ctrl Samp		
24	Vial 24	QAP 17036 #1	SIMALC1	1	Sample		

~~17036~~ *AG*
 4/13/17

17036

AG

Line	Location	SampleName	Method	Inj	SampleType	InjVolume	DataFile
25	Vial 25	QAP 17036 #2	SIMALC1	1	Sample		
26	Vial 26	QAP 17036 #3	SIMALC1	1	Sample		
27	Vial 27	QAP 17036 #4	SIMALC1	1	Sample		
28	Vial 28	QAP 17036 #5	SIMALC1	1	Sample		
29	Vial 29	0.10 CTRL	SIMALC1	1	Ctrl Samp		
30	Vial 30	NEG CTRL	SIMALC1	1	Ctrl Samp		
31	Vial 31	QAP 17037 #1	SIMALC1	1	Sample		
32	Vial 32	QAP 17037 #2	SIMALC1	1	Sample		
33	Vial 33	QAP 17037 #3	SIMALC1	1	Sample		
34	Vial 34	QAP 17037 #4	SIMALC1	1	Sample		
35	Vial 35	QAP 17037 #5	SIMALC1	1	Sample		
36	Vial 36	0.10 CTRL	SIMALC1	1	Ctrl Samp		
37	Vial 37	NEG CTRL	SIMALC1	1	Ctrl Samp		
38	Vial 38	QAP 17038 #1	SIMALC1	1	Sample		
39	Vial 39	QAP 17038 #2	SIMALC1	1	Sample		
40	Vial 40	QAP 17038 #3	SIMALC1	1	Sample		
41	Vial 41	QAP 17038 #4	SIMALC1	1	Sample		
42	Vial 42	QAP 17038 #5	SIMALC1	1	Sample		
43	Vial 43	0.10 CTRL	SIMALC1	1	Ctrl Samp		
44	Vial 44	NEG CTRL	SIMALC1	1	Ctrl Samp		

Calibration Part:

Line	Location	SampleName	Method	CalLev	Update	RF	Update	RT	Interval
2	Vial 2	0.079 CAL 1	SIMALC1	1	Replace		Replace		
3	Vial 3	0.158 CAL 2	SIMALC1	2	Replace		Replace		
4	Vial 4	0.316 CAL 3	SIMALC1	3	Replace		Replace		

Sequence Table (Back Injector):

No entries - empty table!

17036 BT
4/19/17

AG

Washington State Patrol Toxicology Laboratory
2203 Airport Way S Seattle, WA 98134

Inj. Date: 4/13/2017 3:47:05 PM

Sample Name: QAP 17036 #1

Instrument: HSGC#1


Operator: Andrew Gingras

Column: DB-ALC1


Location: Vial 24

Method: C:\HPCHEM\1\METHODS\SIMALC1.M

Sample Info:




#	Compound	Peak Area	RT (min)
1	Ethanol	1532	1.086
2	n-Propanol	2850	1.765



Ethanol 0.125 g/100mL

AK



n-Propanol 0.012 g/100mL

AK


Washington State Patrol Toxicology Laboratory
 2203 Airport Way S Seattle, WA 98134

Inj. Date: 4/13/2017 3:50:18 PM
 Instrument: HSGC#1
 Column: DB-ALC1


Sample Name: QAP 17036 #2
 Operator: Andrew Gingras
 Location: Vial 25

Method: C:\HPCHEM\1\METHODS\SIMALC1.M

Sample Info:




#	Compound	Peak Area	RT (min)
1	Ethanol	1585	1.085
2	n-Propanol	2961	1.764



Ethanol 0.125 g/100mL

Handwritten initials



n-Propanol 0.012 g/100mL

Handwritten initials

Washington State Patrol Toxicology Laboratory
2203 Airport Way S Seattle, WA 98134

Inj. Date: 4/13/2017 3:53:32 PM

Sample Name: QAP 17036 #3

Instrument: HSGC#1


Operator: Andrew Gingras

Column: DB-ALC1


Location: Vial 26

Method: C:\HPCHEM\1\METHODS\SIMALC1.M

Sample Info:




#	Compound	Peak Area	RT (min)
1	Ethanol	1557	1.084
2	n-Propanol	2914	1.763



Ethanol 0.125 g/100mL

pat



n-Propanol 0.012 g/100mL


JB

Washington State Patrol Toxicology Laboratory
2203 Airport Way S Seattle, WA 98134


Inj. Date: 4/13/2017 3:56:45 PM
Instrument: HSGC#1
Column: DB-ALC1
Method: C:\HPCHEM\1\METHODS\SIMALC1.M

Sample Name: QAP 17036 #4
Operator: Andrew Gingras
Location: Vial 27


Sample Info:



#	Compound	Peak Area	RT (min)
1	Ethanol	1498	1.084
2	n-Propanol	2789	1.764



Ethanol 0.125 g/100mL *MA*



n-Propanol 0.012 g/100mL


AG

Washington State Patrol Toxicology Laboratory
 2203 Airport Way S Seattle, WA 98134


Inj. Date: 4/13/2017 3:59:58 PM
 Instrument: HSGC#1
 Column: DB-ALC1
 Method: C:\HPCHEM\1\METHODS\SIMALC1.M

Sample Name: QAP 17036 #5
 Operator: Andrew Gingras
 Location: Vial 28

Sample Info:




#	Compound	Peak Area	RT (min)
1	Ethanol	1500	1.084
2	n-Propanol	2794	1.764



Ethanol 0.125 g/100mL

MT



n-Propanol 0.012 g/100mL

AG


Washington State Patrol Toxicology Laboratory
 2203 Airport Way S Seattle, WA 98134

Inj. Date: 4/13/2017 4:03:12 PM
 Instrument: HSGC#1


Sample Name: 0.10 CTRL
 Operator: Andrew Gingras
 Location: Vial 29

Column: DB-ALC1
 Method: C:\HPCHEM\1\METHODS\SIMALC1.M

Sample Info: 17036




#	Compound	Peak Area	RT (min)
1	Ethanol	1283	1.085
2	n-Propanol	2941	1.764



Ethanol 0.102 g/100mL

AG



n-Propanol 0.012 g/100mL

AG

Washington State Patrol Toxicology Laboratory
2203 Airport Way S Seattle, WA 98134

Inj. Date: 4/13/2017 4:06:25 PM

Sample Name: NEG CTRL

Instrument: HSGC#1


Operator: Andrew Gingras

Column: DB-ALC1


Location: Vial 30

Method: C:\HPCHEM\1\METHODS\SIMALC1.M

Sample Info: 17036




#	Compound	Peak Area	RT (min)
1	Ethanol	0	0.000
2	n-Propanol	2906	1.764



Ethanol 0.000 g/100mL

Handwritten mark



n-Propanol 0.012 g/100mL

Handwritten signature

Sequence Parameters:

Operator: Christie Mitchell-Mata

Data File Naming: Prefix/Counter

Signal 1 Prefix: SIG1

Counter: 0001

Signal 2 Prefix: SIG2

Counter: 0001

Data Directory: C:\HPCHEM\1\DATA\

Data Subdirectory: 170417CM

Part of Methods to run: According to Runtime Checklist

Barcode Reader: not used

Shutdown Cmd/Macro: none

Sequence Comment:

Ethanol Calibrator 1 0.079 g/100 mL, E0217-01 - Exp. 08/21/17

Ethanol Calibrator 2 0.158 g/100 mL, E0217-02 - Exp. 08/21/17

Ethanol Calibrator 3 0.316 g/100 mL, E0217-03 - Exp. 08/21/17

0.04 Control - Lot #FN12181501 - Exp. 12/2020

0.10 Control - Lot #FN08051301 - Exp. 10/2018

0.20 Control - Lot #FN08101505 - Exp. 02/2021

ISTD Lot#P0317 - Exp. 06/13/2017

Dilutor #1

Calibration 1-9 filed with 17036

Sequence Table (Front Injector):

Method and Injection Info Part:

Line	Location	SampleName	Method	Inj	SampleType	InjVolume	DataFile
1	Vial 1	BLANK	SIMALC1	1	Sample		
2	Vial 2	0.079 CAL 1	SIMALC1	1	Calib		
3	Vial 3	0.158 CAL 2	SIMALC1	1	Calib		
4	Vial 4	0.316 CAL 3	SIMALC1	1	Calib		
5	Vial 5	Negative CTRL	SIMALC1	1	Ctrl Samp		
6	Vial 6	0.04 CTRL	SIMALC1	1	Ctrl Samp		
7	Vial 7	0.10 CTRL	SIMALC1	1	Ctrl Samp		
8	Vial 8	0.20 CTRL	SIMALC1	1	Ctrl Samp		
9	Vial 9	Negative CTRL	SIMALC1	1	Ctrl Samp		
10	Vial 10	17036 #1	SIMALC1	1	Sample		
11	Vial 11	17036 #2	SIMALC1	1	Sample		
12	Vial 12	17036 #3	SIMALC1	1	Sample		
13	Vial 13	17036 #4	SIMALC1	1	Sample		
14	Vial 14	17036 #5	SIMALC1	1	Sample		
15	Vial 15	0.10 CTRL	SIMALC1	1	Ctrl Samp		
16	Vial 16	Negative CTRL	SIMALC1	1	Ctrl Samp		

17036
BT
Hala

Calibration Part:

Line	Location	SampleName	Method	CalLev	Update	RF	Update	RT	Interval
2	Vial 2	0.079 CAL 1	SIMALC1	1	Replace		Replace		
3	Vial 3	0.158 CAL 2	SIMALC1	2	Replace		Replace		
4	Vial 4	0.316 CAL 3	SIMALC1	3	Replace		Replace		

Handwritten signature

Sequence: C:\HPCHEM\1\SEQUENCE\CMQAP2.S

Sequence Table (Back Injector):

No entries - empty table!

17036

BT
4/19/17

W

=====
Calibration Table
=====

Calib. Data Modified : Monday, April 17, 2017 8:10:03 AM
Calculate : Internal Standard
Based on : Peak Area
Rel. Reference Window : 5.000 %
Abs. Reference Window : 0.050 min
Rel. Non-ref. Window : 5.000 %
Abs. Non-ref. Window : 0.050 min
Multiplier : 1.0000
Dilution : 1.0000
Sample Amount : 0.00000
Use Multiplier & Dilution Factor with ISTDs
Uncalibrated Peaks : not reported
Partial Calibration : No recalibration if peaks missing
Curve Type : Linear
Origin : Included
Weight : Equal
Recalibration Settings:
Average Response : No Update
Average Retention Time: No Update

Calibration Report Options :
Printout of recalibrations within a sequence:
Normal Report after Recalibration

Sample ISTD Information:

ISTD ISTD Amount Name
[g/100mL]
-----|-----|-----
1 1.20000e-2 n-Propanol

Signal 1: FID1 A,

RetTime [min]	Lvl Sig	Amount [g/100mL]	Area	Amt/Area	Ref Grp Name
1.084	1 1	7.91500e-2	979.27411	8.08252e-5	1 Ethanol
	2	1.58300e-1	2076.62354	7.62295e-5	
	3	3.19520e-1	3773.46948	8.46754e-5	
1.764	1 1	1.20000e-2	2665.24585	4.50240e-6	I1 n-Propanol
	2	1.20000e-2	2714.42358	4.42083e-6	
	3	1.20000e-2	2601.23340	4.61320e-6	


=====
Peak Sum Table
=====

No Entries in table
=====


17036
BT
4/17/17

AM

=====
Calibration Curves
=====



Ethanol at exp. RT: 1.084
FID1 A,
Correlation: 0.99936
Residual Std. Dev.: 0.02711
Formula: $y = mx + b$
m: 5.45664e-2
b: 1.26122e-2
x: Amount Ratio
y: Area Ratio



n-Propanol at exp. RT: 1.764
FID1 A,
Correlation: 1.00000
Residual Std. Dev.: 0.00000
Formula: $y = mx + b$
m: 1.00000
b: 0.00000
x: Amount Ratio
y: Area Ratio

17036

BT
4/19/17

W

Washington State Patrol Toxicology Laboratory
2203 Airport Way S Seattle, WA 98134

Inj. Date: 4/17/2017 7:57:57 AM

Sample Name: BLANK

Instrument: HSGC#1


Operator: Christie Mitchell-Mata

Column: DB-ALC1


Location: Vial 1

Method: C:\HPCHEM\1\METHODS\SIMALC1.M


Sample Info: 17036



#	Compound	Peak Area	RT (min)
1	Ethanol	0	0.000
2	n-Propanol	0	0.000



Ethanol 0.000 g/100mL *BT*



n-Propanol 0.000 g/100mL

BT

Washington State Patrol Toxicology Laboratory
 2203 Airport Way S Seattle, WA 98134

Inj. Date: 4/17/2017 8:01:16 AM

Sample Name: 0.079 CAL 1

Instrument: HSGC#1


Operator: Christie Mitchell-Mata

Column: DB-ALC1


Location: Vial 2

Method: C:\HPCHEM\1\METHODS\SIMALC1.M

Sample Info: 17036




#	Compound	Peak Area	RT (min)
1	Ethanol	979	1.085
2	n-Propanol	2665	1.765



Ethanol 0.078 g/100mL

Handwritten mark



n-Propanol 0.012 g/100mL

Handwritten mark

Washington State Patrol Toxicology Laboratory
 2203 Airport Way S Seattle, WA 98134

Inj. Date: 4/17/2017 8:04:32 AM

Sample Name: 0.158 CAL 2

Instrument: HSGC#1


Operator: Christie Mitchell-Mata

Column: DB-ALC1


Location: Vial 3

Method: C:\HPCHEM\1\METHODS\SIMALC1.M

Sample Info: 17036




#	Compound	Peak Area	RT (min)
1	Ethanol	2077	1.087
2	n-Propanol	2714	1.766



Ethanol 0.165 g/100mL

BA




n-Propanol 0.012 g/100mL


AM

Washington State Patrol Toxicology Laboratory
2203 Airport Way S Seattle, WA 98134

Inj. Date: 4/17/2017 8:07:49 AM Sample Name: 0.316 CAL 3
Instrument: HSGC#1 Operator: Christie Mitchell-Mata
Column: DB-ALC1 Location: Vial 4
Method: C:\HPCHEM\1\METHODS\SIMALC1.M
Sample Info: 17036




#	Compound	Peak Area	RT (min)
1	Ethanol	3773	1.084
2	n-Propanol	2601	1.764



Ethanol 0.316 g/100mL

PK



n-Propanol 0.012 g/100mL

W

Washington State Patrol Toxicology Laboratory
2203 Airport Way S Seattle, WA 98134

Inj. Date: 4/17/2017 8:11:04 AM

Sample Name: Negative CTRL

Instrument: HSGC#1


Operator: Christie Mitchell-Mata

Column: DB-ALC1


Location: Vial 5

Method: C:\HPCHEM\1\METHODS\SIMALC1.M

Sample Info: 17036




#	Compound	Peak Area	RT (min)
1	Ethanol	0	0.000
2	n-Propanol	2643	1.763



Ethanol 0.000 g/100mL

WAT



n-Propanol 0.012 g/100mL

W

Washington State Patrol Toxicology Laboratory
 2203 Airport Way S Seattle, WA 98134

Inj. Date: 4/17/2017 8:14:16 AM

Sample Name: 0.04 CTRL

Instrument: HSGC#1


Operator: Christie Mitchell-Mata

Column: DB-ALC1


Location: Vial 6

Method: C:\HPCHEM\1\METHODS\SIMALC1.M

Sample Info: 17036




#	Compound	Peak Area	RT (min)
1	Ethanol	496	1.084
2	n-Propanol	2666	1.763



Ethanol 0.038 g/100mL

BT



n-Propanol 0.012 g/100mL

LM

Washington State Patrol Toxicology Laboratory
 2203 Airport Way S Seattle, WA 98134

Inj. Date: 4/17/2017 8:17:29 AM

Sample Name: 0.10 CTRL

Instrument: HSGC#1


Operator: Christie Mitchell-Mata

Column: DB-ALC1


Location: Vial 7

Method: C:\HPCHEM\1\METHODS\SIMALC1.M

Sample Info: 17036




#	Compound	Peak Area	RT (min)
1	Ethanol	1225	1.085
2	n-Propanol	2648	1.765



Ethanol 0.099 g/100mL

Handwritten mark



n-Propanol 0.012 g/100mL

Handwritten mark

Washington State Patrol Toxicology Laboratory
 2203 Airport Way S Seattle, WA 98134

Inj. Date: 4/17/2017 8:20:42 AM

Sample Name: 0.20 CTRL

Instrument: HSGC#1


Operator: Christie Mitchell-Mata

Column: DB-ALC1


Location: Vial 8

Method: C:\HPCHEM\1\METHODS\SIMALC1.M

Sample Info: 17036




#	Compound	Peak Area	RT (min)
1	Ethanol	2327	1.083
2	n-Propanol	2560	1.762



Ethanol 0.197 g/100mL

MM



n-Propanol 0.012 g/100mL

MM

Washington State Patrol Toxicology Laboratory
 2203 Airport Way S Seattle, WA 98134

Inj. Date: 4/17/2017 8:23:56 AM

Sample Name: Negative CTRL

Instrument: HSGC#1


Operator: Christie Mitchell-Mata

Column: DB-ALC1


Location: Vial 9

Method: C:\HPCHEM\1\METHODS\SIMALC1.M

Sample Info: 17036




#	Compound	Peak Area	RT (min)
1	Ethanol	0	0.000
2	n-Propanol	2652	1.762



Ethanol 0.000 g/100mL

PM



n-Propanol 0.012 g/100mL

W


Washington State Patrol Toxicology Laboratory
2203 Airport Way S Seattle, WA 98134

Inj. Date: 4/17/2017 8:27:09 AM
Instrument: HSGC#1
Column: DB-ALC1


Sample Name: 17036 #1
Operator: Christie Mitchell-Mata
Location: Vial 10

Method: C:\HPCHEM\1\METHODS\SIMALC1.M

Sample Info:




#	Compound	Peak Area	RT (min)
1	Ethanol	1664	1.089
2	n-Propanol	2830	1.769



Ethanol 0.127 g/100mL

mt



n-Propanol 0.012 g/100mL

cm

Washington State Patrol Toxicology Laboratory
 2203 Airport Way S Seattle, WA 98134

Inj. Date: 4/17/2017 8:30:22 AM

Sample Name: 17036 #2

Instrument: HSGC#1


Operator: Christie Mitchell-Mata

Column: DB-ALC1


Location: Vial 11

Method: C:\HPCHEM\1\METHODS\SIMALC1.M

Sample Info:




#	Compound	Peak Area	RT (min)
1	Ethanol	1500	1.084
2	n-Propanol	2626	1.764



Ethanol 0.123 g/100mL

MT



n-Propanol 0.012 g/100mL

CM

Washington State Patrol Toxicology Laboratory
2203 Airport Way S Seattle, WA 98134

Inj. Date: 4/17/2017 8:33:35 AM

Sample Name: 17036 #3

Instrument: HSGC#1


Operator: Christie Mitchell-Mata

Column: DB-ALC1


Location: Vial 12

Method: C:\HPCHEM\1\METHODS\SIMALC1.M

Sample Info:




#	Compound	Peak Area	RT (min)
1	Ethanol	1489	1.084
2	n-Propanol	2624	1.763



Ethanol 0.122 g/100mL

pat



n-Propanol 0.012 g/100mL

cm


Washington State Patrol Toxicology Laboratory
2203 Airport Way S Seattle, WA 98134

Inj. Date: 4/17/2017 8:36:49 AM
Instrument: HSGC#1
Column: DB-ALC1


Sample Name: 17036 #4
Operator: Christie Mitchell-Mata
Location: Vial 13

Method: C:\HPCHEM\1\METHODS\SIMALC1.M

Sample Info:




#	Compound	Peak Area	RT (min)
1	Ethanol	1514	1.084
2	n-Propanol	2640	1.763



Ethanol 0.123 g/100mL

Handwritten initials



n-Propanol 0.012 g/100mL

Handwritten initials

Washington State Patrol Toxicology Laboratory
2203 Airport Way S Seattle, WA 98134

Inj. Date: 4/17/2017 8:40:02 AM

Sample Name: 17036 #5

Instrument: HSGC#1


Operator: Christie Mitchell-Mata

Column: DB-ALC1


Location: Vial 14

Method: C:\HPCHEM\1\METHODS\SIMALC1.M

Sample Info:




#	Compound	Peak Area	RT (min)
1	Ethanol	1563	1.085
2	n-Propanol	2665	1.765



Ethanol 0.126 g/100mL

Handwritten signature



n-Propanol 0.012 g/100mL

Handwritten signature

Washington State Patrol Toxicology Laboratory
 2203 Airport Way S Seattle, WA 98134

Inj. Date: 4/17/2017 8:43:15 AM

Sample Name: 0.10 CTRL

Instrument: HSGC#1


Operator: Christie Mitchell-Mata

Column: DB-ALC1


Location: Vial 15

Method: C:\HPCHEM\1\METHODS\SIMALC1.M

Sample Info: 17036




#	Compound	Peak Area	RT (min)
1	Ethanol	1230	1.084
2	n-Propanol	2666	1.764



Ethanol 0.099 g/100mL

Handwritten mark



n-Propanol 0.012 g/100mL

Handwritten mark

Washington State Patrol Toxicology Laboratory
2203 Airport Way S Seattle, WA 98134

Inj. Date: 4/17/2017 8:46:28 AM

Sample Name: Negative CTRL

Instrument: HSGC#1


Operator: Christie Mitchell-Mata

Column: DB-ALC1


Location: Vial 16

Method: C:\HPCHEM\1\METHODS\SIMALC1.M

Sample Info: 17036




#	Compound	Peak Area	RT (min)
1	Ethanol	0	0.000
2	n-Propanol	2760	1.768



Ethanol 0.000 g/100mL

BM



n-Propanol 0.012 g/100mL

BM